
Labber

www.keysight.com/find/labber

Notices

Copyright Notice

© Keysight Technologies 2022

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Keysight Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

9018-18184

Edition

November 2022

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Declaration of Conformity

Declarations of Conformity for this product and for other Keysight products may be downloaded from the Web. Go to <http://www.keysight.com/go/conformity> and click on “Declarations of Conformity.” You can then search by product number to find the latest Declaration of Conformity.

U.S. Government Rights

The Software is “commercial computer software,” as defined by Federal Acquisition Regulation (“FAR”) 2.101. Pursuant to FAR 12.212 and 27.405-3 and Department of Defense FAR Supplement (“DFARS”) 227.7202, the U.S. government acquires commercial computer software under the same terms by which the software is customarily provided to the public. Accordingly, Keysight provides the Software to U.S. government customers under its standard commercial license, which is embodied in its End User License Agreement (EULA), a copy of which can be found at <http://www.keysight.com/find/sweula>. The license set forth in the EULA represents the exclusive authority by which the U.S. government may use, modify, distribute, or disclose the Software. The EULA and the license set forth therein, does not require or permit, among other things, that Keysight: (1) Furnish technical information related to commercial computer software or commercial computer software documentation that is not customarily provided to the public; or (2) Relinquish to, or otherwise provide, the government rights in excess of these rights customarily provided to the public to use, modify, reproduce, release, perform, display, or disclose commercial computer software or commercial computer software documentation. No additional government requirements beyond those set forth in the EULA shall apply, except to the extent that those terms, rights, or licenses are explicitly required from all providers of commercial computer software pursuant to the FAR and the DFARS and are set forth specifically in writing elsewhere in the EULA. Keysight shall be under no obligation to update, revise or otherwise modify the Software. With respect to any technical data as defined by FAR 2.101, pursuant to FAR 12.211 and 27.404.2 and DFARS 227.7102, the U.S. government acquires no greater than Limited Rights as defined in FAR 27.401 or DFARS 227.7103-5 (c), as applicable in any technical data.

Warranty

THE MATERIAL CONTAINED IN THIS DOCUMENT IS PROVIDED “AS IS,” AND IS SUBJECT TO BEING CHANGED, WITHOUT NOTICE, IN FUTURE EDITIONS. FURTHER, TO THE MAXIMUM EXTENT PERMITTED BY APPLICABLE LAW, KEYSIGHT DISCLAIMS ALL WARRANTIES, EITHER EXPRESS OR IMPLIED, WITH REGARD TO THIS MANUAL AND ANY INFORMATION CONTAINED HEREIN, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. KEYSIGHT SHALL NOT BE LIABLE FOR ERRORS OR FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH THE FURNISHING, USE, OR PERFORMANCE OF THIS DOCUMENT OR OF ANY INFORMATION CONTAINED HEREIN. SHOULD KEYSIGHT AND THE USER HAVE A SEPARATE WRITTEN AGREEMENT WITH WARRANTY TERMS COVERING THE MATERIAL IN THIS DOCUMENT THAT CONFLICT WITH THESE TERMS, THE WARRANTY TERMS IN THE SEPARATE AGREEMENT SHALL CONTROL.

Safety Information

CAUTION

A CAUTION notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a CAUTION notice until the indicated conditions are fully understood and met.

WARNING

A WARNING notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

Table of Contents

1	Introduction	8
2	Installation	9
2.1	Installation - Microsoft Windows.....	9
2.1.1	Microsoft Windows - Troubleshooting	9
2.1.2	Microsoft Windows - Defender SmartScreen warnings	10
2.2	Installing a Labber License	10
2.2.1	Installing a Local License	10
2.2.3	Installing a Floating License	11
2.3	VISA distribution	13
2.4	Network and firewall settings.....	13
2.5	Program folders	13
2.5.1	Data folder	13
2.5.2	Instrument drivers	13
2.5.3	Scripting	14
3	Instrument Server	15
3.1	Program startup	15
3.2	Server window	16
3.3	Adding instruments.....	16
3.4	Configuring instruments	17
3.5	Instruments with vector-valued quantities	19
3.6	Keeping track of open client connections.....	19
3.7	Troubleshooting - Timing statistics	20
3.8	Troubleshooting - Instrument and Network logs.....	20
4	Controller	21
4.1	Controller operation	21
4.2	Improving controller performance	22
5	Scheduler	24
5.1	Scheduling measurements.....	24
5.2	Scheduler settings.....	26
6	Measurement program.....	27
6.1	Measurement configuration	27
6.2	Adding channels.....	28

6.2.1	Instrument configuration - locks	29
6.3	Sending and retrieving values from instruments.....	30
6.4	Defining step sequences	30
6.4.1	Step setup - Basic settings	30
6.4.2	Step setup - Advanced settings.....	31
6.4.3	Step setup - Sweep mode	33
6.4.4	Step setup - Channel relations	34
6.5	Log channels.....	35
6.5.1	Log channels limits	35
6.6	Timing.....	36
6.7	Log name, Project and User tags, Comments	36
6.8	Tags	37
6.9	Starting a measurement	37
6.10	Signal connections	40
6.11	Hardware timing and synchronization.....	43
6.11.1	Arm/trig mode	44
6.11.2	Hardware looping.....	44
6.12	File locks.....	45
6.13	Measurement settings	46
6.13.1	General.....	46
6.13.2	Optimizer	46
6.14	Comparing Measurement configurations.....	46
7	Optimizer	48
7.1	Optimizer operation.....	48
7.1.1	Cost function	48
7.1.2	Termination and convergence criteria.....	49
7.1.3	Running an optimizer measurement	50
7.2	Optimizer settings.....	50
7.2.1	General optimizer settings.....	50
7.2.2	Individual parameter settings	51
7.3	Custom optimizers	51
7.3.1	Defining custom optimizers	51
7.3.2	Defining optimizers settings	52

7.3.3	Using custom optimizers.....	54
8	Log Browser	55
8.1	Database	55
8.2	Log browser dialog.....	55
8.2.1	Database hierarchy view	56
8.2.2	Log list	56
8.2.3	Graph / Log info	57
8.2.4	Tool bar	57
9	Log Viewer.....	59
9.1	Plot config	59
9.1.1	Equations	62
9.1.2	Physical vs. Instrument units	63
9.2	Entry list	63
9.3	Tool bar	63
9.4	Multi-panel graph mode	64
9.5	Image mode	65
9.6	Views	67
9.7	Exporting data	67
9.7.1	Exporting to Image	68
9.7.2	Exporting to Text	68
9.7.3	Exporting to Matlab	69
9.7.4	Custom Export.....	69
10	Preferences	71
10.1	Folders.....	71
10.2	Server	72
10.3	Measurement.....	74
10.4	Log Viewer.....	75
10.5	Logger.....	76
10.6	Advanced.....	76
11	Scripting	77
11.1	Console options.....	78
11.2	Scripting using Python.....	79
12	Instrument drivers	80

12.1	Driver definition files.....	81
12.1.1	Signal Generators and Signal Analyzers	82
12.1.2	General settings	82
12.1.3	Model and options	84
12.1.4	VISA Settings	85
12.2	Quantities.....	88
12.3	Custom drivers - Python code.....	92
12.3.1	Creating the driver definition file.....	93
12.3.2	Implementing the <i>Python</i> code.....	94
12.3.3	Helper functions for <code>quant</code> objects.....	97
12.3.4	Helper functions for <code>driver</code> objects	98
12.3.5	Testing the driver	102
12.4	Subclassing the <i>VISA</i> driver	103
12.4.1	Helper functions for drivers subclassing the <i>VISA</i> driver	104
12.5	Support for sweeping.....	105
12.5.1	1 Sweeping - Driver definition file.....	105
12.5.2	Sweeping - Python code.....	107
12.6	Controller drivers	108
12.7	Hardware arming and triggering.....	108
12.8	Hardware looping	108
12.8.1	Hardware looping - outputting values	109
12.8.2	Hardware looping - reading values	109
12.9	Python distribution	109
12.9.1	Python distribution, 32-bit version	110
12.9.2	External Python distribution	110
12.9.3	Troubleshooting, external Python distribution.....	111
Appendix A:	Python API	113
A.1	Installation	113
A.1.1	Installing the API with pip	113
A.1.2	Requirements.....	114
A.1.3	Testing the API	114
A.1.4	Upgrading from earlier versions	114
A.2	Instrument server	115

A.2.1	Labber client.....	115
A.2.2	Scheduling measurements.....	115
A.2.3	Connecting to instruments	116
A.2.4	Blocking vs. non-blocking clients	117
A.2.5	Function definitions	119
A.2.6	Class definitions	120
A.3	Log files	133
A.3.1	Reading data from Labber	133
A.3.2	Creating Labber log files.....	134
A.3.3	Function definitions	135
A.3.4	LogFile class.....	138
A.4	Script tools	146
A.4.1	Initialization.....	146
A.4.2	Example.....	146
A.4.3	Function definitions	147
A.4.4	MeasurementObject class	148
A.5	Configurations.....	152
A.5.1	Example.....	153
A.5.2	Scenario class	156
A.5.3	Scenario module	162
A.5.4	Instrument module	166
A.5.5	Step module	169
A.5.6	Lookup module	174
B	Appendix B: Labber Quantum User Guide.....	177
B.1	System requirements	177
B.1.1	Recommended Standard Hardware Configurations.....	177
B.1.2	Labber Software Driver Packages and Dependencies.....	182
B.1.3	M5400xxxA Software/Gateware Compatibility and Dependencies	183
B.1.4	Common Software Compatibility.....	183
B.1.5	Instrument Hardware and Firmware	184
B.2	Installation Guide and Troubleshooting	184
B.2.1	Quantum IP Library	184
B.2.2	Configuring a Python Environment for Labber Drivers.....	185

B.2.3	Multi-chassis configuration	187
B.2.4	PXI Trigger Reservations	189
B.2.5	Common Problems When Getting Started with TSE/Quantum IP Library	190
B.3	Keysight PXI Labber Drivers	197
B.3.1	Keysight PXI AWG.....	197
B.3.2	Keysight PXI PWTSE Trigger	201
B.3.3	Keysight PXI Digitizer Demod	207
B.3.4	Multi-Qubit Pulse Generator for Agile sequences	239
B.3.5	Keysight PXI Agile AWG.....	258
B.3.6	Keysight PXI SMU	263

1 Introduction

The software package consists of three separate programs. The *Instrument Server* handles the communication with the instruments, the *Measurement* program allows instrument values to be controlled and measured in user-defined sequences, while the *Log Browser* is used to organize and analyze the acquired data. The relation between the parts is visualized in Fig. 1.1.




Fig 1.1. Overview and structure of the components in the Labber software package.

In a typical experimental setup, the *Instrument Server* keeps track of and communicates with all the instruments and equipment available in the setup. The communication can be over GPIB, serial, USB, TCPIP, or any other interface. During an experiment, the *Measurement* program will connect to an *Instrument Server* to output values to one specific instrument, or to read data from another one. Note that the *Measurement* program only talks to the *Instrument Server*, and not directly with the instruments. This modular approach allows the same generic procedure to be used for setting/reading values, regardless of the instrument type or the communication interface.

The *Measurement* program saves the experimental configuration, the instrument settings and the acquired data into a central log database. The *Log Browser* provides a fast and efficient method for browsing, visualizing and organizing the measured data. Finally,

the *Log Viewer* provides functionality for data analysis and for generating high-quality plots and figures.

In addition to the *Instrument Server*, *Measurement* and *Log Browser* programs, there is a Python API which allows all functionality to be accessed programmatically for scripting purposes or for writing custom applications.

2 Installation

2.1 Installation - Microsoft Windows

The setup file will install *Labber* to the default Microsoft Windows installation directory, as well as create folders for storing data and local driver files in the user's home directory. The default directories for local files can be set in the *Preferences* window, see Section [PrefsFolder](#)).

After installation, the *Instrument Server*, *Log Browser* and *Measurement* programs can be started by clicking the corresponding file from the Windows start menu. Note that the *Log Browser* and *Measurement* programs can be opened from within the *Instrument Server*, so it is usually sufficient to start just the server program.

The Windows installer will also attempt to install the Labber Python API. See section [A1. Installation](#).

Note that as of Labber 1.8 and greater, the default install location has changed to "C:\Program Files\Keysight\Labber". Unless you want to maintain side-by-side installations, it is recommended to uninstall the older version in "C:\Program Files\Labber". Shortcuts pinned to the Task Bar will also need to be updated to point to the new installation location. Start Menu shortcuts will be updated automatically.

2.1.1 Microsoft Windows - Troubleshooting

Depending on security settings, some virus scanners may prohibit *Labber* from being installed or run on your computer. If you're experiencing difficulties installing or running the program, try to temporarily disable the virus scanner.

Some Microsoft Windows distributions lack a few support files needed by the program to run correctly. If the program won't start, download and install the redistributable support files for Microsoft Visual C++ from <http://www.microsoft.com/en-us/download/details.aspx?id=26368>. Click on "Download" and select the file [vc_redist_x86.exe](#).

2.1.2 Microsoft Windows - Defender SmartScreen warnings

On certain Windows distributions, a dialog may pop up when installing *Labber* stating that the application is unrecognized and hasn't been screened by Microsoft. To install *Labber*, you need to override the dialog by clicking on "More info", and then click on the "Run anyway"-button.

2.2 Installing a Labber License

Labber is now licensed through the Keysight standard licensing system: PathWave License Manager. In order to register your license, you must download PathWave License Manager [here](#).


There are currently two options to install a license in your system:

- Install a local license (node locked, transportable, trial and USB portable)
- Use a floating license from a remote license server

Information on obtaining a license can be found [here](#).

2.2.1 Installing a Local License


- Once your license file is obtained, open PathWave License Manager and click on "+ Add a License File"



- Browse to the license file and select it. You should see a progress bar followed by a success message.
- The new license that was just installed should now appear in the corresponding tab (Local) in the **License Status** panel in PathWave License Manager. Note that the **Local Licenses** tab will list USB Portable and node-locked licenses.

2.2.2 Installing a Floating License

Floating licenses (network licenses) reside on a license server (a separate computer) and are checked out for use by Keysight products (in our case the Instrument Server), then returned (checked in) when no longer needed so that they can be used on another computer or instrument.




Setting up a license server


The process of setting up a license server is explained in more detail under Keysight PathWave License Manager Administrator's Guide, under *Floating License Setup: First Time*.

If a License Server is already set up and running, Client Machines can remotely connect to the license servers:

- Open PathWave License Manager and click on "Specify a Remote License Server"



- Choose "quantum" in Step 1 and enter the license server host name or IP address in Step 2



- The new license should now appear in the corresponding tab (Floating) in the **License Status** panel in PathWave License Manager.
- When the Instrument Server is launched, it should now check out a license from your specified floating license server. The status of the license can be checked by navigating to **Help > About Labber**. Here the IP address of the floating license server will be listed.

2.3 VISA distribution

To communicate with instruments over the VISA protocol, a VISA distribution needs to be installed on the computer. A VISA distribution can be downloaded from Keysight Technologies, www.keysight.com/find/iosuite.

2.4 Network and firewall settings

The *Instrument Server*, *Log Browser* and *Measurement* programs communicate using TCP/IP, which makes it possible to perform measurements involving instruments connected to different computers, even on different networks. The default settings assign TCP port 9406 for server/client communication and TCP port 9407 for sending internal notifications between the program parts. If you want to perform measurements in a multi-computer network and firewall is enabled on your system, the firewall must be configured to allow traffic on these ports. The port numbers can be changed in case they are occupied on your system (see Chapter [Prefs](#)).

2.5 Program folders

In addition to the folders with executables the program uses a few extra folder locations, as listed below.

2.5.1 Data folder

The program needs a folder for saving the measured data. By default, this folder is set to "<User home directory>/Labber/Data", but it can be changed at any time from the *Preferences* window (see Section [PrefsFolder](#)).

2.5.2 Instrument drivers

The program has two separate folders for storing instrument drivers, one main folder (set to be the "Drivers" subfolder under your installation) and one local folder (called "*Local drivers*" in the *Preferences*). The main driver folder resides in the same folder location as the executables, and should not be altered in a typical setup. The local driver folder is set to "<User home directory>/Labber/Drivers", but its location can be changed in the *Preferences*.

When creating a new instrument driver, the driver definition file should always be placed in the "*Local drivers*" folder. This allows the user's own drivers to be kept separately from the drivers provided with *Labber*, and it also prevents drivers written by users from being

deleted when updating the *Labber* program to a newer version. See Section [Drivers](#) for more information on creating instrument drivers.

2.5.3 Scripting

The Python API that contains scripting helper functions are located in the *Script* folder of the main program directory. See Section [scriptPython](#) for more information on scripting.

3 Instrument Server

3.1 Program startup

When starting the *Instrument Server*, the program will create a tray icon and a tray menu for controlling the server, showing preferences and launching the *Measurement* and the *Log Browser* programs (see Fig. 3.1). The tray menu also shows the status of the network server. Whenever the network server is running, clients are allowed to connect to the server to communicate with the instruments. Note that the network server keeps running in the background even after the server window has been closed. To stop the server, either select “*Stop Network Server*” from the tray menu or quit the server by selecting the “*Quit Server*” menu item.




Fig. 3.1 System tray menu for the *Instrument Server* program. In addition to controlling the server settings, the menu provides options for starting the *Measurement* and the *Log Browser* programs.




Fig. 3.2 The main Instrument Server window.

3.2 Server window

The main server window contains a list with all instruments defined in the setup (see Fig. 3.2). The standard procedure of *Instrument Server* is to populate this list with the instruments that can be controlled by the computer. Once the instruments are defined and properly configured, they are ready to be used by the *Measurement* program.

3.3 Adding instruments

To add an instrument, click the “Add” button or select “Edit/Add...” from the pull-down menu. The program will scan the global and local *Instrument driver*-folders (defined in the preferences, see Section [PrefsFolder](#)), and bring up a list with available drivers. Select the instrument type to be added and define the communication interface and address. The instrument can also be given a unique name, which is convenient if many instruments of the same type are present in the setup.

3.4 Configuring instruments

Once an instrument has been added to the server, it needs to be configured to perform the desired operation. Select the instrument to be configured in the instrument list and click the “Config” button (or just double-click the instrument name). This will bring up a window with instrument configuration settings (see Fig. 3.3 for an example of a driver for a DC source). The window contains a list with (at least) two sections with controls:




Fig. 3.3 Driver configuration window for a current source. In addition to controls defining the instrument configuration, there are buttons for sending and retrieving the configuration from the hardware. Some quantities, like “Voltage” in the figure, have additional controls for defining sweep rates.

Communication:

This section contains communication controls that define the interface type and address. In addition, if the driver supports multiple instrument models with different installed options, the model type and available options will be shown here.

Settings:

This section (and all other sections, if present) contain instrument-specific configuration settings.

The toolbar at the top of the window provides the following buttons and controls for communicating with the hardware:

Set cfg:

Send the configuration defined in the dialog to the instrument hardware. This requires the communication interface and address to be properly defined.

Get cfg:

Read the configuration from the instrument hardware and update the driver dialog.

At startup:

This controls defines the operation to be performed directly after the instrument driver has started. The default is “*Set config*”, which will configure the instrument hardware according to the settings in the driver dialog. Other options are “*Get config*”, which will read the configuration from the instrument hardware and update the *Labber* driver configuration, or “*Do nothing*”, in which case neither the hardware configuration nor the *Labber* driver configuration are updated.

Start:

When clicking this button, the *Instrument Server* will connect to the instrument and perform the operation defined by the “*At startup*”-control. After successfully performing these tasks, the instrument will be in the *Active* state (marked by an indicator in the lower-right hand corner of the driver window and in the *Instrument server* window). Note that once the driver is active, any subsequent changes made to any of the controls will directly be sent to the instrument hardware. If an instrument is controlled by a client, it is no longer possible to change the configuration from the driver window (all controls will be grayed out). Values can still be sent to or read from the instrument, but only by using the “*Set Value*” or “*Get Value*” buttons in the server window, or if a client asks a value to be measured/updated. The “grayed out”-behavior can also be turned on by default from the “*Server*” section of the *Preferences* dialog (see Section [PrefsServer](#)).

Stop:

This will take the instrument driver out of the *Active* state, stop any eventual instrument operation and close the communication interface.

3.5 Instruments with vector-valued quantities

Some instruments like oscilloscopes, network analyzers and digitizers measure not only scalar values but also traces containing vector values. Drivers for such instruments contain a few extra controls (see Fig. 3.4. for an example). If the “Show trace” checkbox is enabled, the user can acquire and plot the current instrument data by selecting a trace to show and clicking the “Get trace” button. The “Save trace...” button allows the trace currently visible to be saved to the log database. Note that instrument driver must be started and in the active state to acquire and show data traces.




Fig. 3.4. Example of an instrument driver that returns vector-valued data.

3.6 Keeping track of open client connections

Once all instruments are defined, clients can connect to the server to control and measure instruments quantities. The server-client model of *Labber* is very flexible: The *Measurement* program can setup experiments that involve instruments connected to different servers on different computers, and a single server can handle simultaneous calls from multiple measurement programs. This flexibility also brings potential complications, like situations where two clients simultaneously try to access the same instrument. To avoid these complications, the server provides a way for clients to exclusively lock an instrument, thereby preventing other clients from accessing it. The locks are described in more detail in Section [MeasDriverCfg](#).

To keep track of open connections and locked instruments, the *Instrument Server* program features a dialog that lists open client connections and the instruments those clients are using. The dialog is shown by selecting “*Server/Show Open Connections...*” from the menu bar.

3.7 Troubleshooting - Timing statistics

The *Instrument Server* program keeps track of the time each instrument driver needs to perform operations, which can be useful information when benchmarking instrument communication. To turn on the timing statistics, select “*Tools/Show Timing Statistics*” from the *Instrument Server* menu bar. This will add two columns to the main *Instrument Server* window, one displaying the total number of calls performed to a specific instrument quantity, and one displaying average the time per call. The timing statistics can be reset for all instrument by selecting “*Tools/Reset Statistics*”, or for individual quantities by right-clicking the item and selecting “*Reset Timing Statistics*”.

3.8 Troubleshooting - Instrument and Network logs

The *Instrument Server* program keeps logs of recent activities, both for instrument and network communication. The log files are useful if problems arise with instrument communication or if clients have difficulties connecting to the server. To inspect the log files, select “*Log/View Instrument Log*” or “*Log/View Network Log*” from the *Instrument Server* menu bar. The log files provide dated entries with the data strings sent to or received from instruments or from clients.

The amount of logging detail can be controlled in the preferences dialog (see Section [PrefsServer](#)); select “*Debug*” for the most detailed information. However, once the problems have been resolved and the instruments and networks are working as expected, it is recommended to reduce the logging detail level to minimize overhead.

4 Controller

In addition to standard instruments, *Labber* also provides special controller instruments for implementing functionality such as PID controller loops. The controller instruments work by reading an input value from a separate instrument such as a thermometer, applying a controller logic to regulate temperature (for example), and then sending the controller output value to another instruments such as a heater.




Fig. 4.1. The user interface for the PID Controller driver.

4.1 Controller operation

To use the controller functionality, start by adding a controller instrument to the *Instrument* server by clicking the *Add Instrument*-button in the *Instrument* server toolbar. *Labber* provides a built-in PID controller, and additional custom controllers can be created as described in Section [ControllerDriver](#). In addition to the usual sections and settings specific to the particular driver, a *controller* driver also have a number of extra settings related to running the controller loop. An example of the built-in *PID Controller* driver dialog is shown in Fig. 4.1, with the controller settings seen in the right-hand side of the figure. The dialog contains the following settings:

Enabled:

If checked, the controller loop will run in the background and call the input/output instruments at a fixed interval set by the **Period**-setting.

Period:

Intended controller loop period, in seconds.

Measured period:

Actual controller period, which may be different than the set period depending on the time it takes to read/write the input/output values from/to the instruments.

Input signal:

Input signal for the optimizer.

Input value:

Current input value.

Output signal:

Output signal for the optimizer.

Output value:

Current output value.

To set up the controller, first start the instrument driver by clicking the *Start*-button. Next, select the proper *Input* and *Output signals* from the pull-down controls. Finally, set the intended controller period, make sure that both the *Input* and *Output* instruments are running, and then press the *Enabled* checkbox to start the controller. The controller loop will now run in the background and call the input/output instruments at a fixed interval set by the **Period**-setting.

4.2 Improving controller performance

If the controller needs to run at a high repetition rate, set the **Period** control to **0.0** to run the controller loop continuously without gaps. The actual controller loop period will not be zero due to the time it takes to read/write values from the instruments.

Note that the updating the user interface introduces a slightly delay, so for the fastest operation it is advised to run the controller with its dialog window closed. The measured controller loop period can be probed even if the controller window is closed by expanding the `PID Controller/Controller settings` items in the main *Instrument* server dialog.

5 Scheduler

The *Labber Instrument server* contains a scheduler that allows the user to define a queue of experiments to run, as well as functionality for repeating a specific measurement at fixed intervals, for example once per day. The scheduler automatically launches and executes the measurement program whenever an experiment is due.




Fig. 5.1. The user interface for scheduling a measurement.

5.1 Scheduling measurements

Measurements can be scheduled from the user interface, or from the Python API. To schedule an experiment from the user interface, open the *Instrument server* program and select “Scheduler/Schedule Measurement” from the main pull-down menu. This will open a dialog (see Fig. 5.1) with the following settings:

Name:

An optional name for the measurement

Path:

Path to *Labber* measurement configuration to run, in `.labber`, `.json` or `.hdf5` format.

Priority:

Checkbox for setting priority in scheduling system. If a prioritized and non-prioritized measurement are both ready for execution at a specific time, the prioritized one will run first.

Schedule time:

Scheduled time for measurement to run. If the date is in the past, the measurement will execute as soon as the dialog is closed.

Repeat periodically:

If checked, the experiment will be repeated at a fixed interval. If unchecked, the measurement will run only once.

Repeat period:

Repeat interval, in hours.

When closing the dialog, the measurement configuration will be added to the queue of experiments to execute. If there are no other experiments in the queue, and if the “*Schedule time*” is right now or in the past, the measurement will start as soon as the dialog is closed.

To view a list of scheduled measurements from the user interface, select “*Scheduler/Schedule Measurement*” from the *Instrument server* pull-down menu. In addition to displaying the measurements currently scheduled in the queue, the dialog has an option to remove a scheduled measurements (Fig. 5.2).




Fig. 5.2. The user interface for displaying the list of scheduled measurements.

5.2 Scheduler settings

By default, the scheduled experiments will run in a separate process from the one used by the main *Measurement* program. This allows a queued experiment to execute at the same time as one launched from the *Measurement* program. However, this may cause issues if both experiments are trying to access the same resource, for example a specific instrument in the *Instrument server*.

This can be avoided by unchecking the setting “*Run queued experiments in separate process*” under the section “*Measurement/Advanced*” in the *Labber* preferences. If unchecked, an experiment started from the *Measurement* user interface will not start immediately upon pressing *Start* in the dialog, but rather be added to the scheduler queue and execute when the other experiments in the queue have finished.

Note that a restart of both the *Instrument server* and the *Measurement* program may be required for the changes to fully go into effect.

6 Measurement program

The *Measurement* program allows instrument quantities to be measured as a function of other parameters. The program is highly flexible, allowing multi-dimensional sweeps involving any instrument quantity defined in the *Instrument Server*. The *Measurement* is started by from the system tray menu (“*Show Measurement Editor*”) or by selecting “*Window/Show Measurement Editor*” from the main *Instrument Server* window.




Fig. 6.1. The main Measurement configuration window.

6.1 Measurement configuration

The main measurement configuration window is shown in Fig. 6.1. The left-hand panel contains a list with instrument quantities (or *Channels*) involved in the experiment, the top-right section defines the sequence of *Channels* to sweep, while the lower-right panel shows a list of channels to measure. Measurements are easily configured by dragging *Channels* between the lists.

6.2 Adding channels

The first step for setting up a measurement is to define the *Channels* involved in the experiment. A *Channel* represents an instrument quantity on an *Instrument Server*, together with additional properties like name, unit, conversion factors and limits. The easiest way to define channels is to add instruments already present on an *Instrument Server*, which is done by clicking “Add Instruments from Server” in the main *Measurement* configuration window. This will bring up a dialog with options for connecting to an *Instrument Server*, and a list of instruments that can be added to the measurement.

There is also an option for adding channels without having the corresponding instrument previously defined on a server. Choosing the “Edit/Add Instruments...” from the menu bar will bring up a dialog where the user can select which instrument to use and how to communicate with it. In this case, the user needs to specify both the communication protocol of the instrument as well as the server address, so that the new instrument can be created on the *Instrument Server* when starting the measurement.

By default, the program will add *channels* for every quantity active in the instrument configuration. To minimize clutter and allowing an easy overview of the measurement setup, it's advisable to remove channels that will not be controlled from an experiment. This is done by selecting a quantity and pressing the “Remove”-button below the list. Quantities can always be re-introduced later by clicking the “Add”-button. In addition, the value of any instrument quantity can be controlled by opening the *Instrument driver* configuration window (either by double-clicking the instrument name in the *Channels* list or by selecting an instrument and pressing “Show cfg...”).

Note that the *Instrument driver* configuration window serves different purposes in the *Measurement* program and in the *Instrument Server*. In the *Instrument Server*, the instrument configuration dialog is used to *directly* control the hardware settings, meaning that any changes to the dialog will directly affect the state of the hardware. In contrast, in the *Measurement* program the dialog is used to define a configuration that will be used in a specific *Measurement*, but no changes are made to the hardware until the measurement is started. To avoid confusion, *Instrument driver* configuration windows have a different background color when opened within the *Measurement* program and in the *Instrument Server*.




Fig. 6.2. The Channel configuration window allows the user to define properties of the physical quantity measured by an instrument.

Channels also contain properties for describing the physical quantity measured by an instrument. The properties are set in the *Channel configuration* window (see Fig. 6.2), which is brought up by selecting a channel and clicking the “Edit...”-button. The dialog allows the user to set channel max/min limits, and to define conversion factors between the physical quantity investigated in the experiment and the quantity measured by the instrument. An example where such a conversion is useful is when current biasing a circuit by applying a voltage over a large resistor in series with the circuit. In this case, the physical quantity would be current (with units Ampere), while the instrument quantity would be voltage. The equations for converting between physical and instrument units are defined in the text box next to the *Conversion factors*-controls.

6.2.1 Instrument configuration - locks

After the channels have been added to the *Measurement* configuration, the corresponding *Instrument driver* configuration window can be shown by double-clicking the instrument name or selecting the instrument and clicking “Show Config...”. In addition to the settings listed when describing the *Instrument Server* (see Section [ConfigInstr](#)), the dialog contains an extra checkbox (“Lock instrument on server”, under the section “Communication”) for determining whether the instrument will be used exclusively by the current experiment. If the control is checked (default behavior), no other clients can connect to or change the instrument values during the duration of the measurement. See Section [OpenClients](#) for more information about locks.

6.3 Sending and retrieving values from instruments

The “Set Value...”- and “Get Value”-buttons below the channel list allow the user to quickly set or retrieve the current instrument value. Note that the “Set value...”-operation will immediately send the new value to the instrument hardware.

Instrument values can also be controlled from the *Instrument driver* configuration window, which is opened by double-clicking the instrument name or selecting the instrument and clicking “Show Config...”. However, in contrast to the “Set/Get Value”-buttons, changing the value of a control in the *Instrument driver* window will only update the local value kept in the *Labber* configuration. The actual instrument hardware is not updated until the user clicks “Set Cfg” in the driver window, or when the measurement is started (provided that the “At Measurement Start”-option in the driver window is set to “Set config”).

6.4 Defining step sequences

A measurement consists of a list of *Step sequences* that output values to instruments in a specified order. To define a *Step sequence*, drag the channel to sweep from the *Channels* list on the left to the *Step sequence* list on the top right of the main *Measurement* configuration window. This will bring up the *Basic settings*-dialog for defining the range of values to output. Once defined, the step items can be re-ordered by dragging the entries within the list.

6.4.1 Step setup - Basic settings

The basic settings dialog allows the user to define single-point step values or basic ranges, either by defining start-stop or center-span values. The step size is specified either by setting a fixed step size, or by defining the total number of points in the step range, see Fig. 6.3. When defining the number of points in the range, the user can set the interpolation to be linear or logarithmic. Note that logarithmic interpolation only works if all values in the step range are positive.




Fig. 6.3. Basic dialog for defining step sequences.

If the channel is sweepable or if there are scaling factors defined between physical/instrument units, the dialog contains a few extra controls as described in Section [AdvancedStepSetup](#).

6.4.2 Step setup - Advanced settings

The advanced settings dialog contains a few extra controls to provide better control of the step parameters. First, there is a list with step ranges, allowing multiple ranges to be defined with different step sizes (see Fig. 6.3 for an example). Use the “Add...”, “Edit...” and “Remove”-buttons to add/edit ranges, and drag the entries in the list to make the values appear in the right order. The graph in the upper-right corner of the dialog shows a visual representation of the step output values, with the step number on the y-axis and the corresponding output value on the x-axis.




Fig. 6.4. Advanced dialog for defining step sequences.

In addition, the upper-left part of the dialog contain the following controls for fine-tuning the step sequence:

Step units:

The control sets whether the step values are given in instrument or physical units of the channel (see Section [Channels](#)). The step values are updated to reflect the unit settings whenever the control is updated. The control is only visible if physical/instrument unit conversion factors have been defined in the *Channel* setup dialog.

Wait after each step:

Time to wait after a step value has changed. Note that the actual time to wait will be the maximum of this time and the delay time between step and measure as set in the main *Measurement* configuration window (see Section [Timing](#)).

Alternate step direction:

If checked, the execution of the step sequence in multi-dimensional experiments will alternate between forward and reversed direction, eliminating the need to go back to the first step point between loops. This feature is also useful when looking for hysteresis when sweeping a field up and down.

After last step:

Defines the operation to perform after the step sequence has completed. Possible values are *"Goto first point"*, *"Stay at final"* or *"Goto specified value"*. Default behavior is *"Goto first point"*.

6.4.3 Step setup - Sweep mode

If the step channel is sweepable, the sweep mode controls provide a few extra user interface elements for controlling the sweep settings. For more information on how to write drivers that supports sweeping, see Section [SweepDriver](#).

Sweep mode:

The program supports three different sweep modes:

Sweep mode - Off:

Sweep mode is off, step values are set directly.

Sweep mode - Between points:

The instrument is swept between step points, but the output is held constant while acquiring data for the log channels.

Sweep mode - Continuous:

In this mode, the instrument is configured to continuously sweep from the first to the last value in the step list. The log channels are being measured at the points defined in the step list, but the program will not stop sweeping the output channel while acquiring data for the log channels.

Rate:

Sweep rate, in units s^{-1} . The sweep rate is also shown in units min^{-1} . Note that the sweep rate in this dialog will overwrite the sweep rate defined in the configuration window of the corresponding instrument driver.

Time between points:

If sweep mode is *Continuous*, this text shows the typical interval between measurement points for the given sweep rate and step list.

Use different sweep rate outside loops:

If checked, a numerical control will appear below the checkbox, allowing the user to define a separate sweep rate for sweeping to the init/final values and for sweeping between loops. If unchecked, the program will use the common sweep rate defined in the control above when setting init/final/between loop values.

6.4.4 Step setup - Channel relations

One useful feature of the *Advanced step configuration* is the ability to define relations between channels. For instance, imagine a situation where we want to sweep two voltages V1 and V2 in a way that V2 is always exactly 1.5 V higher than V1. To implement this, we first define the step sequence for V1 as usual. Next, we create a step configuration for V2 and switch to the advanced settings. Clicking the “*Enable channel relations*” will allow us to enter an equation relating the output of V2 to other channels. The values of other channels are accessible through parameters, shown in the list on the right-hand side of the dialog. The “*Add...*”, “*Edit...*” and “*Remove*”-buttons below the parameter list are used to edit the parameter names. The parameter “x” refers to the step values as defined in the step list in the upper part of the dialog.

For this particular example, we would enter $p1 + 1.5$ in the equation box for channel V2 (assuming that the parameter $p1$ is linked to channel V1). The equation string can involve basic mathematical functions like $\cos(x)$, $\sin(x)$, $\sqrt{x}$, $\exp(x)$, etc... Also, note that the raised operator ($\wedge$) is given by two multiplication signs ($**$).

Before starting a measurement, it is good practice to check that the relation equation produces the intended output. By default, the graph in the top-right corner of the dialog

shows the step values generated by the step list, but it can also be configured to visualize the output of the relations as a function of any other channel in the measurement. The graph contents are set by the “Plot to show”-control to the left of the figure.

6.5 Log channels

The *Log channels* list defines the channels to measure at each step point. To add log channels, simply drag a channel entry from the main *Channels* list on the left to the *Log channels* list. Note that it is not possible to add the same channel to both the step and the log list.

If the checkbox “Log in parallel” is checked, the program will try to measure all channels simultaneously at each step point. If the “Log in parallel” is unchecked, the channels will instead be measured sequentially, starting with the top-most one in the list. The log channels can be reordered by dragging the items within the list.

6.5.1 Log channels limits

Each log channel has an associated range limit, which is defined by double-clicking the log channel or clicking the “Edit...” button below the log channel list. If the measured value falls outside the defined limits during a measurement, one of the following actions will be taken:

Nothing

No action is taken, the measurement continues as usual.

Continue to next step item

The measurement program stops execution of the innermost step sequence, and continues to the next item of the second step sequence.

Stop, stay at current values

Stop the measurement, hold all instruments at the current values.

Stop, go to init/final configuration

Stop the measurement, go to final values as defined in the *Advanced step setup dialog* (see Section [AdvancedStepSetup](#)). The default is to go to the initial values.




Fig. 6.5. Limit options for log channels.

6.6 Timing

The *Timing* section in the lower-right corner of the main *Measurement* configuration window allows the user to set a time to wait between outputting new values to the step channels and measuring the log channels. In addition, the section gives an estimate for how long time it'll take to run the measurement. The estimate is based on the duration of the previous experiment; the expected time needed per point can be adjusted manually, if required.

6.7 Log name, Project and User tags, Comments

When running a measurement, both the measurement configuration and the obtained data are saved into a single file in the database folder (see Section [LogDatabase](#) for a discussion of the folder hierarchy and the structure of the log database). The log file name is defined by the *Log name* control in the toolbar in the upper-right part of the main *Measurement* configuration window. When starting a measurement with a file name that already exists, a dialog will pop up presenting the user with the following options:

Create New

The new measurement will be save into a new log file, with a modified file name ("_2", "_3", "_4", etc, will be appended to the log name).

Append Data

The new measurement will append new data to the old log. This option is only available if the measurement is one-dimensional (that is, if only one of the step

sequences contains more than a single value), and if the previously existing log has the same structure as the new one.

Overwrite

The old log is deleted before starting the new measurement.

The comment field allows experiment-specific descriptions to be added to the measurement configuration file. Note that there is no need to type any information related to instrument settings here; all the instrument configurations are automatically saved into the configuration file.

6.8 Tags

The *Project*, *User* and *Tags* controls provide ways of keeping the log database organized. The controls are shown by clicking the “*Show Tags*”-button in the dialog toolbar.

The *Project* field supports a hierarchy structure, with subprojects separated by a forward slash (“/”). For example, entering “*Sample2/DeviceA/IV-curves*” will put the log in the subproject “*IV-curves*” of subproject “*DeviceA*”, which is located in the project “*Sample2*”.

The *Project* tag can be entered directly into the text field, or by clicking the folder icon next to the control to bring up a hierarchy tree with all projects defined in the database.

The *User* name can be entered directly into the text field, or by clicking the user icon to bring up a list with users already present in the log database. A log file can only belong to a single *Project* and *User*.

Contrary to the *Project* and *User* fields, a log can contain multiple *Tags*. The *Tags* are added/removed by clicking the plus/minus signs next to the tag list. Similar to the *Project* field, the *Tags* support a hierarchy tree of tags and subtags.

6.9 Starting a measurement

Once the step sequences, log channels and the log name have been defined, the measurement is ready for execution. When clicking the “*Start measurement*”-button in the upper-right corner of the *Measurement* dialog, the program will perform the following sequence:

1. The program will go through all step sequences to make sure that all of the step values are valid and within the min/max ranges allowed by the corresponding channels.
2. Next, connections will be established to the *Instruments Servers* of all instruments in the *Channels* list.
3. For every instrument, the program will either set or read the current the hardware configuration, depending on the value of the “*At Measurement Start*”-control of each instrument driver (see Section *AtStartup*). The recommended setting is “*Set cfg*”, since this will assure that instrument is hardware in the same state every time the measurement is performed. Note that the “*At Measurement Start*”-operation will be performed for all instruments defined in the measurement window, even for channels that aren’t used in step sequences or as log channels.
4. The program will perform the measurement by stepping through all the values defined in the *Step sequences*. The order of the step sequences defines the order in which the values are outputted, starting with the values in the top-most step sequence. For example, consider a situation with two defined step sequences: one for “*Channel 1*” with values {1,2,3} and one for “*Channel 2*” with values {10,20}. There are a total of $3 \times 2 = 6$ step points. If “*Channel 1*” occurs before “*Channel 2*” in the sequence list, the program will set the values in the following order:

Step No.	Channel 1	Channel 2
1	1	10
2	2	10
3	3	10
4	1	20
5	2	20
6	3	20

5. On the other hand, if “*Channel 2*” occurs before “*Channel 1*”, the step order will be:

Step No.	Channel 1	Channel 2
1	1	10
2	1	20
3	2	10
4	2	20
5	3	10
6	3	20

6. The values of all log channels are measured at each step point.
7. When the measurement is finished, the program will close the connections to all instruments and all *Instrument Servers* and return to the main *Measurement* configuration window. The new log will be available for viewing in the *Log Browser* window (see Chapter [BrowserDlg](#))

Figure 6.6 depicts the dialog shown when a measurement is running. The list on the left contains a list of the step and log channels defined in the measurement, together with current values and progress indicators (for step channels). The green light indicate that a value is currently being sent/received from an instrument. The graph on the right visualize the measurement progress for the step channel selected in the channel list on the left. Alternatively, if a log channel is selected, the graph will show the currently measured trace for that channel.




Fig. 6.6. The Measurement window.

The tool bar at the top of the dialog contains buttons for showing the measured data in real-time, either as a line plot or as an image map. In addition, there are buttons for skipping traces, pausing and stopping the experiment. The *Skip* button will stop execution of the innermost step sequence, save the current trace, and then continue to the next step item.

6.10 Signal connections

Many experiments involve sending or reading waveforms from instruments like arbitrary waveform generators, digitizers or digital oscilloscopes. For example, imagine an experiment where we want to control the amplitude of a sine signal outputted using an arbitrary waveform generator. The process can be divided into two tasks: The first task is to numerically calculate a waveform with the correct amplitude, the second task is to send that waveform to the output of the arbitrary waveform generator. Another example would be to measure a waveform with an oscilloscope, and then apply some function to extract the signal's amplitude or frequency, which would allow us to record only a single or a few values characterizing the signal instead of saving the whole waveform to the log file.




Fig. 6.7. A Measurement configuration with signal connections for pulse generation and signal demodulation.

The *Measurement* program provides *Signal Connections* for managing situations like this. The idea is to separate the signal generation or signal analyzing from the instrument communication, to make it possible to develop generic *Signal Generator* or *Signal Analyzer* drivers for creating or analyzing waveforms, which operate independently of the specific hardware that is used to output or measure the waveforms. In this way, the waveform generation/analyzing functions can be used interchangeably with instruments from different vendors. For more information on how to create your own *Signal Generator* or *Signal Analyzer* drivers, see Section [DriverINI](#).

The *Signal Connections* button in the toolbar of the main *Measurement* configuration window is used to show/hide a list of signal connections defined in the current setup. The button is only enabled when the setup contains instruments that allow waveform generation/analyzing. The *Signal Connections* control contains a list with all the instrument quantities that can output or analyze a waveform. To make a connection, double-click one of the outputs and select the source signal from the dialog that pops

up, or simply drag a channel that represents a signal source from the main *Channels* list onto the correct output in the *Signal Connections* list. Figure 6.7 shows an example of a *Measurement* configuration with a few signal connections.

Note that signal connections are also possible for scalar-valued channels. To make a signal connection between two scalar-valued channels, click the “*Show scalar-valued signals*”-checkbox below the signal connection list. Scalar-valued channels are listed in italics in the signal connection list, to distinguish them from the waveform signal connections.

As mentioned earlier in this section, two types of signal connections can be made: The first type is when a *Signal Generator* driver is used to generate waveforms that will be sent to the output of an arbitrary waveform generator, for example. The second type of connection is when a waveform that is acquired using an instrument such as a digitizer or an oscilloscope is sent to a *Signal Analyzer* driver. The example in Fig. 6.7 illustrates both examples: The signals generated by the “Pulse Generator” *Signal Generator* driver are configured to be sent to various output channels of a Tektronix Arbitrary Waveform Generator (labelled “AWG” in the figure), whereas the waveforms acquired by the “Acqiris U1084A Digitizer” will be sent to the of the “Signal demodulation” *Signal Analyzer* driver that will extract the amplitude of the waveform at a specific frequency.

When running a *Measurement* that contains *Signal connections*, for each step in the step sequence the program will perform the following sequence:

1. Update step values: Update the values of channels defined in the *Step sequence*.
2. Generate and output signals: If present, calculate signals with *Signal Generator* drivers and send the resulting waveforms to the corresponding instruments outputs.
3. Wait: Wait for the time specified in the *Timing* section in the lower-right corner of the main *Measurement* configuration window.
4. Acquire and analyze signals: If present, measure instrument channels that acquire waveforms, and send the acquired waveforms to the corresponding *Signal Analyzer* drivers.
5. Log results: Measure the channels specified in the *Log channels* list and save them to disk. If a *Signal Analyzer* driver is in use, the *Log channels* list is where the user defines what quantities to store in the log file.

In the example of Fig. 6.7, the step sequence will update a few parameters of the “Pulse Generator” *Signal Generator* driver. Once all parameters have been updated, the “Pulse Generator” driver will calculate new waveforms that will be sent to the “AWG” output channels. After that, the program will wait for 0.1 second to give the sample time to settle, before acquiring two waveforms (“Ch1 - Data” and “Ch2 - Data”) with the “Acqiris U1840A Digitizer”. The measured waveforms will be sent to the “Signal demodulation” *Signal Analyzer* driver, which will analyze the waveforms and return the result in the channel named “Value”, which will be stored in the log file.

When running an experiment that contains *Signal connections*, it is possible to look at the measured waveforms in real-time as they are being acquired and processed. In the window that is visible when an experiment is running, mark the channel to investigate in the step sequence list in the left-hand part of the dialog (see Fig. 6.8).




Fig. 6.8. The measurement window during an experiment with signal connections, showing a measured waveform.

6.11 Hardware timing and synchronization

In standard operation mode, *Labber* handles instrument synchronization by waiting for an instrument to report that all step channels values have been outputted before reading

the log channels. Since this requires communicating with the instruments over the computer, the time synchronization may not be precise enough for certain applications. For such applications, *Labber* supports operating in arm/trig mode and hardware looping for enhancing the synchronization and timing precision.

6.11.1 Arm/trig mode

In Arm/trig mode, log instrument will be armed to wait for an external trigger before starting to acquire data. To turn on arm/trig mode, click the *Arm/trig mode* checkbox below the Step sequence configuration list in the *Measurement* setup window. The trigger channel must be an instrument channel represented by boolean or a button, and it is defined by the pull-down menu next to the check box. The instrument used to generate the trigger must also be represented in the step configuration list.

In Arm/trig mode, the following operations are performed at each point of the measurement sequence:

1. Set output values of the step channels, but instruct the instrument to wait for a trigger before outputting any signals.
2. Arm the log channels to get ready to acquire data.
3. Wait for the time specified in the *Timing* section in the lower-right corner of the main *Measurement* configuration window. This time can be zero.
4. Generate the trigger signal. The output of the trigger should be physically connected to the step and log instruments, so that the step instrument can start outputting signals and the log instruments can start acquiring data.
5. Read out the acquired data.

6.11.2 Hardware looping

Some instruments can perform looping of values within the instrument hardware. This allows for implementing more efficient looping, since there will be no need for the computer to send new values to the instrument at each step value. This mode requires that the instrument outputting values support hardware looping, that the instrument reading values supports both hardware looping and hardware arming, and that there is a trigger defined for instrument synchronization. If the instruments used in the *Measurement* configuration fulfill these requirements, hardware looping can be

activated by clicking the *Hardware loop* checkbox next to the arm/trig mode trigger controls.

In hardware loop mode, the top-most step item in the *Step sequence* configuration list of the *Measurement* setup window will be controlled by the instrument hardware. Instead of setting and getting values point-by-point using the computer, the looping of the top-most step item will be handled in the following way:

1. Calculate the number of values *n* to step in the top-most step item.
2. Send all *n* values of the top-most step item to the output instrument, but instruct the instrument to wait for a trigger before outputting any signals.
3. Arm the log channels to get ready to acquire data. The log instruments will be configured to acquire *n* values.
4. Wait for the time specified in the *Timing* section in the lower-right corner of the main *Measurement* configuration window. This time can be zero.
5. Generate the trigger. The output of the trigger signal should be physically connected to the step and log instruments, so that the step instrument can start outputting signals and the log instruments can start acquiring data. The step instrument will output the *n* values in the prescribed order, and the log instrument will acquire *n* values.
6. Read out the acquired *n* values.

Since multiple values will be read out at once in hardware looping mode, the progress bar and the time estimate shown in the *Measurement* window during experiments may not update often enough to be accurate.

6.12 File locks

The software uses a file locking system to prevent the *Log Browser*, *Log Viewer* or *Measurement Editor* from making changes to a file while a measurement is running. The locking mechanism works by creating an empty file with the same name as the log file that is being measured, but with the ending *.lock*. Under normal operation, the *.lock*-file will be removed when the measurement is completed, but if the experiment was unexpectedly interrupted (for example, if the computer suddenly lost power), the *.lock*-file will persist and prevent any changes from being made to the file. To

manually remove the locking, go to the folder where the log file is stored, locate the *.lock*-file and delete it.

6.13 Measurement settings

In addition to the general preferences described in Section [Prefs](#), there are a number of settings that can be uniquely defined for each specific *Measurement*. These specific *Measurement settings* can be accessed by clicking the “*Show Settings*” toolbar icon in the top-left corner of the *Measurement* dialog. The settings are described in detail below.

6.13.1 General

Send values in parallel:

When outputting multiple values in a measurement step sequence, define if data should be sent to all instruments in parallel, or sequentially one after each other. Default value is *True*.

Only send signal if source instrument has been updated:

If checked, Labber will only perform a signal connection if the source instrument has been updated since last call.

Data compression:

The value ranges from 0 (no compression) to 9 (maximum compression). Higher compression reduces the log file size, but may slightly increase time for loading/saving data.

6.13.2 Optimizer

The optimizer functionality and the corresponding settings are described in more detail in Section [Optimizer](#) below.

6.14 Comparing Measurement configurations

For complex measurement scenarios containing a large number of instruments, it is sometimes difficult to keep track of all setting and parameters involved in the experiment. For these cases, *Labber* provides a convenient feature to compare and highlight differences between the current scenario and a previous measurement saved in

the Log database. The function is accessed by selecting *Tools/Compare Configurations* in the pull-down menu and selecting the measurement configuration to compare the current scenario to.

7 Optimizer

In *Optimizer* mode, the *Measurement* program will try to minimize the value of an expression based on the measured channels instead of looping through the step channels in a pre-determined sequence. This can be useful when the goal of the experiment is to minimize a certain quantity as opposed to mapping out the value of the quantity over the full parameter space.

7.1 Optimizer operation

To enable the optimizer, simply double-click on one of the *Step items* in the *Step sequence* list in the *Measurement* program, switch to “*Basic settings*” (if not already in that mode), then click the “*Optimize...*”-button to convert the step item to an optimizer parameter. Instead of sweeping over the parameter, *Labber* will try to optimize the cost function (see below) by varying the parameter over the range specified by the “*Min value*” and “*Max value*” controls in the dialog. The various options in the dialog are described in more detail in Section [ParameterSettings](#) below.

7.1.1 Cost function

The next step is to define the cost function and the general settings of the optimizer. These options are available by clicking the “*Show Settings*” toolbar icon in the top-left corner of the *Measurement* dialog, and clicking *Optimizer* in the section list in the left part of the dialog. The most important setting is the optimizer cost function, which is defined by the expression in the “*Minimization function*”-control. The cost function takes the latest measured values of the log channels as inputs and must return a single scalar value. The optimizer algorithm will then try to minimize the value of the cost function by iteratively varying the various optimizer parameters.

The inputs available to the cost function are the latest values of the measured log channels, provided in the numpy list `y`. Each element in the list corresponds to a channel, and the order of the elements is the same as the order at which the log channels appear in the *Measurement Editor*. If you are using a single log channel, its value can be accessed by `y[0]`. However, note that `y[0]` may be scalar or vector-valued, depending on if the particular log channel returns a trace of a single value. For the optimizer to work, the cost function must always return a scalar, so if your log channel is

vector-valued you need to apply some operation to convert the vector to a scalar. For example, `mean(y[0])` would optimize with respect to the mean of the measured trace. In addition to `y`, the vector `x` with the latest values of the optimizer parameters is also available as an input to the minimization function. You can use any Python and numpy expression when defining the cost function.

7.1.2 Termination and convergence criteria

There are three possible criteria for defining when the optimizer should terminate the optimization process.

Absolute target reached:

If the value of the cost function is less than the *Target value*, the optimizer will terminate.

Relative tolerance reached:

If the change in the cost function between calls is smaller than the value given by “*Relative tolerance*”-setting, AND if the change in the optimizer parameter values between calls are smaller than the “*Precision*” setting of each parameter, the optimizer will terminate. Note that both criteria need to be fulfilled for termination.

Max number of evaluations reached:

The optimizer will automatically terminate after performing the number of measurements specified by “*Max evaluations*”.

By default, the “*Target value*” is set to minus infinity, which means that it will never terminate the optimizer. In addition, the “*Relative tolerance*” is set to infinity by default, which means that only the “*Precision*” of the individual optimizer parameters matter for relative convergence.

Note that the termination/convergence criteria may differ for different optimizer algorithms, the description above only refers to the default *Nelder-Mead* optimizer provided by *Labber*.

7.1.3 Running an optimizer measurement

When running a measurement with the optimizer enabled, *Labber* automatically will add a step item named “*Optimizer iteration*” that handles the optimizer loop. Note that it is possible to run an experiment with a mix of optimized and non-optimized parameters, where the optimizer will execute to find the optimal value of one parameter while stepping over different values of another parameter.

7.2 Optimizer settings

In order to use the optimizer, both the general optimization protocol and the individual optimization parameters must be configured. The various settings are described below.

7.2.1 General optimizer settings

These settings define the cost function and the algorithm-specific settings of the optimizer, and can be accessed by clicking the “*Show Settings*” toolbar icon in the top-left corner of the *Measurement* dialog. The settings are described in detail below.

Method:

Algorithm used for optimization.

Max evaluations:

Maximum number of function evaluations/measurements performed before terminating the optimization.

Minimization function:

Function for optimizer to minimize. The measured channels are available in the variable `y`, which is a list of log channel values. Each list item may be a number or a numpy array, depending on the channel datatype. Default is `min(y[0])`, which will minimize the value of the first log channel.

Target value:

Absolute value of minimization function at which the optimization will terminate. Default value is `-inf`, which will prevent the optimizer to terminate until the other optimization goals are met.

Relative tolerance:

Change in minimization function between iterations that is acceptable for convergence. Default value is `inf`, which will make the optimizer run until the `Precision`-value of all involved parameters are met.

7.2.2 Individual parameter settings

These settings are individual to each optimization parameter and can be accessed by double-clicking a channel in the *Step sequence* list and going to “Optimize...”-mode.

Start value:

Initial value for parameter.

Initial step size:

Initial step size for the parameter.

Min value:

Lowest parameter value allowed during the optimization procedure.

Max value:

Highest parameter value allowed during the optimization procedure.

Precision:

Target precision for optimizer that will trigger optimizer termination.

7.3 Custom optimizers

It is possible to create custom optimizer modules to implement a specific optimization protocol. The sections below describe how to define and test a custom optimizer algorithm.

7.3.1 Defining custom optimizers

It is recommended to use one of the already present optimizer configuration files as a template. The custom optimizers should be contained in a single python `.py` file, which must contain a function called `optimize` that takes exactly two parameters:

config:

Python dict with optimizer settings. The keys have the same names as the labels of the optimizer settings in the *Measurement* program. The individual parameter settings are stored as a list in the same dictionary, with key `optimizer_channels`.

minimize_function:

Python callable that takes exactly one argument (`x`). The function will run the Labber measurement for the provided parameter values `x`, where each value in the vector `x` corresponds to an optimizer parameter. The function is typically passed directly to the `scipy` optimizer, see the provided optimizer `Nelder-Mead` for an example.

The function must return a Python dictionary with results from optimizer, using `scipy`'s `OptimizeResult` format. The only necessary key is "`x`", containing the final optimizer parameters.

When creating a new optimizer, the python file should be given a unique name and placed in the *local* optimizer folder (the folder named "*Local optimizers*" in the *Preferences* window), instead of the global one (within your installation directory). This allows the user's own optimizers to be kept separately from the optimizers provided by *Labber*, and it also prevents optimizers written by the user from being deleted when updating the *Labber* program to a newer version.

Note that even when making additions/changes to an existing optimizers from the global folder, the best practice is to copy that optimizer file from the global folder to the local folder, and only make changes to the optimizer version. If optimizers with the same names exist in both the local and the global optimizer folders, *Labber* will always use the optimizer in the "*Local optimizer*"-folder.

7.3.2 Defining optimizers settings

For custom optimizers, it is possible to define optimizer-specific configuration parameters in addition to the general settings in Section [OptimizerSettings](#) above. The optimizer-specific settings are defined by adding a function `define_optimizer_settings()` to the same python `.py` file that contain the optimizer code. The function should return a list of python dicts, where each dict represents a specific setting. The settings are

defined in a similar way to quantities of an instrument driver (see Section [Quantities](#)), with the difference that the settings are specified in a python function instead of a `.ini` configuration file. Each setting must define the `name` and `datatype` parameter, all other parameters are optional.

The custom settings will show up in the *Optimizer*-section of the *Settings*-pane of the *Labber Measurement* dialog, allowing the user to change their values prior to running a measurement. The values of the custom parameters will then be accessible as entries in the `config` input in the `optimize` function defined above.

As an example, the code below will define custom settings with three parameters for the *Bayesian-Gaussian-Process* optimizer.

```
def define_optimizer_settings():
    """Define extra settings for optimizer

    Returns
    -----
    optimizer_cfg : list of dict
        List of configuration items for optimizer, each item is a dict.
        Necessary keys are "name" and "datatype".

    """
    # Bayesian optimization settings
    optimizer_cfg = [
        dict(name='Acquisition function',
            datatype='COMBO',
            combo_defs=['LCB', 'EI', 'PI', 'gp_hedge'],
            def_value='gp_hedge',
            tooltip=('See https://scikit-optimize.github.io/ for more info'),
        ),
        dict(name='kappa',
            datatype='DOUBLE',
            def_value=1.96,
            state_item='Acquisition function',
            state_values=['LCB', 'gp_hedge'],
            tooltip=('Controls how much of the variance in the predicted ' +
                'values should be taken into account. Higher value ' +
                'favours exploration over exploitation and vice versa'),
        ),
        dict(name='xi',
            datatype='DOUBLE',
            def_value=0.1,
            state_item='Acquisition function',
            state_values=['EI', 'PI', 'gp_hedge'],
            tooltip=('Controls how much improvement one wants over the ' +
                'previous best values. Higher value ' +
                'favours exploration over exploitation and vice versa'),
        ),
    ]
```

```
    ),  
]  
return optimizer_cfg
```

7.3.3 Using custom optimizers

To make the new optimizer available to *Labber*, place it in the local optimizer folder and click the menu alternative “*Tools/Reload Optimizers*” in the *Measurement Setup* dialog. This will scan the optimizer folders and update the “*Method*” control in the general optimizer settings.

It is highly recommended to first test the optimizer in a pure Python environment. For an example of how to test the optimizer, see the code at the end of the file `NeIded-mead.py` provided in the global optimizer folder.

8 Log Browser

8.1 Database

The log database consist of a set of *Labber* log files organized in a special folder structure. When running a measurement, the program will save the configuration and data in a single file in the folder “<Database folder>/xxxx/yy/Data_yyzz”, where <Database folder> is the base database folder as set in the *Preferences* window (see Section [PrefsFolder](#)), xxxx is the current year (four digits), yy is the current month (two digits) and yyzz is the current month+day (total four digits). As an example, the database foLabber_manuallder for logs created on January 29, 2014 would be “<Database folder>/2014/01/Data_0129”.

The user can add additional data like images, scripts or even subfolders to folders within the main database folder. The *Log Browser* and the *Measurement* programs will ignore any additional files when scanning for *Labber* log files.

8.2 Log browser dialog

The *Log browser* is used to browse through the measured data and give a quick overview of the individual log files. The *Log Browser* is started by from the system tray menu (“*Show Log Browser*”) or by selecting “*Window/Show Log Browser*” from the main *Instrument Server* window. When starting, the program will scan through the default database folder for log files.




Fig. 8.1. The Log Browser window.

Figure 8.1 shows the main *Log Browser* window. The dialog consists of a tool bar at the top, with sections on the left showing the database structure, a list in the center with log files, and finally a graph and info controls on the right giving a preview of the selected log file. The function of the various controls are described below.

8.2.1 Database hierarchy view

The four fields *Project*, *Tags*, *User* and *Date* on the left-hand side of the *Log Browser* dialog give an overview of the database hierarchy and allow the user to limit the selection of logs visible in the main log list. The *Project*, *Tags* and *Date* entries can be expanded to reveal subfolder selection. The filtering process is exclusive, meaning that only log files that fulfill the constraints of all the four fields are shown in the log list. Selecting the “- All entries -” item in one of the controls will disable the filtering for that field. The *Project*, *Tags*, *User* and *Date* fields can be hidden from the *View*-menu in the main menubar.

8.2.2 Log list

The list in the center of the dialog shows the logs in the database that fulfill the constraints set by the hierarchy fields on the left. Logs that contain particularly

important data can be starred, which will make them easier to spot in a large selection of logs. To star a log, either right-click the entry and select “Star” from the pop-up menu, or select a log and use the “Star”-button in the tool bar at the top of the window, or just press the space bar after a log has been selected.

The controls in the toolbar above the log list provide options for showing only starred log and for filtering logs by name. The log list can be sorted by *Name*, *Creation date* or *Sweep dimensions* by clicking the corresponding list header.

8.2.3 Graph / Log info

The graph in the top-right corner gives a quick preview of the contents of the selected log file. By default, the graph will show an image map if the data is two-dimensional in nature, and otherwise a line plot containing the first few entries in the log. The preview can be changed by creating a *View* in the *Log Viewer*, see Section [Views](#) for more information on *Views*. Below the graph, there are controls showing the *Step sequence*, the *Log channels* and the comment (if present) of the currently selected log.

8.2.4 Tool bar

The toolbar contains the following buttons:

Open Database:

Open another database than the default one used by the *Measurement* program. The dialog that opens should be pointed to the folder containing the *Year*-folders of the log database.

Reload Database:

Reload the current database and scan through all the log files. This is needed if log files have been manually added to the database folder.

Star:

Star/unstar the currently selected log.

Plot Traces:

This will open the currently selected log in the *Log Viewer*, see Section [LogViewer](#) for more information about the *Log Viewer* dialog.

Plot Image Map:

If the selected log contains 2-dimensional data, this button will open the data as an image plot in the *Log Viewer*.

Show config:

The button will show/hide an additional side bar with all the instrument configurations in the currently selected log. At the bottom of the sidebar, there will also be a checkbox “*Show all quantities*”; if checked, all instrument quantities and values are shown in the list, otherwise only the quantities present as channels in the *Measurement* configuration are displayed.

9 Log Viewer

The *Log Viewer* provides an environment for plotting and analyzing log files. The viewer is started by double-clicking a log in the *Log Browser*, which will bring up the main *Log Viewer* window (see Fig. 9.1). The dialog contains a list with all entries in the log, a plot showing the currently selected entries, and a set of controls on the left for configuring the plot. The data can be plotted as individual *Traces* or as an *Image*, and the user can quickly switch between the two modes using the *Traces/Image* buttons in the tool bar. The *Log Viewer* also provide options for saving a *View* of the log data, which allows plot settings and analysis configurations to be easily restored.




Fig. 9.1. The Log Viewer dialog in Trace mode.

9.1 Plot config

The plot configuration tools on the left-hand side of the dialog determines what is plotted in the graph. The controls are:

Y-channel:

The top-most control sets the channel to plot in the graph. If the channel contains complex values, the user can choose to plot the *Real*, *Imaginary*, *Magnitude* or *Phase* of the signal.

X-channel:

The next control defines the x-axis in the plot.

Rotation/x (complex only):

This introduces a phase rotation per x-unit to a complex trace, which has the same effect as compensating for electrical delay when plotting signal transmission versus frequency. The “*Calculate*”-button next to the control estimates the rotation value that will best compensate such delays. The controls are only visible if the plotted quantity is complex.

Fixed Rotation (complex only):

This introduces a fixed phase rotation to a complex trace. The “*Calculate*”-button next to the control estimates the rotation value that will maximize the signal in the real component, while minimizing the signal in the imaginary one. The controls are only visible if the plotted quantity is complex.

Plot in dB (complex only):

If plotting magnitude, this converts the value to dB, using the formula $20 \cdot \log_{10}(y)$. The control is only visible if the plotted quantity is complex.

Unwrap angle (complex only):

If plotting the phase, this will unwrap phase jumps around $\pm 180^\circ$. The control is only visible if the plotted quantity is complex.

Auto-rotate (complex only):

If checked, the program will automatically rotate the phase of each trace to maximize the real part of the signal. This is the same as the “*Fixed Rotation/Calculate*”-button, but the algorithm is applied for each trace individually,

which means that the individual traces will generally be rotated by different amounts.

Operation:

Applies an operation to each selected trace. The following operations are available:

$y - \langle y \rangle$:	Subtract the average value from each trace.
Normalize:	The traces are normalized using the formula $(y - \langle y \rangle) / \text{std}(y)$.
d/dx:	Calculate the numerical derivative.
FFT:	Numerical Fourier transform. Only positive frequency are shown in the resulting plot.
Histogram:	Bin the data into a histogram, the number of bins is set by the "Bins" control.
Histogram-2D:	Bin complex data into a 2D histogram, with the x/y-axes given by the real and imaginary parts of the data. The function only works for complex values.

Smooth:

Smooth the trace by taking a running average over the number of data points specified in the control.

Traces:

Applies an operation to the collection of all selected traces. The following trace operations are available:

Show individual:	Standard operation, plot the individual traces as they are.
Subtract first:	Subtract the first selected trace from the following traces.

Subtract previous: Subtract the previous trace. The plot will contain $N/2$ elements, with data (Trace 2 - Trace 1), (Trace 4 - Trace 3), (Trace 6 - Trace 5), etc...

Average: Plots the average of all selected traces.

Standard deviation: Plots the standard deviation of all selected traces.

9.1.1 Equations

If the “Enable equations”-control is checked, the x - and y -values in the plot are modified according to the equations given in the two text controls. The equation supports most standard mathematical functions like `cos(x)`, `sin(x)`, `sqrt(x)`, `exp(x)`, etc... Note that the raised operator ($\wedge$) is implemented as two multiplication signs (`**`).

In addition to the variables x and y that represents the input data, the following parameters can be used in the equations:

$p\#$:

Value of other channels in the measurement. The channels are accessed by the parameter $p\#$, where $\#$ is a number that represents the channel shown in the list below the equation controls. Note that the value will be complex if the channel represents a complex quantity; use `real(p#)`, `imag(p#)` or `abs(p#)` to get real, imaginary or the magnitude of the data.

n :

A vector with values $\{1, 2, 3, \dots, n_{tot}\}$, where n_{tot} is the number of elements in the trace.

m :

Trace number, starting with **1** for the first measured trace, which is the same as the $\#$ -parameter in the log entry list.

$m0$:

Trace number, starting with **1** for the first selected trace.

9.1.2 Physical vs. Instrument units

Select “*Tools/Plot Data in Instrument Units*” to show the data in instrument units instead of physical units. The default units (physical or instrument) can be set in the *Preferences* dialog, under “*Measurements/Units*”.

9.2 Entry list

The entry list shows the content of the log file, with each entry representing a one-dimensional trace of data. For multi-dimensional logs, the *Log Viewer* supports two different modes, which are controlled by the “*Table*”-control to the right of the list. The two modes are:

Log list:

This is the default mode, where the list contains the entries in the order that they were measured, and where the selected entries are shown directly in the graph.

Multi-column:

In this mode, the list becomes multi-dimensional, with each column representing a step dimension in the *Measurement* configuration file. The mode supports data slicing along different dimensions. The slice directions is set by the “*Slice parameter*”-control directly above the log list.

9.3 Tool bar

The toolbar contains the following buttons:

Open Log:

Open another log file.

Reload Log:

Reload the current log file, which is useful if the measurement is ongoing and new data has been added to the log.

Export:

Show the *Export Figure* dialog, for exporting the currently selected data to an image file.

Views/Save View:

Select/save the current view. See Section [Views](#) for more information.

Traces/Image:

Switch between *Trace* and *Image* plot mode.

Plot controls:

The plot controls contain tools for zooming/panning the graph, and for enabling/disabling the cursors. If in *Image* mode, there are a few extra buttons for transposing the data and enabling cross-sections and contrast controls.

Show config:

The button will show/hide an additional side bar with all the instrument configurations in the currently selected log. Below the list of quantities, there is a checkbox “*Show all quantities*”; if checked, all instrument quantities and values are shown in the list, otherwise only the quantities present as channels in the *Measurement* configuration are displayed. The “*Project*” and “*User*” controls at the bottom of the dialog allow the *Project* and *User* tags to be modified.

9.4 Multi-panel graph mode

The multi-panel graph allows multiple channels from a single log entry to be plotted in one or multiple graphs, as shown in Fig. 9.2. To enable multi-panel graph mode, select “*Views/Show multiple Graphs*” from the pull-down menu. The multi-panel mode is enabled by default if the log file contains more than one log channel. When the multi-panel graph mode is enabled, the toolbar at the top of the window contains an extra sub-menu for selecting number of figures to be shown, and for controlling whether the x- and y-axes of the figures should be synchronized or not.




Fig. 9.2. The Log Viewer dialog in multi-panel graph mode.

In multi-panel graph mode, the Y-channel control in the *Plot config* group in the upper-left corner of the window is replaced by a list of all measured channels. The channels can be assigned to one or multiple graphs by right-clicking the channel name and selecting a graph, by right-clicking one of the graphs and selecting a channel, or by dragging a channel entry onto one of the graphs. The default multi-panel graph configuration can be set in the *Preferences* dialog.

9.5 Image mode

In *Image* mode, the graph with the individual traces is replaced by an image map, as depicted in Fig. 9.3. The third-dimension data is specified by the controls in the "Third dimension"-group in the left-hand side of the dialog. In addition to specifying the data source, there are controls for performing basic signal operations along the third dimension, similar to the trace operations described in Section [PlotConfig](#). The "Third dimension"-group also contains the following buttons:

Show trace list:

If checked, a trace list is shown allowing the user to select which traces to include in the image map.

Ignore x-data:

If checked, the program will take the x-data from the first data trace and ignore the x-values of subsequent traces. This is useful for confining the representation to a square image plot for data where the x-values are changing from trace to trace.




Fig. 9.3. The Log Viewer dialog in Image mode.

If the log files contains more than two dimension, there will be a “Data selection”-list appearing under the group of “Third dimension” controls. The list allows the user to select which subset of data to show in the image plot. In image mode, there are a few extra buttons in the toolbar at the top of the window.

Transpose:

Switch X- and Y-axes in the image map.

Contrast:

Show the contrast controls, allowing the user to set the contrast range of the image. The range can be manually controlled by shifting the range cursors in the spectrum plot. The “Auto range”-button will optimize the contrast by removing outliers, while the “Full range”-button will return to full range.

X/Y cross sections:

Show the X/Y-cross sections. The position of the cross section is controlled by moving the cursors around.

9.6 Views

Views provide a way to save the current plot settings and selection of log entries, so that the current view can be easily restored. The most recently defined View will also be the preview of the log that is shown in the preview graph of the *Log browser*.

To save the current view, click the “Save view”-button in the tool bar or select “Views/Save view...” in the menu, define a name of the View and click the “OK”-button.

To restore a previously saved view, select the View to show from the “View”-control in the toolbar. Views can be renamed or deleted from the *Edit View*-dialog, which is accessed by selecting “Views/Edit views...” in the menu bar.

9.7 Exporting data

Data can be exported to a text file, to a *Matlab* “.mat” file, or as an image. The export options are available from the “File”-menu.

9.7.1 Exporting to Image




Fig. 9.3. The Export Image dialog makes it easy to generate publication-quality figures.

The currently selected plot can be saved as an image file by either selecting “Tools/Save Screenshot”, or saved to the clipboard by selecting “Tools/Copy Graph”. To modify the labels and the style of the image, select “Export/Export as Image...” or click the *Export* button in the toolbar. This will open the *Export Figure* dialog, which makes it possible to modify the figure axes and labels to render publication-quality figures (see Fig. 9.3). The resulting image can be copied to the clipboard or saved in JPEG, PNG, SVG or Adobe PDF format.

9.7.2 Exporting to Text

When exporting the text, a dialog will open allowing the user to define what to export. The following options are available:

Data to export:

Determines whether to export all or only the traces currently selected in the *Log Viewer*.

Include header with log information:

If checked, a header with log info will be included at the top of the text file.

Include separate x-data with each trace:

If checked, each data entry will contain two rows of data, one for x and one for y .
If unchecked, the first row is x -data, and all the following rows are y -data.

Include data for third dimension:

If checked, data for a third dimension is added to the end of the text file. The channel for which the third dimension data is taken is defined in the control below the check box.

9.7.3 Exporting to Matlab

The *Matlab* export will export the selected traces to a single “.mat” file. The data is saved into a *Matlab* struct; the data structure is shown when opening the file in the *Matlab* workspace browser.

9.7.4 Custom Export

The *Custom Export*-function allows the user to export data into a custom format. Note that the custom export function always exports the raw data, without applying any operations such as smoothing, FFT, etc. The following options are available when exporting custom data:

Data to export:

Determines whether to export all or only the traces currently selected in the *Log Viewer*.

Custom script:

Path to python file containing the custom `exportData` function.

The custom export functionality needs to be implemented in a python function called `exportData`, which should be located in a separate python file (.py). An example of a custom export script can be found in the file *ExportScript.py* in the *Script* folder of the main program directory (see Section [Folders](#) for an overview of folder locations). The function definition of the `exportData`-function must have the following format:

```
def exportData(file_name, step_data, log_data, step_name, log_name,
```

```
step_unit, log_unit, comment='')
```

file_name:

Output path for the exported data.

step_data:

Data for stepped channels. The data is defined in a nested python list of 1-d *numpy* arrays (one for each trace). The first index is the step channel number as defined in the *Measurement* dialog step list, the second index is the trace number. For example, to access the data for the innermost step channel and the third trace, use `step_data[0][2]`.

log_data:

Data for log channels, defined in the same way as the `step_data`.

step_name:

List of strings defining the step channel names.

log_name:

List of strings defining the log channel names.

step_unit:

List of strings defining the step channel units.

log_unit:

List of strings defining the log channel units.

comment:

Log comment.

10 Preferences

To access the Preferences dialog, select *Preferences...* from the *Instrument Server* tray icon menu, or “*Edit/Preferences...*” from the pull-down menu. The dialog has the following sections and settings:




Fig. 9.4. Dialog for setting preferences.

10.1 Folders

The *Folders* section defines the location of various folders used by the program.

Database folder:

Main database folder for saving data from the *Measurement* program. Default value is “<User home directory>/Labber/Data”.

Local drivers:

Folder containing user-defined instrument drivers.

Local optimizers:

Folder containing user-defined optimizer functions.

10.2 Server

The *Server* section contains settings related to the *Instrument Server*.

Start server on program startup:

If True, the server starts listening for incoming connections when the program starts up. If False, the user has to start the server manually. Default value is True.

TCP port:

TCP port used for communication between the *Instrument Server* and the clients. Default port is 9406.

Notification TCP port:

TCP port used for sending notifications between the various program parts. The communication only occurs on the local computer. The default port is 9407.

Data transfer format:

Format use for data transfer over the network. Binary is faster, whereas text is human-readable and better for debugging purposes. Default is Binary.

Server timeout:

Maximum waiting time before the server returns an error. This value should be reasonable long, in case an instrument takes a long time to perform an operation. Default value is 1,000,000 seconds.

Restrict client IP addresses:

Restrict allowed clients according to the list defined below. Default is True.

Allowed clients:

List of IP numbers of allowed clients. Request from computers with IP numbers outside the list will be rejected. Note that it is possible to define wild cards, for example `"192.168.*"` will allow connections from any client with IP starting with 192.168. Default value is *localhost*, which only allows connections from the same computer that is running the server.

Allow instruments to be controlled from driver configuration window:

If True, instrument settings can be updated directly from the driver configuration window while a driver is running. If False, the controls are grayed out once the driver is started, which makes it less likely that incorrect/too large instrument values are outputted by mistake. Default is True.

Keep instrument drivers running after the measurement ends:

Starting up a driver may take a few seconds, depending on system. Therefore, stopping and starting the driver between measurements may slow down experiments, which can be avoided by keeping the driver running after a measurement ends. Default is True.

Change background color for active instruments:

Change background color of driver dialog for active instruments, to highlight that changes to any parameter of the instrument driver window will directly update the instrument hardware. Default is True.

Change background color for instruments in Measurement dialog:

Use different background color for instrument configuration dialogs in the Measurement program than in the Instrument server, to make it clear which program the dialog belongs to. In the *Instrument Server*, the instrument configuration dialog is used to *directly* control the hardware settings, meaning that any changes to the dialog will directly affect the state of the hardware. In contrast, in the *Measurement* program the dialog is used to set up a configuration that will be used in a specific *Measurement*, but no changes are made to the hardware until the measurement is started. To avoid confusion, if this setting is True the *Instrument driver* configuration windows have a different background color when opened within the *Measurement* program and in the *Instrument Server*. Default is True.

Instrument log level:

Amount of information to log when performing instrument communication. The log can be viewed by selecting “Log/View Instrument Log...” in the *Instrument Server* menu bar. Default value is Basic.

Network log level:

Amount of information to log when performing network communication. The log can be viewed by selecting “Log/View Network Log...” in the *Instrument Server* menu bar. Default value is Basic.

10.3 Measurement

This section contains settings related to the *Measurement* program.

Sort step items before starting Measurement:

If checked, step items are sorted according to instrument type before starting a Measurement.

Default units, step sequences:

Default units when defining a new step sequence in the *Measurement Setup* dialog.

Default units, viewing data:

Default units when viewing data in the *LogBrowser* and the *LogViewer*.

Default sweep units:

Set if sweep rates should be defined in terms of rate per second or rate per minute in the *Instrument Server* and *Measurement Editor* programs. If set to *Instrument default*, the program will use the default units defined in the settings of each instrument driver (see Section [SweepDriver](#)). Note that this setting only affects the sweep units shown in the dialog windows, the sweep units used within a particular instrument driver implementation is always set by the configuration file of the driver (see Section [SweepDriver](#)).

Graph refresh interval:

Refresh interval for graph. Use larger values if the user interface becomes unresponsive. Default value is 80 ms.

Default live colormap:

Default colormap for viewing image data in the live graph shown during measurements.

10.4 Log Viewer

Default colormap:

Default colormap when viewing data as images in the *LogBrowser* and the *LogViewer* dialogs.

Default cursor type:

Default cursor type in all graphs.

Default complex representation:

Default format for representing complex scalar data.

Default complex representation, vector:

Default format for representing complex vectors, typically from instruments such as spectrum analyzers and vector analyzers.

Default panel configuration, 2 channel:

Default multi-panel graph configuration for showing two log channels.

Default panel configuration, 3 channel:

Default multi-panel graph configuration for showing three log channels.

Default panel configuration, 4 channel:

Default multi-panel graph configuration for showing four log channels.

Save current view when closing Log Viewer:

Automatically save current view when closing the Log Viewer.

10.5 Logger

Logger folder:

Database folder for saving logging data from the Logger program.

Number of points in Acquire graph:

Number of points shown in the live logger graph.

Alarm de-activation range:

Range at which an out-of-range alarm de-activates.

Dark mode:

If checked, the visualizer will plot data on a dark background.

Refresh interval in Logger Visualize:

Data refresh interval in Logger visualize.

10.6 Advanced

Python distribution:

Path to custom Python executable. The Python distribution must be running Python 3.6 or later. Leave blank to use the built-in Labber Python distribution. For Windows, pick the executable `pythonw.exe` instead of `python.exe` to avoid creating a console window for each driver process. For more information, see Section [PythonDistExternal](#).

Temporary items:

Folder for storing settings and temporary items. Do not alter this item unless having good reasons for doing so.

Show error if setting the value of an inactive quantity:

If unchecked, the program will not show an error if trying to set the value of an inactive quantity.

Send status updates to clients:

Send status updates from Instrument server to log clients during slow operations such as sweeping.

Interval for checking swept instruments:

Time interval between checks when testing if a swept instrument has reached the final value.

VISA library:

Path to VISA library. Leave blank to use default library.

Delay for wait dialog:

Shortest delay time for showing the wait dialog. Default value is 2 seconds.

Show error dialog in script mode:

If unchecked, no error dialog will be shown if an error occurs during a scripted Measurement. This can be useful if no user interactions is required to handle errors.

Run queued experiments in separate process:

If checked, queued measurements will run in a separate instance of the Measurement program. This may cause conflicts if queued experiments and measurement from the user interface are started at the same time.

11 Scripting

This chapter describes how to write scripts to perform sequence of experiments. Scripting is useful for running multiple experiments after each other, or for defining sequences where some properties of a measurement is updated depending on the result of a previous measurement.

In a typical scripting setup, the user would first create a number of *Measurement* configurations using the standard *Measurement* configuration dialog, and then save those configurations files to a folder on disk. The script would then be programmed to execute those configurations, either as they are or by first updating one or multiple parameters of the *Measurement* configurations.

The most basic way of implementing scripting is to call the *Measurement* program with command-line arguments, as described in Section [console](#) below. The advantage of this method is that one can use any programming language that supports calling an external program for writing the script, the disadvantage is that the function calls can become rather long and difficult to read. If you plan to script experiments using the programming language *Python*, there are a number of helper functions that will simplify the procedure. These helper functions are described in Section [scriptPython](#).

11.1 Console options

In addition to the user-interface based *Measurement* configuration dialog, it is also possible to start experiment from the command line. The program is called `Measurement-Console.exe`, and is located in the main program folder (see Section [Installation](#) for the folder structure). The command-line arguments are:

```
Measurement-Console.exe [-h] -i INPUT_PATH [-u CHANNEL VALUE TYPE]
                        [-m CHANNEL] [-o OUTPUT_PATH] [-e EXPORT_PATH] [-r CHANNEL]
```

-h, -help:

Show a help message and exit

-i INPUT_PATH:

Path to the measurement configuration file to execute or rearrange.

-u CHANNEL VALUE TYPE:

Update the step item `CHANNEL` with a new value. The `TYPE`-argument defines what property of the step item to update, and must be one of `SINGLE`, `START`, `STOP`, `CENTER`, `SPAN`, `STEP`, `N_PTS`. Note that scripted measurements do not raise an error if updating an inactive step item. Instead, the step item is automatically switched to the new step type. For example, if the original step type is `SINGLE`, and the user updates the `START` value, the step type is changed to `START-STOP`. Note that it is up to the user to ensure that all other relevant quantities are updated as well (in the example, the `STOP` value and the `STEP` or `N_PTS` value).

-m CHANNEL:

Specifies the primary channel name. Values of all other updated channels will be defined by look-up tables relative to the primary channel values.

-o OUTPUT_PATH:

Specifies the path of the output log file. If not given, the data will be saved to the input measurement configuration file.

-e EXPORT_PATH:

After completed measurement, export the last trace to the specified text file. Any previous contents in that file will be overwritten. This is useful for creating scripts where future measurements depend on the results of previous measurements.

-r CHANNEL:

Re-arrange a log with N 1-dimensional entries of length M to a 2-dimensional log with dimensions (N, M). The **CHANNEL** determines which data to use when defining the second dimension. It is also possible to rearrange into a multi-dimensional log by specifying multiple channels, but if so lists of step values for each dimension need to be specified as well. For example, to rearrange a log with 6 entries into a multi-dimensional log with 3*2 entries, use **-r "Channel 1" "1.0, 2.0, 3.0" "Channel 2" "1.0, 2.0"**. Note that the internal order of the new dimensions is defined by the order in which they appear in the step list of the original *Measurement* configuration file, *not* by the order they are listed after the **-r** command. Also, note that no measurement will be performed when running the program with this option.

If no arguments are given, the program will open the standard user interface window for configuring the experiment.

11.2 Scripting using Python

The Python scripting helper functions are part of the *Labber API*, located in the *Script* folder of the main program directory (see Section [Folders](#) for an overview of folder locations). The *Labber API* can be installed into a local Python environment by using pip. For details on how to do this, see section [A1.1 PythonAPI](#).

The helper functions in the *ScriptTools* module are designed for repeatedly performing a number of *Measurements* that each contain one-dimensional sweeps, and where one or multiple parameters of the *Measurement* configurations are updated between each measurement. The functions are best explained by an example, which we'll take from the domain of superconducting qubits. For the purpose of this example, we can view the qubit as a slightly anharmonic oscillators whose frequency tunes with applied magnetic flux. The qubit is read out by coupling it to microwave resonators, and the coupling is arranged in a way that changing the qubit frequency will cause a slight shift of the resonator frequency.

Now, say that we want to probe the qubit frequency as a function of applied flux. The difficulty is that the changing the flux will affect both the qubit and the resonator frequencies, which means that we can not use a fixed-frequency read-out tone. Instead, we need to implement the following procedure:

1. **Set new magnetic flux value**
2. **Measure resonator**
3. **Find resonance frequency f_0 of resonator**
4. **Measure qubit, while keeping the resonator at f_0**
5. **Repeat for all values of magnetic flux**

The file *ExampleScript.py* in the *Script* folder contain an example script for performing the sequence described above. The script assumes that the user has created two *Measurement* configurations, one for measuring the resonator, and one for measuring the qubit, and that both *Measurement* configurations have a single-valued step item called 'Flux bias' that control the magnetic flux.

12 Instrument drivers

This chapter describes the definition of instrument drivers and instructions for how to create custom drivers. The general driver structure is visualized in Fig. 12.1.

The *Communication* part describes the interface and address used for communication, and is normally handled by the *Labber Instrument server*. The *Model and options* part provides a way to enable/disable certain features of a driver depending on the instrument model/installed options. Finally, the list of *Quantities* defines all properties and settings available on the instrument.




Fig. 12.1. Structure of an Instrument driver. The Communication part is handled by the Instrument server, while the Model and Options and the Quantities are defined in the driver configuration file.

12.1 Driver definition files

The driver definition file is a file specifying the instrument name, vendor, model and options as well as a list of quantities. For basic instrument drivers, where the value of each *quantity* can be set or read using a single text command over GPIB, serial, USB or ethernet using the VISA protocol, all information about the instrument and the communication is contained in the definition file. For more advanced drivers, for example network analyzers that capture vector data, the driver definition file needs to be complemented with *Python* source code for implementing the more advanced instrument operations (see Section [PythonDriver](#)). The Python code must be Python 3 compatible.

The driver definition files provided by *Labber* are located in the “*Drivers*” subfolder under your local installation. The definition files are plain text files using the *INI* file format, which consists of a number of “sections”, each containing a list of “properties”. The driver file requires implementing sections for *General settings*, *Model and options* and *VISA Settings*, see below for more information about each section. It is recommended to use one of the already present driver configuration files as a template. For an example of creating an instrument driver from scratch, see Section [PythonDriver](#).

When creating a new driver, the definition file should be placed in the *local* driver folder (the folder named “*Local drivers*” in the *Preferences* window), instead of the global one

(the installation location). This allows the user's own drivers to be kept separately from the drivers provided by *Labber*, and it also prevents drivers written by the user from being deleted when updating the *Labber* program to a newer version.

The *INI* configuration file can be placed directly in the “*Drivers*” folder, or within a subfolder of that directory. Using a subfolder is the recommended approach, since it gives a natural place to store extra files related to the driver.

Note that even when making additions/changes to an existing driver from the global folder, the best practice is to copy that driver file from the global folder to the local folder, and only make changes to the local version. If drivers with the same names exist in both the local and the global driver folders, *Labber* will always use the driver in the “*Local drivers*”-folder.

12.1.1 Signal Generators and Signal Analyzers

Signal Generators and *Signal Analyzers* are drivers that are used to generate or analyze waveforms. The drivers do not perform any instrument communication, which means that the *Model and options* and the *VISA Settings* parts of the *INI* file do not need to be defined. To define that a driver is a *Signal Generators* or a *Signal Analyzers*, set the corresponding item in the *General settings*-part of the *INI* file as described below. See Section [Signals](#) for more information about how *Signal Generators* or a *Signal Analyzers* are used in an experiment.

12.1.2 General settings

The *General settings*-section define name and version of the driver. Note that it is the *name* property in this section that sets the driver name, not the name of the driver definition *INI* file.

name:

The name is shown in all the configuration windows.

version:

The version string should be updated whenever changes are made to this config file.

driver_path:

Name of folder containing the code defining a custom driver. Do not define this item or leave it blank for any standard driver based on the built-in VISA interface.

interface:

Pre-defined communication interface for instrument, default is `GPIB`.

Valid values are `GPIB`, `TCPIP`, `USB`, `PXI`, `Serial`, `VISA`, `Other`, `None`.

address:

Pre-defined address for instrument, default is an empty string.

startup:

Pre-defined startup option for instrument, default is `Set config`.

signal_generator:

Set to *True* if driver is a *Signal Generator*. Default is *False*.

signal_analyzer:

Set to *True* if driver is a *Signal Analyzer*. Default is *False*.

controller:

Set to *True* if driver is a *Controller*. Default is *False*. For more information, see Section. [ControllerDriver](#) below.

support_hardware_loop:

Set to *True* if driver supports hardware looping. Default is *False*.

support_arm:

Set to *True* if driver supports hardware arming. Default is *False*.

use_32bit_mode:

Set to *True* if driver should run in a 32-bit Python environment. Default is *False* (run in 64-bit). For more information, see Section [PythonDist](#) below.

12.1.3 Model and options

The *Model and options*-section provides a way to enable/disable certain features of a driver depending on the instrument model/installed options.

`model_str_1, model_str_2, etc:`

List of models supported by the driver.

`check_model:`

If *True*, the driver checks the instrument model id at startup (*True* or *False*). The model is checked by sending the `model_cmd` command (see below) over the VISA interface. Default is *False*.

`model_cmd:`

Command used to check the instrument model. Default command is `*IDN?`.

`model_id_1, model_id_2, etc:`

Model strings expected to be returned by the instrument by the `*IDN?` call. If not defined, the program assumes `model_str_1, model_str_2, etc` as default values

`option_str_1, option_str_2, etc:`

List of available instruments options. The options are shown as checkbox controls in the driver configuration window.

`check_options:`

If *True*, the driver checks the installed instrument options at startup (*True* or *False*). The option is checked by sending the `option_cmd` command (defined below). Default is *False*.

`option_cmd:`

If `check_options` is set to *True*, define command for getting the options from the instrument.

option_id_1, option_id_2, etc:

If `check_options` is set to *True*, supply valid id option strings that the instrument returns when sending the `option_cmd`. The list of `option_id` should match the elements in the list `option_str`.

12.1.4 VISA Settings

This section contains configuration of the VISA protocol. The VISA protocol enables text-based communication with instruments over GPIB, USB, serial and ethernet interfaces.

use_visa:

Enable or disable communication over the VISA protocol (*True* or *False*). If *False*, the driver will not perform any instrument operations (unless there is a custom *Python* driver, see Section [PythonDriver](#)).

reset:

Reset the interface (not the instrument) at startup (*True* or *False*). Default is *False*.

query_instr_errors:

Query instrument errors (*True* or *False*). If *True*, every command sent to the device will be followed by an error query. This is useful when testing instruments, but may degrade performance by slowing down the instrument communication.

error_bit_mask:

If `query_instr_errors` is *True*, set bit mask for checking status byte errors (default is 255, include all errors). The bits signal the following errors:

- 0: Operation
- 1: Request control
- 2: Query error
- 3: Device error
- 4: Execution error

5: Command error

6: User request

7: Power on

error_cmd:

Command string to be sent to instrument when querying for instrument error messages.

init:

Initialization commands are sent to the instrument when starting the driver. `*RST` will reset the device, `*CLS` clears the interface.

final:

Final commands sent to the instrument when closing the driver.

str_true: _

String used for sending boolean *True* to the instrument, default is `1`.

str_false:

String used for sending boolean *False* to the instrument, default is `0`.

str_value_out:

Conversion string used for converting value to string to be sent to the instrument. Default is `%.9e`, which creates 9-digit string using exponential notation. To create strings with floating-point notation, use `%.9f` instead.

str_value_strip_start:

Number of characters to strip from the beginning of the string returned from the instrument, before trying to convert to a number. Default is `0`.

str_value_strip_end:

Number of characters to strip from the end of the string returned from the instrument, before trying to convert to a number. Default is `0`.

always_read_after_write:

If **True**, the program will automatically read the response from the instrument after each write command. Useful for instruments that always reply to all commands. Default is **False**.

The following entries are optional, they provide detailed settings for the communication interface. Note that the values provided in the *INI* file will be the default setting for the driver, but the user can always change the settings by going to the *Communication* settings of the instrument driver user interface and clicking “*Show advanced interface settings*”.

timeout:

Time (in seconds) before the timing out while waiting for an instrument response. Default is 10 seconds.

term_char:

Termination character used by the instrument, valid values are **Auto**, **None**, **CR**, **LF**, **CR+LF**.

send_end_on_write:

Assert end during transfer of last byte of the buffer

suppress_end_on_read:

Suppress end bit termination on read

baud_rate:

Communication speed for serial communication. Default is 9600.

data_bits:

Number of data bits for serial communication. Default is 8.

stop_bits:

Number of stop bits for serial communication. Default is 1, possible values are 1, 1.5 and 2

parity:

Parity used for serial communication, possible values are `No parity`, `Odd parity`, `Even parity`.

gpiib_board:

GPIB board number. Default is 0.

gpiib_go_to_local:

Make GPIB instrument automatically go to local after closing. Default is `False`.

tcipip_specify_port:

Use specific TCPIP socket port. Default is `False`.

tcipip_port:

TCPIP socket port. Only relevant if `tcipip_specify_port` is `True`.

12.2 Quantities

All quantities are defined in separate sections, with the name of the quantity given by the section header. The properties of a quantity are defined by a number of keywords, see below for a list the possible options. Only the `datatype` keyword is mandatory, the other ones are optional.

datatype:

The data type should be one of `DOUBLE`, `BOOLEAN`, `COMBO`, `STRING`, `COMPLEX`, `VECTOR`, `VECTOR_COMPLEX`, `PATH` or `BUTTON`. Only `DOUBLE`, `BOOLEAN` and `COMBO` datatypes can be stepped in a measurement. The `BUTTON` datatype does not have an associated value, and can therefore not be controlled from the *Measurement* program. It is typically used to manually force an instrument to perform a certain task.

label:

Label shown next to control in user interface. If not specified, the label defaults to the name of the quantity.

unit:

Unit for the quantity.

def_value:

Default value.

tooltip:

Tool tip shown when hovering the mouse over the control in the driver GUI.

low_lim:

Lowest allowable value. Defaults to `-INF`.

high_lim:

Highest allowable values. Defaults to `+INF`.

x_name:

X-axis label for a vector data. Only valid if `datatype` is `VECTOR` or `VECTOR_COMPLEX`.

x_unit:

X-axis unit for a vector data. Only valid if `datatype` is `VECTOR` or `VECTOR_COMPLEX`.

combo_def_1, combo_def_2, ...:

Options for a pull-down combo box. Only used when `datatype` is `COMBO`.

group:

Name of the group where the control belongs.

section:

Name of the section where the control belongs.

state_quant:

Quantity that determines this control's visibility.

second_state_quant:

A second quantity that determines this control's visibility. This is an 'AND' operation: if a control has a state_quant and a second_state_quant, both need to be True for it to appear

state_value_1,state_value_2, ...:

Values of "state_quant" for which the control is visible.

second_state_value_1,second_state_value_2, ...:

Values of "second_state_quant" for which the control is visible, if "state_quant" is also True

model_value_1, model_value_2, ...:

Values of "model" for which the control is visible. The value must match one of the models defined in the *Model and Options*-section described above.

option_value_1, option_value_2, ...:

Values of "option" for which the control is visible. The value must match one of the options defined in the *Model and Options*-section described above.

permission:

Sets read/writability, options are BOTH, READ, WRITE or NONE. Default is BOTH.

show_in_measurement_dlg:

This setting is optional. If True, the quantity will be automatically shown when adding the instrument to a *Measurement* configuration. This is useful for instrument that contain a lot of quantities, but where most are not likely to be stepped in a measurement.

set_cmd:

Command used to send data to the instrument. Put "<*>" where the value should appear. If "<*>" does not occur in the string, the value will be added after the command.

get_cmd:

Command used to get the data from the instrument. Default is `set_cmd?`.

cmd_def_1, cmd_def_2, ...:

List of strings that define what is sent to/read from an instrument for a quantity that is defined as a list of multiple options. Only used when `datatype` is `COMBO`.

See Section [SweepDriver](#) for a list of extra properties that need to be defined for instruments that support sweeping.

The *Instrument Server* uses the list of quantities to create the controls in the driver dialog window, as shown in Fig. 12.2.




Fig. 12.2. An Instrument driver dialog, shown together with the corresponding instrument definition file.

12.3 Custom drivers - Python code

Custom drivers are required when single-line command strings `get_cmd` and `set_cmd` as defined in the instrument definition file are too simple to read or write a value to an instrument. This is often the case for instruments like network analyzers or oscilloscopes, which contained vector-valued quantities that depend in complicated ways on other settings of the instrument.

The process of creating a custom driver is best described by an example. We are going to create a driver that generates a sinusoid, but without doing any actual instrument communication (the driver will be a *Signal Generator*, as described in Section [Signals](#)). For

an example involving instrument communication, see the drivers for one of the network analyzers or oscilloscopes in the *Instrument Drivers* folder.

12.3.1 Creating the driver definition file

Every driver, even the custom ones, require a definition file. We start with the *General settings*-section:

```
[General settings]

# The name is shown in all the configuration windows
name: Simple Signal Generator

# The version string
version: 1.0

# Name of folder containing the code defining a custom driver
driver_path: SimpleSignalGenerator

# Define that the driver is a Signal Generator
signal_generator: True
```

Note that we define the `driver_path`: this signals that there is a custom driver available for this instrument. When starting the driver, the *Instrument Server* will look for the Python file “SimpleSignalGenerator/SimpleSignalGenerator.py” in the *Instrument Drivers* folder, or for the file “SimpleSignalGenerator.py” in the folder where the *INI* configuration file is located. See Section [PythonCode](#) below for more information on how to implement the code for custom drivers.

This particular instrument driver does not do any instrument communication and therefore does not have any model or option definitions, so we can skip the *Model and options*-section and the *VISA settings*-section.

Next, we need to define the quantities of the driver. For this example, we want to be able to define the amplitude, frequency and phase of the signal to be generated. In addition, we want to add the option of adding white noise to the signal.

```
[Frequency]
datatype: DOUBLE
unit: Hz
def_value: 10.0

[Amplitude]
```

```

datatype: DOUBLE
unit: V
def_value: 1.0

[Phase]
datatype: DOUBLE
unit: deg
def_value: 0.0

[Add noise]
datatype: BOOLEAN
def_value: False

[Noise amplitude]
datatype: DOUBLE
unit: V
def_value: 0.1
state_quant: Add noise
state_value_1: True

[Signal]
datatype: VECTOR
permission: READ
x_name: Time
x_unit: s

```

Note that the *Noise amplitude* quantity will only be visible if *Add noise* is *True*. The last quantity (“*Signal*”) represents the signal we want to generate in the *Python* code. The `permission` of this quantity is set to `READ`, to indicate that this quantity can only be read, not written.

For your convenience, this example driver *INI* definition file and the corresponding *Python* code are available under *Examples* in the *Instrument Drivers* folder.

12.3.2 Implementing the *Python* code

Once the *INI* file has been created, we need to implement the *Python* code that generates the signal. The code should define a subclass of either the `InstrumentDriver.InstrumentWorker` or the `VISA_Driver` class, depending on if the driver will use the *VISA* protocol for communication or not. The `VISA_Driver` class is a subclass of `InstrumentDriver.InstrumentWorker`, and details for how to subclass the `VISA_Driver` is described in Section [SubClassVISA](#) below. *Labber* is running *Python* 3 for all instrument drivers, make sure that all code is *Python* 3 compatible. See Section [PythonDist](#) below for more information about the *Python* distribution.

The new class should re-implement the four functions `performOpen`, `performClose`, `performSetValue` and `performGetValue`, which are called when an instrument is started, stopped, and called for setting or getting an instrument value, respectively. To describe the procedure, we create a *Python* class for the *Simple Signal Generator*-example shown above:

```
import InstrumentDriver
import numpy as np

class Driver(InstrumentDriver.InstrumentWorker):
    """ This class implements a simple signal generator driver """

    def performOpen(self, options={}):
        """Perform the operation of opening the instrument connection"""
        pass

    def performClose(self, bError=False, options={}):
        """Perform the close instrument connection operation"""
        pass
```

As described in the previous section, the *Simple Signal Generator* is only for demonstration purposes and will not involve any actual instrument communication, so we subclass `InstrumentDriver.InstrumentWorker` instead of `VISA_Driver`. In this example, the functions `performOpen` and `performClose` don't do anything.

The code for the more interesting functions `performSetValue` and `performGetValue` follow below:

```
def performSetValue(self, quant, value, sweepRate=0.0, options={}):
    """Perform the Set Value instrument operation. This function should
    return the actual value set by the instrument"""
    # just return the value
    return value

def performGetValue(self, quant, options={}):
    """Perform the Get Value instrument operation"""
    # proceed depending on quantity
    if quant.name == 'Signal':
        # if asking for signal, start with getting values of other controls
        amp = self.getValue('Amplitude')
        freq = self.getValue('Frequency')
        phase = self.getValue('Phase')
        add_noise = self.getValue('Add noise')
        # calculate time vector from 0 to 1 with 1000 elements
        time = np.linspace(0,1,1000)
        signal = amp * np.sin(freq*time*2*np.pi + phase*np.pi/180.0)
        # add noise
        if add_noise:
```

```

        noise_amp = self.getValue('Noise amplitude')
        signal += noise_amp * np.random.randn(len(signal))
        # create trace object that contains timing info
        trace = quant.getTraceDict(signal, t0=0.0, dt=time[1]-time[0])
        # finally, return the trace object
        return trace
    else:
        # for other quantities, just return current value of control
        return quant.getValue()

```

The functions `performSetValue` and `performGetValue` take a `quant` object as a first parameter. The object represents the quantity to be read/set, and all properties of the quantity (as defined in the *INI* configuration file) can be accessed from the object's data members. This is used in the `performGetValue`-function, where the object variable `quant.name` is accessed to find out which quantity to read. See Section `quantObj` below for more info about the `quant` objects.

The `options` variable present in both `performSetValue` and `performGetValue`-definitions is a *Python* dictionary that contains additional options for setting/getting a value. It is used to provide a way to determine if a driver is called multiple times within a single step of a *Measurement* (see functions `isFirstCall` and `isFinalCall` in the list of driver helper functions in Section `driverObj` below).

When implementing `performSetValue`, note that the value of the quantity in the configuration is set *after* this method returns `value`. Typically `performSetValue` should return the argument `value` without modification, but in certain cases it may be useful to perform some processing on `value` and return a different value. As a result of this behavior, it is possible to access the *previous* value of the quantity stored in the configuration via the helper functions `quant.getValue()` or `self.getValue(quant)` (see below). If it is required that these methods return the *new* value of the quantity, the configuration can be updated within `performSetValue` via the helper function `quant.setValue(value)` (see below). This consideration is particularly important when performing operations only after multiple quantities have been set (using the `isFinalCall` method), as by default all quantities will have been updated in the configuration, *with the exception of the last one*.

12.3.3 Helper functions for `quant` objects

The `quant` object represents an instrument quantity, and it provides a few helper functions that are useful when writing drivers:

`quant.getValue():`

The function returns the current value of the quantity. Note that it will just return the local value stored in the driver, no instrument communication is performed when calling this function.

`quant.getValueIndex(value=None):`

The function returns the value as an index number, only useful for quantities with `datatype=COMBO`. If `value=None`, the function will return the local value stored the driver. Note that no instrument communication is performed when calling this function.

`quant.setValue(value, rate=None):`

The function sets the current value of the quantity. Note that it will just update the local value stored in the driver, no instrument communication is performed when calling this function.

`quant.getTraceDict(y, x0=0.0, dx=1.0, x1=None, x=None, logX=False):`

Returns a python dictionary containing the numpy array `y`, together with additional x-scale info. The x-scale information can be supplied either as start value and step size `(x0, dx)`, as start and stop values `(x0, x1)`, or as a full vector (input parameter `x`, must have same length as `y`). If using the start/stop notation `(x0, x1)`, it is possible to set `logX` to `True` to create a trace with logarithmic interpolation between the start/stop values. These dictionaries are used to pass waveform data between drivers with vector-valued quantities, like *Signal Generator* and *Signal Analyzers*.

`quant.getCmdStringFromValue(value=None):`

Convert the input value to a string formatted for sending to the instrument. If the input parameter `value` is `None`, the current value is used.

quant.getValueFromCmdString(sValue):

Inspect the input string `sValue` coming from the instrument and return a numerical value.

12.3.4 Helper functions for `driver` objects

The base driver object `InstrumentDriver.InstrumentWorker` provides the following helper functions. The `options` variable present in the functions `isFirstCall` and `isFinalCall` is a *Python* dictionary with additional options that is passed to the `performSetValue` and `performGetValue`-functions when calling the driver from outside.

getName():

Return name of instrument, as defined in the user-interface dialog.

getInterface():

Return instrument interface, as defined in the dialog. The interface type is one of `GPIOB`, `TCPIP`, `USB`, `Serial`, `VISA`, `Other`, `None`.

getAddress():

Return address of instrument, as defined in the user-interface dialog. This is function can be used to determine when opening communication to an instrument.

getCommunicationCfg():

Return communication configuration as a dictionary, with the following keys: `Timeout`, `Term. character`, `Send end on write`, `Suppress end bit termination on read`, `Baud rate`, `Data bits`, `Stop bits` and `Parity`. The configuration items are described in Section [CommunicationCfg](#) above.

getValue(quant_name):

The function is used to access the current local value of any quantity of the driver. The function is used repeatedly in the example above for getting the amplitude, frequency and phase when creating the sinusoid. Note that if you access a quantity by `self.getValue` or `self.getValueArray` within `performSetValue`,

it will return the previous value of the quantity. This is by design, in case access to the old state is required before setting the new state. To update to the current value include `quant.setValue(value)` at the start of `performSetValue`.

`getValueArray(quant_name):`

Same as above, but will return current value as a numpy array instead if the quantity is vector-valued. Otherwise, it'll return an empty numpy array. Note that if you access a quantity by `self.getValue` or `self.getValueArray` within `performSetValue`, it will return the previous value of the quantity. This is by design, in case access to the old state is required before setting the new state. To update to the current value include `quant.setValue(value)` at the start of `performSetValue`.

`getValueIndex(quant_name):`

Get value of quantity as numerical index. Only useful for quantities with `datatype=COMBO`.

`getCmdStringFromValue(quant_name):`

Get command string for current value of quantity with name `quant_name`. See `quant.getCmdStringFromValue` in section above for more info.

`setValue(quant_name, value, sweepRate=None):`

The function is used to set the local value of any quantity of the driver. No hardware communication will take place; to actually set the instrument value, use the function `sendValueToOther` defined below.

`readValueFromOther(quant_name, options={}):`

The function will read the value of another quantity from the instrument. In contrast to the `getValue` mentioned above, this function will perform actual hardware communication to retrieve the current value from the instrument. The function will return the updated value.

`sendValueToOther(quant_name, value, sweep_rate=0.0, options={}):`

The function will communicate with the hardware to update the value to another quantity of the instrument. The function will return the updated value.

getModel():

Get model string.

setModel(model_name):

Set model string.

getOptions():

Get list of strings describing installed options.

setInstalledOptions(list_of_options):

Set list of strings describing installed options.

isConfigUpdated(bReset=True):

Returns true if any non-read-only quantity of the instrument has been updated since the last call to this function where `bReset` was `True`.

isFirstCall(options):

If a driver is used in a *Measurement* and there are multiple quantities of that driver that will be updated within a single step, it can be advantageous to delay outputting data to an instrument until all local driver quantities have been updated. This function returns `True` if the current call is the first one within the current measurement step.

isFinalCall(options):

Same as above, but returns `True` if the current call is the last one.

isStopped():

Return `True` if the user stopped the measurement. If the instrument communication is expected to take a long time, it's recommended to periodically call this function to ensure that the driver remains responsive to user interaction.

isHardwareTrig(options):

Return `True` if the caller is in hardware trig mode.

isHardwareLoop(options):

Return `True` if the caller is in hardware loop mode.

getHardwareLoopIndex(options):

Get the current hardware loop index. The function returns a tuple (`index`, `n_pts`), where `index` is the index for the current call, and `n_pts` is the total number of points of the hardware loop.

log(message, level=20):

Log a message to the instrument logger. The log level is an integer ranging from 30 (warning, always shown) to 10 (debug, only shown in debug mode).

wait(wait_time=0.05):

Pause execution and put the process to sleep for the given time (in seconds).

getValueFromUserDialog(value=None, text='Enter value:', title='User input'):

Show user interface dialog to ask the user for an input value. The function returns the value entered by the user.

reportStatus(message):

Report status update to the *Instrument Server* and connected clients. The argument `message` should be a string.

reportProgress(quant, progress):

Report progress update when setting/getting the value of the quantity `quant` to the *Instrument Server* and connected clients. The argument `progress` should be a floating point value between 0.0 and 1.0. The function is used to provide feedback to the user when performing slow instrument operations, for example when sweeping a magnetic field.

reportCurrentValue(quant, value):

Report current value of the quantity `quant` to the *Instrument Server* and connected clients. The function is used to provide feedback to the user when performing slow instrument operations, for example when sweeping a magnetic field.

12.3.5 Testing the driver

As stated previously, this example driver *INI* definition file and the corresponding *Python* code are available under *Examples* in the *Instrument Drivers* folder. To test the driver, move the *INI* file and the folder with the *Python* code to reside directly in the *Instrument Drivers* folder. Next, start the *Instrument Server* and add a new instrument. If the driver is defined properly, the *Simple Signal Generator* should show up in the instrument driver list. Select the new driver and although the driver doesn't perform any instrument communication, we still need to provide an address. Select "Other" under "Interface" and type any string in the "Address" text box. This is to ensure that every instrument has a unique address so that the *Instrument Server* can access multiple instances of the *Simple Signal Generator*, if needed. Finally, click "OK" to close the dialog.




Fig. 12.3. The example instrument driver *Simple Signal Generator*.

The new instrument should appear in the main *Instrument Server* list. Double-click the instrument name will bring up the driver configuration window, as shown in Fig. 12.3. To test the code, start the driver with the "Start" button and make sure the "Trace"-checkbox is checked to view the sinusoid. To control parameters while the driver is running, either just update one of the controls, or go to the server window, expand

the “Simple Signal Generator”-item in the instrument list and use the “Set Value”-button to set a new value. If the “Update continuously”-control in the driver dialog is checked, you will see the trace change in real time as parameters are modified.

12.4 Subclassing the VISA driver

The previous example was a little bit unusual, since no actual instrument communication was performed in the driver. A more common situation would be where most communication can be handled with simple text-based commands as defined in the driver *INI* configuration file, but where a few advanced quantities need special functionality. For these cases, the easiest way to proceed is to subclass the `VISA_Driver` and only re-implement code for the special cases. The code below shows an example of what such a driver would look like:

```
from VISA_Driver import VISA_Driver

class Driver(VISA_Driver):
    """ This class re-implements the VISA driver """

    def performOpen(self, options={}):
        """Perform the operation of opening the instrument connection"""
        # calling the generic VISA open to make sure we have a connection
        VISA_Driver.performOpen(self, options=options)
        # do additional initialization code here...
        pass

    def performClose(self, bError=False, options={}):
        """Perform the close instrument connection operation"""
        # calling the generic VISA class to close communication
        VISA_Driver.performClose(self, bError, options=options)
        # do additional cleaning up code here...
        pass

    def performSetValue(self, quant, value, sweepRate=0.0, options={}):
        """Perform the Set Value instrument operation. This function should
        return the actual value set by the instrument"""
        # check quantity name
        if quant.name == 'Some_Special_Operation':
            # special case, perform special code to set value
            pass
        else:
            # otherwise, call standard VISA case
            value = VISA_Driver.performSetValue(self, quant, value, sweepRate,
options)
        return value

    def performGetValue(self, quant, options={}):
        """Perform the Get Value instrument operation"""
        # check quantity name
```

```

if quant.name == 'Some_Special_Operation':
    # special case, perform special code to get value
    value = 0.0
else:
    # for all other cases, call generic VISA driver
    value = VISA_Driver.performGetValue(self, quant, options)
return value

```

Note that both the `performOpen` and `performClose` functions have to call the generic VISA class, to make sure that the communication is properly initiated.

12.4.1 Helper functions for drivers subclassing the VISA driver

In addition to the helper functions provided by the generic driver object (described in Section `driverObj` above), the `VISA_Driver` provides the following helper function:

`writeAndLog(sCmd, bCheckError=True):`

The function will send the command string `sCmd` to the instrument.
If `bCheckError` is `True`, an error check is performed after the command has been sent.

`write(sCmd, bCheckError=True):`

Same as above, but no entry will be created in the *Instrument Log*, regardless of the log level. See Section `logs` for more information about the logging feature.

`write_raw(data):`

Write raw data bytes to the instrument, without interpreting termination characters, etc. No logging is performed.

`reply = askAndLog(sCmd, bCheckError=True):`

The function will send the command string `sCmd` to the instrument and wait for a reply. The reply is returned as a text string. If `bCheckError` is `True`, an error check is performed after the command has been sent and data has been received.

`reply = ask(sCmd, bCheckError=True):`

Same as above, but no entry will be created in the *Instrument Log*, regardless of the log level. See Section `logs` for more information about the logging feature.

reply = read(n_bytes=None, ignore_termination=False):

Read a total of `n_bytes` from the device, ignoring any termination characters.

If `n_bytes` is `None`, the complete buffer is read. If `ignore_termination` is set to `True`, the program will not check for or remove termination characters.

queryErrors():

Check for instrument errors by checking the event status register (`*ESR?`). An exception is raised if the instrument reports an error. The check only takes place if the item `query_instr_errors` in the VISA settings of the *INI* file is set to `True`.

12.5 Support for sweeping

Some quantities, for example the B-field of a magnet, require the output to be changed with a well-defined sweep rate whenever the value is updated. *Labber* provides support for swept quantities, but such drivers require a few extra configuration settings compared to standard drivers. The extra settings are described in the subsection below. For more information on how swept experiments are implemented in the *Measurement Setup* dialog, see Section [SweepModeSetup](#).

12.5.1 1 Sweeping - Driver definition file

In addition to the properties listed in Section [Quantities](#), the driver *INI*-file of an instrument that supports sweeping needs to define the following properties for a sweepable quantity:

sweep_cmd:

Command used to sweep data. Use "<sr>" for sweep rate or "<st>" for sweep time, and "<*>" for the value. Note that sweep rate will be defined in terms of change per second or change per minute, as set by the `sweep_minute`-setting defined below. If the instrument does not have a built-in command for sweeping, a similar effect can be achieved by repeatedly using the `set_cmd` to incrementally change the instrument value. To enable this feature, set `sweep_cmd` to `****REPEAT SET****`, followed by the time interval between setting values (in seconds). If no time interval is defined, default is 0.1 seconds.

sweep_check_cmd:

Command used to check if the instrument is currently in sweep mode. The instrument should return `True` or `1` if the instrument is sweeping towards a value. If `sweep_check_cmd` is not defined, the program will determine if an instrument is in sweep mode by continuously reading the current value and comparing it against the target value with resolution `sweep_res`, as defined below.

`sweep_res:`

Attainable resolution when sweeping an instrument, in absolute units. Default value is 10^{-10} , to avoid float rounding errors. This parameter is not used if the `sweep_check_cmd` is defined.

`stop_cmd:`

Command used to stop a sweep.

`sweep_rate:`

Default sweep rate, in rate per second or rate per minute (as set by the `sweep_minute` parameter defined below). If this value is non-zero, sweeping will be turned on automatically for this quantity. Default value is `0`.

`sweep_minute:`

If `True`, sweep rates are defined in terms of value rate per minute, otherwise in rate per second. Default is `False` (rate per second).

`sweep_rate_low:`

Minimal sweep rate. Default is `0`.

`sweep_rate_high:`

Maximal sweep rate. Default is `+Inf`

Note that the existence of the `sweep_cmd`-parameter defines whether a quantity is sweepable or not. If a quantity is sweepable, the *Instrument Server*, *Instrument Driver* and the *Measurement Setup* configuration dialogs will contain a few extra options for controlling the sweep rates. If the `sweep_cmd`-parameter is defined but the `set_cmd`-parameter is not, the driver will not allow direct setting of output values. This is useful

for instruments like magnets, whose output currents must always be swept at a certain rate.

12.5.2 Sweeping - Python code

If the sweeping functionality of a driver cannot be implemented using the built-in functionality based on the parameters in the driver definition file listed above, it is possible to write custom *Python* code for carrying out the sweeping. To begin with, the `performSetValue`-function for setting an instrument value (described in

Section [PythonCode](#) and Section [SubClassVISA](#) above) needs to be implemented to support sweeping:

```
def performSetValue(self, quant, value, sweepRate=0.0, options={}):
```

When re-implementing the `performSetValue`-function for a swept quantity, it is important that the code inspects the `sweepRate` parameter to see if the user wants to set the value directly (`sweepRate=0.0`), or perform sweeping (`sweepRate>0.0`). Note that in sweep mode (`sweepRate>0.0`), the function should not wait for the sweep to finish, since the sweep checking/waiting is handled by the *Instrument Server*. The `sweepRate` parameter is defined in terms of change per second or change per minute, as set by the `sweep_minute` configuration parameter defined in the section above.

In addition to the four standard functions `performOpen`, `performClose`, `performSetValue` and `performGetValue` described in Section [PythonCode](#), drivers that support sweeping may also re-implement the following functions:

```
def checkIfSweeping(self, quant, options={}):
```

The function should return `True` if the instrument is currently sweeping to the target value. The standard implementation will either send the `sweep_check_cmd` to the instrument or continuously read the current value and compare to the target, as described in Section [SweepDriver](#) above.

```
def performStopSweep(self, quant, options={}):
```

This function should stop the current sweep. The default implementation will send the `stop_cmd` to the instrument, as described in Section [SweepDriver](#) above.

12.6 Controller drivers

To make a controller driver, start by setting the key `controller` in the driver configuration file to `True`. This will make *Labber* automatically add a few quantities such as *Period* and *Input/output* signals for handling the controller operation. Note that these quantities will be added automatically, and shall not be included in the driver `.ini` file.

For the controller to operate properly, the driver `.py` file must implement the `performGetValue` for the *Output value*-quantity. The function should typically read the value of the *Input value* control, and then apply the proper control logic to generate the output value. For an example of controller driver, see the *PID Controller* driver provided with *Labber*.

12.7 Hardware arming and triggering

In hardware trigger mode, log instrument will be armed to wait for a hardware trigger before starting to acquire data. The function `isHardwareTrig(options)` can be used by both instruments outputting and instruments reading values to check if the measurement is in hardware trig mode. Instruments that supports hardware arming need to define the `support_arm` parameter in the *General settings* of the driver definition file (see Section [DriverINIGeneral](#) above), and implement the following function:

```
def performArm(quant_names, options={}):
```

The function should arm the instrument, to make it ready to acquire values for the list of quantities defined by `quant_names`.

The function `performArm` is called before issuing the trigger starting the measurement. See Section [MeasProg](#) for more information on how hardware triggering is configured in the *Measurement* program.

12.8 Hardware looping

Some instruments can perform looping of values within the instrument hardware. This allows for implementing more efficient looping than with a computer, since there will be no need for the computer to send new values to the instrument at each step value. For a

more detailed description of how hardware looping works and how it is configured in the *Measurement* program, see Section [MeasProg](#).

Hardware looping requires that both the instrument outputting and the instrument reading values support hardware looping, and that the instrument reading values supports hardware arming, as defined by the `support_hardware_loop` and `support_arm` parameters in the *General settings* of the driver definition file (see Section [DriverINIGeneral](#) above).

12.8.1 Hardware looping - outputting values

In hardware looping mode, the `performSetValue` function will be called n times for a step sequence containing n points. At the final call, the instrument should be configured to start outputting values when a trigger is issued. The function `isHardwareLoop(options)` can be used to check if the measurement is in hardware loop mode, and the function `(index, n_pts) = getHardwareLoopIndex(options)` can be used to get the current hardware loop index and the total number of point `n_pts`.

12.8.2 Hardware looping - reading values

In addition to defining the `support_hardware_loop` and the `support_arm` parameters in the driver definition file, the driver *Python* file needs to implement the `performArm` function for arming the instrument to acquire multiple values. The number of values to expect is given by the output `n_pts` of the function `(index, n_pts) = getHardwareLoopIndex(options)`. After the instrument has been armed and a trigger has been sent, the function `performGetValue` will be called multiple times to acquire the results.

In the same way as for instruments outputting values, the functions `isHardwareLoop(options)` and `(index, n_pts) = getHardwareLoopIndex(options)` can be used to check if the measurement is in hardware loop mode, and to get the current hardware loop index and the total number of point `n_pts`, respectively.

12.9 Python distribution

When starting an instrument driver, *Labber* will launch a dedicated driver process and execute its Python code in the new process. By default, *Labber* will use the default, built-in Python distribution, which currently is a 64-bit version of Python 3.7. However, this may change to a newer Python version in a future release of *Labber*.

12.9.1 Python distribution, 32-bit version

For compatibility reasons, *Labber* is also shipped with a 32-bit version of the same Python distribution (Windows only), to allow control of older instruments for which only 32-bit Windows DLLs/drivers are available.

To activate the 32-bit Python version for a specific driver, open the driver's configuration window in the *Instrument Server*, go to the "Communication"-section, click "Show advanced interface settings", and check the "Run in 32-bit mode"-box prior to starting the instrument driver. Note that each instance of an instrument driver is running in its own, separate process, which makes it is possible to have some drivers run in 32-bit mode, while others are running in 64-bit mode. To set an instrument's default setting to run in 32-bit mode, use the `use_32bit_mode`-flag described in Section [DriverINIGeneral](#) above.

12.9.2 Custom external Python distribution

Sometimes it is convenient to use an external Python distribution instead of *Labber's* built-in one. For example, a custom driver may be relying on a number of external Python packages, and it is convenient to install/update those packages using a Python package manager instead of manually copying them into the folder location of the *Labber* driver.

When setting up a custom external Python distribution, it is recommended to use Python 3.7, as well as using a python environment manager, such as *venv* or *Anaconda*. The drives depend on certain pip packages installed in this environment. A list of packages can be found in the installation directory at *python-labber/requirements.txt*. To install these packages, run the following from your python environment.

```
pip install -U -r requirements.txt
```

Once your environment is set up, you can configure *Labber* to use it for your custom drivers. Open *Labber's Preferences* dialog, go to the "Advanced"-section and point the "Python distribution" control to the file representing the `python` executable for a given Python environment. For *Anaconda/miniconda* distributions, the `python` executable is located directly in the root of each environment folder ("`pythonw.exe`", Windows). Note that the external Python distribution will be used for only instrument drivers located in the "Local Drivers" preference.

For more information about the *Anaconda/miniconda* Python distributions, see <https://www.continuum.io/>.


12.9.3 Troubleshooting, external Python distribution

If a driver process terminates immediately upon starting or if a dialog pops up with a “*Broken pipe*”-message, there are most likely missing packages in the external Python distribution. To find out which package is missing, quit the *Instrument Server* and restart it from a terminal window to get access to the standard error output.

- For Windows, open a terminal window, `cd` to the location of the *Instrument Server* application, then run the application `InstrumentServer-Console.exe` from the command line. In addition, the “*Python distribution*” variable mentioned in Section [PythonDistExternal](#) above should point to the file `"python.exe"` instead of `"pythonw.exe"`.
- On Windows, there have been reports of incompatibilities with certain versions of Anaconda. Anaconda3 v. 4.4.0 is the most recent working version to have been tested to work.

Labber

APPENDIX A: PYTHON API



Appendix A: Python API

The *Labber* Python API (application program interface) provides Python classes and functions for controlling instruments in the *Labber* Instrument Server, for reading and writing *Labber* log files, and for scripting *Labber* measurements.




Fig A1. Overview and structure of the components in the *Labber* software package, including the Python API

A.1 Installation

The API is included when installing *Labber*, and the files are located in the "Script" folder of the main program directory.

A.1.1 Installing the API with pip

The *Labber* API can be installed into a local Python environment by using `pip`. It is recommended to use the *Anaconda/miniconda* package manager for configuring the environment. To Install the *Labber* API, activate this environment, navigate into the Installation directory, navigate down one level into the "Script" directory, and execute `pip install .`. This will use the `setup.py` file to install the API.

The Windows installer will install the *Labber* API by default if it detects a Python distribution on the system path. Whether or not a Python distribution is detected, the installer will prompt with the option to select a specific Python environment into which to install the API. To do this, specify the top-level directory of the Python environment

(the one which includes `python.exe`). This will ensure the API is updated along with the Labber programs.

A.1.2 Requirements

The Labber API requires the following Python packages, and will install them automatically during setup.

- Python 3.6 to 3.9
- NumPy
- h5py
- PyQt5
- qtpy
- msgpack
- sip
- future

A.1.3 Testing the API

To test the Labber API installation, execute the following code in a Python console:

```
import Labber
print (Labber.version)
```

A.1.4 Upgrading from earlier versions

In Labber versions earlier than 1.8, the Python API was not installed via pip, but relied on adding the "Script" folder to the PYTHONPATH, either as an environment variable or with a `.pth` file. This addition to the PYTHONPATH should be removed when updating to Labber 1.8 or later. Since the Labber install location has changed, this will no longer point to the correct location

Note that in Labber version 1.1 and later, the *ScriptTools* module has been moved into the `Labber` module. To make scripts written for older versions of Labber work with version 1.1 and later, replace `import ScriptTools` with `from Labber import ScriptTools`.

A.2 Instrument server

The instrument server API provides functions that allow Python to communicate with the Labber Instrument server. The API provides functionality both for communicate and control instruments as well as scheduling measurements using the Instrument server's built-in scheduler.

A.2.1 Labber client

The instrument control API uses a client object to communicate with the Instrument server. To initialize the client, use the `connectToServer` function. The following example will connect to an instrument server on the local computer and list all available instruments.

```
import Labber
# connect to server
client = Labber.connectToServer('localhost')
# get list of instruments
instruments = client.getListOfInstrumentsString()
for instr in instruments:
    print(instr)
# close connection
client.close()
```

A.2.2 Scheduling measurements

The function `schedule_measurement` is used to schedule a measurement using the Instrument server's queueing system. For measurement that are scheduled to run immediately, the function will wait until the measurement has finished, and return the full path of the final measurement file. For measurement that are scheduled for the future, the function will return directly without waiting.

Note that the server timeout should be set to `None`, to allow the client to wait for a long time for the measurement to finish before timing out. The following example illustrates the procedure.

The queue can also be accessed and modified directly from the API. The function `get_scheduler_queue()` returns an ordered list of scheduled measurements. The function `prioritize_measurement(name)` can be used to move a measurement to the top of the queue.

```
import Labber
# connect to server
client = Labber.connectToServer('localhost', timeout=None)
# schedule experiment, wait to finish
```

```
output_file = client.schedule_measurement(name= 'Test', path='~/Desktop/Test.hdf5')
print('Final output file:', output_file)
# close connection
client.close()
```

A.2.3 Connecting to instruments

To set or read values from an individual instrument, first use the function `connectToInstrument` of the Labber client object to access the instrument. The function will return an instance of a `InstrumentClient` object, from which the entire configuration or individual quantities of the instrument can be set or read.

The following example will first connect to a Labber Instrument server, and then connect to a dc voltage source and a voltmeter using a GPIB interface. Next, it will output voltages with the dc source, and measure the corresponding response with the voltmeter.

```
import Labber
import time, numpy as np
# connect to server
client = Labber.connectToServer('localhost')
# connect to specific instruments
volt_source = client.connectToInstrument('Yokogawa 7651 DC Source',
                                         dict(interface='GPIB', address='6'))
volt_meter = client.connectToInstrument('Agilent 34401 Multimeter',
                                         dict(interface='GPIB', address='1'))

# start drivers
volt_meter.startInstrument()
volt_source.startInstrument()

# put zero voltage to source and turn on the output
volt_source.setValue('Voltage', 0.0)
volt_source.setValue('Output', True)

# do a loop from zero to one volt
for volt in np.linspace(0.0, 1.0, 11):
    # set value to voltage source and wait
    volt_source.setValue('Voltage', volt)
    time.sleep(0.2)
    # read value from voltmeter
    measVolt = volt_meter.getValue('Voltage')
    # print result
    print('Set value: %.2f V, measured value: %.2f V' % (volt, measVolt))

# close client
client.close()
```

See the reference documentation below for more information on functions available in the `InstrumentClient` class.

A.2.4 Blocking vs. non-blocking clients

Labber supports two types of clients. The standard client type discussed so far is a *blocking* client, which will block program execution while waiting for a response from the server. Blocking clients are typically used in scripts where a number of instrument values are set or read in a pre-defined sequential order.

Function calls to a *non-blocking* client, on the other hand, will not wait for a response from the instrument server. Instead, the client uses a system of callback functions to notify the program that a new instrument value has been set or read. The advantage of a non-blocking client is that the main program thread will not be blocked while waiting for a result, and that multiple instrument operations can be performed in parallel. Non-blocking clients are typically used in user-interface driven applications, where the dialogs and user-interface elements need to remain responsive.

To create a non-blocking client, use the same `connectToServer`-function described above but with the argument `wait_for_reply` set to `False`. The following script will perform the same procedure of setting and reading voltages as the blocking-client script shown above, but it will do it in a non-blocking framework. Instead of directly returning values, the `connectToServer`, `connectToInstrument`, `setValue` and `getValue` functions will use *callback* functions to retrieve values.

Note that the event handling of the non-blocking client and the callback functionality is handled by the Qt framework.

```
import Labber
from qtpy.QtCore import QApplication
import functools
import numpy as np

class TestNonBlockClient():
    def __init__(self):
        # keep track of instrument references
        self.dInstr = dict()

    def connect(self):
        # connect to server, call 'connectionOpen' upon completion
        self.client = Labber.connectToServer('localhost', wait_for_reply=False,
                                             callback_open=self.connectionOpen,
                                             callback_network_error=self.on_error,
                                             callback_instrument_error=self.on_error)
```

```

def connectionOpen(self, data):
    # callback after connection has been established
    print('Connection open')
    self.nStart = 0
    # connect to instruments, call 'connected' upon completion
    newCallback = functools.partial(self.connected, 'Agilent 34401')
    self.client.connectToInstrument('Agilent 34401 Multimeter',
                                    dict(address='1', interface='GPIB'),
                                    callback=newCallback)
    newCallback = functools.partial(self.connected, 'Yokogawa 7651')
    self.client.connectToInstrument('Yokogawa 7651 DC Source',
                                    dict(address='3', interface='GPIB'),
                                    callback=newCallback)

def on_error(self, message):
    # print error
    print('Error: %s.\n\n' % message)

def connected(self, sHardware, instr):
    # keep track of instruments
    self.dInstr[sHardware] = instr
    # start driver, call 'started' upon completion
    newCallback = functools.partial(self.started, sHardware)
    instr.startInstrument(callback=newCallback)

def started(self, sHardware, data):
    # check that both instruments have been started
    self.nStart += 1
    if self.nStart==2:
        self.yoko = self.dInstr['Yokogawa 7651']
        self.volt = self.dInstr['Agilent 34401']
        # set yoko loop
        self.lLoop = np.linspace(0.0, 1.0, 11)
        self.nIndex = -1
        # start loop
        self.loop()

def loop(self):
    # check if looping is done
    self.nIndex += 1
    if self.nIndex < len(self.lLoop):
        # keep looping
        self.yoko.setValue('Value', self.lLoop[self.nIndex],
                           callback=self.stepDone)
    else:
        print('Finished!')

def stepDone(self, data):
    # voltage has been set, read response
    self.yoko.getValue('Value', callback=self.logDone)

def logDone(self, measVolt):
    # print result
    volt = self.lLoop[self.nIndex]
    print('Set value: %.2f V, measured: %.2f V' % (volt, measVolt))

```

```

        # keep looping
        self.loop()

if __name__ == '__main__':
    # start Qt event loop
    app = QApplication([])
    # create test object and connect to server
    test = TestNonBlockClient()
    test.connect()

```

A.2.5 Function definitions

Labber.connectToServer(*address='localhost', wait_for_reply=True, port=None, timeout=10, callback_open=None, callback_network_error=None, callback_instrument_error=None, binary_transfer_for_mat=None*)

Connect to Labber Instrument server and return a Labber client object.

There are two version of Labber clients, blocking and non-blocking ones. Blocking clients will wait for the instrument server to send a reply before returning, whereas non-blocking client will return immediately and call callback functions once the values are available.

Parameters

- **address** (*str, optional*) – IP address of Labber Instrument server. Default is localhost.
- **wait_for_reply** (*bool, optional*) – If True, the function will return a blocking client. Default is True.
- **port** (*int, optional*) – Port number for server communication. Default is 9406.
- **timeout** (*int or float, optional*) – Longest time to wait for the server to reply. Default is 10 seconds.
- **callback_open** (*function, optional*) – Callback function called after communication has been established. The function should have a single boolean argument, which will state if the connection was successful or not. Only relevant if *wait_for_reply* is False, ie for non-blocking clients.
- **callback_network_error** (*function, optional*) – Callback function called in case of network error. The function should take a single argument that will contain the error message. Only relevant if *wait_for_reply* is False, ie for non-blocking clients.

- **callback_instrument_error** (*function, optional*) – Callback function called in case of instrument error. The function should take a single argument that will contain the error message. Only relevant if *wait_for_reply* is False, ie for non-blocking clients.
- **binary_transfer_format** (*bool, optional*) – If True, data between the client and the server is sent as binary data instead of text. The value must match the Data transfer format in the Labber Instrument server preferences. Default is True.

Returns

client – Labber client object, either blocking or non-blocking version.

Return type

Client object

Examples

Open connection to server and list connected instruments.

```
>>> import Labber
>>> client = Labber.connectToServer('localhost')
>>> instruments = client.getListOfInstrumentsString()
>>> print(instruments)
>>> client.close()
```

A.2.6 Class definitions

A.2.6.1 Blocking client

The client object should not be initialized directly. Instead, use the `connectToServer()`-function defined above.

```
class Labber.LabberBlockingClient(sAddress='localhost', port=9406, timeout=10, binary_transfer_format=None, convert_to_unicode_if_py2=True)
```

Bases: `object`

Labber client, blocking execution while waiting for server response.

Parameters

- **sAddress** (*str, optional*) – IP address of Labber Instrument server. Default is localhost.
- **port** (*int, optional*) – Port number for server communication. Default is 9406.

- **timeout** (*int or float, optional*) – Longest time to wait for the server to reply. Default is 10 seconds.
- **binary_transfer_format** (*bool, optional*) – If True, data between the client and the server is sent as binary data instead of text. Default is True.

close()

Close the connection to the server.

connectToInstrument(sHardware, dComCfg, bCreateNew=False)

Connect to an instrument object on the instrument server.

Parameters

- **sHardware** (*str*) – Name of instrument hardware to connect to.
- **dComCfg** (*dict*) –

Dictionary describing the communication address of the instrument. Either the *name* key or the *interface* + *address* keys must be defined. The dictionary is defined by the following keys:

Namestr

Name of instrument.

interface{'GPIB', 'TCPIP', 'USB', 'Serial', 'VISA', 'Other', 'None'}

Communication interface.

addressstr

Instrument address string

startup{'Set config', 'Get config', 'Do nothing'}

Operation to perform at instrument startup.

lockbool

If True, instrument will be locked while in use.

- **bCreateNew** (*bool, optional*) – If True, a new instrument will be created if the requested one is not already present. Default is False.

Returns

instr – Object representing an instrument on the Labber instrument server.

Return type

InstrumentClient object

createInstrument(*sHardware*, *dComCfg*)

Create an instrument on the instrument server.

Parameters

- **sHardware** (*str*) – Name of instrument hardware to connect to.
- **dComCfg** (*dict*) –

Dictionary describing the communication address of the instrument.
The dictionary is defined by the following keys:

namestr

Name of instrument.

interface{'GPIB', 'TCPIP', 'USB', 'Serial', 'VISA', 'Other', 'None'}

Communication interface.

addressstr

Instrument address string.

startup{'Set config', 'Get config', 'Do nothing'}

Operation to perform at instrument startup.

lockbool

If True, instrument will be locked while in use.

getListOfInstruments()

Get a list of instruments present on the Labber instrument server.

Returns

instruments – List of instruments on the server. Each element of the list is a two-element tuple (name, comcfg), where *name* is the hardware name and *comcfg* is a dict with communication settings.

Return type

list of tuple

getListOfInstrumentsString()

Get a list of instruments present on the Labber instrument server.

Returns

instruments – List of strings describing instruments on the server.

Return type

list of str

schedule_measurement(*path_to_configuration*, *output_path=None*, *priority=False*, *schedule_d=None*, *period=None*, *name=None*, *command_args=[]*)

Schedule measurement using the instrument server queueing system.

For measurement that are scheduled to run immediately, the function will wait until the measurement has finished, and return the full path of the final measurement file. For measurement that are scheduled for the future, the function will return None directly without waiting.

Parameters

- **path_to_configuration** (*str*) – Path to Labber measurement configuration to run, saved in either .labber, .json or .hdf5 format.
- **output_path** (*str, optional*) – Path for output measurement file. Default is None, in which case the resulting output file is put in the Labber database.
- **priority** (*bool, optional*) – Priority in scheduling system. Default is False.
- **scheduled** (*float, optional*) – Scheduled time for measurement to run, in number of seconds passed since epoch. Default is None, which schedules immediately.
- **period** (*float, optional*) – Periodicity of measurement, measured in seconds. Default is None, in which case the measurement will only run once.
- **name**(*str*) - Optional name for the measurement. Required to reprioritize the measurement from the API
- **command_args** (*list, optional*) – Command-line arguments to pass on to the Measurement engine. Only used if scheduled and period are None.

Returns

Path of output file if scheduled is None and period is None, else None.

Return type

Str

get_scheduler_queue()

Return a list of scheduled measurements in the queue. Each measurement is described by a dictionary including name, path, priority, periodicity, and time scheduled.

Returns

List of dictionaries

Return type

List

`clear_queue()`

Remove all scheduled measurements from the queue

Returns

None

`prioritize_measurement(name)`

Add the priority designation to a measurement. If multiple measurements have the same name, priority will be added to all of them

Parameters

- **name(str)** - Name of the measurement to be prioritized

Returns

None

A.2.6.2 Non-blocking client

The client object should not be initialized directly. Instead, use the `connectToServer()`-function defined above.

```
class Labber.LabberClient(callbackNetworkError, callbackInstrError, sAddress='localhost', port=None, timeout=None, callbackOpen=None, callbackMessage=None, parent=None, convert_to_unicode_if_py2=True, binary_transfer_format=None)
```

Bases: `PyQt5.QtCore.QObject`

Labber client, non-blocking version.

Parameters

- **callbackNetworkError** (*function*) – Callback function called in case of network error. The function should take a single argument that will contain the error message.
- **callbackInstrError** (*function*) – Callback function called in case of instrument error. The function should take a single argument that will contain the error message.
- **sAddress** (*str, optional*) – IP address of Labber Instrument server. Default is localhost.
- **port** (*int, optional*) – Port number for server communication. Default is 9406.
- **timeout** (*int or float, optional*) – Longest time to wait for the server to reply. Default is 10 seconds.
- **callbackOpen** (*function, optional*) – Callback function called after communication has been established. The function should have a single boolean argument, which will state if the connection was successful or not.
- **callbackMessage** (*function, optional*) – Callback function for status updates from the server. The function should have a single string argument with the status.
- **binary_transfer_format** (*bool, optional*) – If True, data between the client and the server is sent as binary data instead of text. Default is True.

close(*bForce=False*)

Close the connection to the server.

Parameters

bForce (*bool, optional*) – If True, the connection is shut down without waiting for it to close. Default is False.

connectToInstrument(*sHardware, dComCfg, callback, bCreateNew=False*)

Connect to an instrument object on the instrument server.

Parameters

- **sHardware** (*str*) – Name of instrument hardware to connect to.
- **dComCfg** (*dict*) –

Dictionary describing the communication address of the instrument. Either the *name* key or the *interface* + *address* keys must be defined. The dictionary is defined by the following keys:

namestr

Name of instrument.

interface{'GPIB', 'TCPIP', 'USB', 'Serial', 'VISA', 'Other', 'None'}

Communication interface.

addressstr

Instrument address string

startup{'Set config', 'Get config', 'Do nothing'}

Operation to perform at instrument startup.

lockbool

If True, instrument will be locked while in use.

- **callback** (*function*) – Callback function called after the instruments has been created. The first argument will be an InstrumentClient object representing the instrument on the Labber instrument server.
- **bCreateNew** (*bool, optional*) – If True, a new instrument will be created if the requested one is not already present. Default is False.

createInstrument(*sHardware, dComCfg, callback=None*)

Create an instrument on the instrument server.

Parameters

- **sHardware** (*str*) – Name of instrument hardware to connect to.
- **dComCfg** (*dict*) –

Dictionary describing the communication address of the instrument. The dictionary is defined by the following keys:

namestr

Name of instrument.

interface{'GPIB', 'TCPIP', 'USB', 'Serial', 'VISA', 'Other', 'None'}

Communication interface.

addressstr

Instrument address string

startup{'Set config', 'Get config', 'Do nothing'}

Operation to perform at instrument startup.

lockbool

If True, instrument will be locked while in use.

- **callback** (*function, optional*) – Callback function called after the instruments has been created. The first argument will be an InstrumentClient object representing the instrument on the Labber instrument server.

firstCallbackStatus(*callbackProgress, callbackCurrentValue, data*)

First callback occurring after the server sends back status updates. The function will handle errors and then call the next callback

getListOfInstruments(*callback*)

Get a list of instruments present on the Labber instrument server.

Parameters

callback (*function*) – Callback function called after the list of instruments has been retrieved. The first argument will be the list of instruments.

Returns

instruments – List of instruments on the server. Each element of the list is a two-element tuple (name, comcfg), where *name* is the hardware name and *comcfg* is a dict with communication settings.

Return type

list of tuple

getListOfInstrumentsString(*callback*)

Get a list of instruments present on the Labber instrument server.

Parameters

callback (*function*) – Callback function called after the list of instruments has been retrieved. The first argument will be the list of instruments.

Returns

instruments – List of strings describing instruments on the server.

Return type

list of str

A.2.6.3 *Instrument client*

The InstrumentClient represents an instrument on the server. Note that the InstrumentClient object should not be initialized directly, but rather created using the `connectToInstrument` or `createInstrument` functions of a *LabberClient* object.

```
class Labber.InstrumentClient(client, instrRef, IdQuant, dOption, block=True)
```

The InstrumentClient is representing an instrument on the server.

abortCurrentOperation(*callback=None*)

Abort current operation, but keep instrument running.

Parameters

callback (*function, optional*) – Callback function called after the instruments has been aborted. Only relevant for non-blocking clients.

arm(*quantities, callback=None, options={}*)

Arm instrument to prepare for later hardware-triggered data acquisition

Parameters

- **quantities** (*list of str*) – Name of quantities that will be acquired when the instrument is triggered.
- **callback** (*function, optional*) – Callback function called after the instrument has been armed. Only relevant for non-blocking clients.

disconnectFromInstr(*callback=None*)

Disconnect from instrument.

Parameters

callback (*function, optional*) – Callback function called after the instruments has been disconnected. Only relevant for non-blocking clients.

getInstrConfig(callback=None)

Get values from the driver.

Parameters

callback (*function, optional*) – Callback function called after the instrument config has been retrieved. Only relevant for non-blocking clients.

Returns

values – Dictionary with instrument values. The keys are names of instrument quantities. Note that only blocking clients will return a value.

Return type

dict

getLocalInitValuesDict()

Get instrument values as recorded at instrument initialization.

Returns

values – Dictionary with instrument values. The dict keys are names of the instrument quantities.

Return type

dict

getLocalOptionsDict()

Get instrument options as recorded at instrument initialization.

Returns

options – Dictionary representing instrument options. The dictionary is defined by the following keys:

modelstr

Instrument model number/name.

optionslist of str

List of strings describing installed options.

Return type

dict

getValue(*sQuant*, *callback=None*, *callbackProgress=None*, *callbackCurrentValue=None*, *options={}*)

Get value of the specified quantity

Parameters

- **sQuant** (*str*) – Name of quantity to set.
- **callback** (*function, optional*) – Callback function called after the instrument value has been retrieved. Only relevant for non-blocking clients.
- **callbackProgress** (*function, optional*) – Callback function for progress updates from the server. The function must take a single argument, which will be a float between 0.0 and 1.0 indicating progress. Only relevant for non-blocking clients.
- **callbackCurrentValue** (*function, optional*) – Callback function for value updates from the server. The function must take a single argument (current value), and is used to show the current value during slow operations, like averaging. Only relevant for non-blocking clients.

Returns

value – Value of the instrument. Note that only blocking clients will return a value.

Return type

float, bool or numpy array.

isRunning(*callback=None*)

Check if instrument driver is running.

Parameters

callback (*function, optional*) – Callback function called after the instruments has been checked. Only relevant for non-blocking clients.

Returns

isRunning – True if instrument is running. Note that only blocking clients will return a value.

Return type

bool

setInstrConfig(dValues={}, callback=None, always_update_all=True)

Send values to the driver.

Parameters

- **dValues** (*dict*) – Dictionary with new values. The keys are names of instrument quantities.
- **callback** (*function, optional*) – Callback function called after the instrument config has been set. Only relevant for non-blocking clients.
- **always_update_all** (*bool, optional*) – If True, the instrument settings are updated even if values have not changed compared to the local settings stored in the driver.

Returns

values – Dictionary with actual values. The keys are names of instrument quantities.

Return type

dict. Note that only blocking clients will return a value

setValue(sQuant, value, rate=0.0, wait_for_sweep=True, callback=None, callbackProgress=None, callbackCurrentValue=None, options={})

Set new value to the specified quantity

Parameters

- **sQuant** (*str*) – Name of quantity to set.
- **value** (*float, bool or numpy array*) – New value.
- **rate** (*float, optional*) – Sweep rate.
- **wait_for_sweep** (*bool, optional*) – If True and *rate* is non-zero, the instrument is waiting for a sweep to finish.
- **callback** (*function, optional*) – Callback function called after the instrument value has been set. Only relevant for non-blocking clients.

- **callbackProgress** (*function, optional*) – Callback function for progress updates from the server. The function must take a single argument, which will be a float between 0.0 and 1.0 indicating progress. Only relevant for non-blocking clients.
- **callbackCurrentValue** (*function, optional*) – Callback function for value updates from the server. The function must take a single argument (current value), and is used to show the current value during slow operations (sweeping). Only relevant for non-blocking clients.

Returns

value – Actual value of the instrument. Note that only blocking clients will return a value.

Return type

float, bool or numpy array.

startInstrument(dOption=None, callback=None)

Start the instrument.

Parameters

- **dOption** (*dict, optional*) –
Dictionary representing instrument options. The dictionary is defined by the following keys:

modelstr

Instrument model number/name.

optionslist of str

List of strings describing installed options.

- **callback** (*function, optional*) – Callback function called after the instruments has been started. Only relevant for non-blocking clients.

stopInstrument(bForceQuit=False, callback=None)

Stop the instrument.

Parameters

- **bForceQuit** (*bool, optional*) – If True, the instrument is shut down without waiting for it to close. Default is False.

- **callback** (*function, optional*) – Callback function called after the instruments has been stopped. Only relevant for non-blocking clients.

```
waitForSweep(sQuant, value=None, callback=None, options={}, callbackCurrentValue=None)
```

Wait for swept instrument to reach final point or certain value.

A.3 Log files

The `LogFile` class provides functionality for reading and writing data from Labber log files.

A.3.1 Reading data from Labber

Labber log files are accessed using the `LogFile` class, which contain a number of functions for reading and writing data (see [class definition](#) below). A `LogFile` object is created by passing the path to a Labber log file as the first argument.

A.3.1.1 Log information

A Labber log file contains both instrument settings and measured data, as well as metadata information from the database such as *User*, *Tags*, *Project* and *Comment*. The following example will print basic information about the log file `TestLog.hdf5`:

```
import Labber

f = Labber.LogFile('TestLog')
print('Number of entries:', f.getNumberOfEntries())

print('Step channels:')
step_channels = f.getStepChannels()
for channel in step_channels:
    print(channel['name'])

print('Log channels:')
log_channels = f.getLogChannels()
for channel in log_channels:
    print(channel['name'])

print('User:', f.getUser())
print('Tags:', f.getTags())
print('Project:', f.getProject())
print('Comment:', f.getComment())
```

The `LogFile` class also contains functions for setting log metadata, see the [class definition](#) below.

A.3.1.2 Log data

A Labber log file contains data from one or multiple *channels*. The data is organized into log *entries*, where each entry contains a one-dimensional vector of values for each channel. The entries correspond to the traces shown in the Labber *Log Viewer* program. The *LogFile* class provides a number of functions for accessing the data, as illustrated in the example below:

```
import Labber

f = Labber.LogFile('TestLog')

# get values of all channels for a specific entry (in this case first entry)
d = f.getEntry(0)
for (channel, value) in d.items():
    print(channel, ":", value)

# get entry as x,y data, let Labber determine which channels to read
(x,y) = f.getTraceXY()

# get data for all entries for a specific channel as a 2D numpy array
data = f.getData('Voltage')

# get last recorded value of a specific channel in the measurement config
# this function also works for channels that are not step items
value = f.getChannelValue('Integration time')

# get last recorded values of all channels
# useful for extracting all instrument settings
value = f.getChannelValuesAsDict()
```

For more information on the various class methods, see the [class definition](#) below.

A.3.2 Creating Labber log files

The API provides functionality for creating Labber log files that can be opened by the Labber *Log Browser* and *Log Viewer* programs. This makes it possible to add custom data to the Labber database, such as simulation results or data acquired outside of the Labber *Measurement* program.

The following lines of Python code will create a log file with sinusoid signals with different frequencies in the Labber database.

```
import Labber
import numpy as np

# create step data
vTime = np.linspace(0,1,501)
```

```

vFreq = np.linspace(1,10,100)
# define step channels
lStep = [dict(name='Time', unit='s', values=vTime),
         dict(name='Frequency', unit='Hz', values=vFreq)]
# define log channels
lLog = [dict(name='Signal', unit='V', vector=False)]

# create log file
f = Labber.createLogFile_ForData('TestSinusoid', lLog, lStep)

# add log entries
for freq in vFreq:
    data = {'Signal': np.sin(2*np.pi*freq*vTime) }
    f.addEntry(data)

```

Note that log files created with the function `createLogFile_ForData` can only be used with the *Log Browser* and *Log Viewer* programs. It is not possible to open or run such a file in the *Labber Measurement* program.

A.3.3 Function definitions

Labber.getTraceDict(value=[], x0=0.0, dx=1.0, x1=None, logX=False, x=None)

Create a dict with metadata for Labber (x,y) traces.

Parameters

- **value** (*list or np.array*) – Vector of y-values for trace data.
- **x0** (*float, optional*) – x-value for first data point in trace vector. Default is 0
- **dx** (*float, optional*) – Step size for x data. If specified, the x-vector starts at “x0”, and every subsequent point is spaced by “dx”. Default is 1
- **x1** (*float, optional*) – x-value for last data point in trace vector. If specified, the “dx” parameter is ignored, and the x-vector will be a linear ramp between “x0” and “x1”.
- **logX** (*bool, optional*) – If True, the values between x0 and x1 are interpolated logarithmically. Only valid if “x0” and “x1” are specified.
- **x** (*list or np.array, optional*) – Vector of x-value to match the y-values. The input must have the same number of elements as the “values” parameter. If specified, the values of “x0”, “dx”, “x1” and “logX” are ignored.

Returns

d – Python dict with values and metadata for describing a (x,y) trace.

Return type

dict

Labber.createLogFile_ForData(*name*, *log_channels*, *step_channels*=[], *use_database*=True)

Create a log file for custom data storage in the Log database.

Parameters

- **name** (*str*) – Name or path of log file.
- **log_channels** (*list of dict*) –

List of dict describing the log channels. The list corresponds to the log channels in the Measurement dialog. The dictionary is defined by the following keys:

namestr

Name of channel.

unitstr, optional

Unit of channel.

complexbool, optional

If True, the channel contains complex data. Default is False.

vectorbool, optional

If True, the channel contains vector data. Default is True.

x_namestr, optional

Label of the x-axis for vector data. Default is "Index".

x_unitstr, optional

Unit of x-values for vector data.

- **step_channels** (*list of dict, optional*) –

List of dict describing the step channels. The list corresponds to the Step sequence in the Measurement dialog. If *step_values* is left undefined, the resulting log file will contain a collection of traces without a uniform pre-defined dimensionality. The dictionary is defined by the following keys:

namestr

Name of channel.

values1D numpy array

Step values for step channels. The length of the vector defines the dimensionality of the data in the resulting log file.

unitstr, optional

Unit of channel.

combo_defs list of str, optional

Enumerator labels for quantity. If specified, Labber will define the channel to be of "COMBO" datatype, and the "values" data must be integer values between 0 and len(combo_defs) - 1.

- **use_database** (*bool, optional*) – If True, the log file is put in the central log database, otherwise the path set by the log name. Default is True.

Returns

log – LogFile object representing the newly created log.

Return type

LogFile object

Examples

Example 1: Create log without step values or fixed dimensions. Note that entries do not need to have the same length.

```
>>> import Labber
>>> lLog = [dict(name='x'), dict(name='y')]
>>> l = Labber.createLogFile_ForData('TestLog', lLog)
```

To add two entries to the log defined above:

```
>>> x = np.linspace(0,1,501)
>>> data = {'x': x, 'y': np.sin(2*np.pi*5*x) }
>>> l.addEntry(data)
>>> x = np.linspace(-1.2,1.2,201)
>>> data = {'x': x, 'y': x**2 }
>>> l.addEntry(data)
```

Example 2: Create log file using pre-defined step values. In this example, the data dimensions are defined by the step channels, and all entries need to have the same length as the first step channel. Note the the presence of `vector=False` for the *Signal* channel, which notifies that the entry size is defined by the first step channel.

```
>>> import Labber
>>> vTime = np.linspace(0,1,501)
>>> vFreq = np.linspace(1,10,10)
>>> chTime = dict(name='Time', unit='s', values=vTime)
>>> chFreq = dict(name='Frequency', unit='Hz', values=vFreq)
>>> chSig = dict(name='Signal', unit='V', vector=False)
```

```
>>> f = Labber.createLogFile_ForData('TestData', [chSig], [chTime, chFreq])
```

To add data to the log defined above:

```
>>> for freq in vFreq:
>>>     data = {'Signal': np.sin(2*np.pi*freq*vTime)}
>>>     f.addEntry(data)
```

Example 3: Create log file using pre-defined step values, but allow individual entries to have different lengths. Compared to *Example 2* above, we use the “getTraceDict” function to define the x-values for the vector-valued data.

```
>>> import Labber
>>> import numpy as np
>>> frequencies = np.linspace(1, 10, 10)
>>> channel_f = dict(name='Frequency', unit='Hz', values=frequencies)
>>> channel_y = dict(name='Signal', unit='V', x_name='Time', x_unit='s')
>>> f = Labber.createLogFile_ForData('TestData', [channel_y], [channel_f])
```

To add data to the log defined above:

```
>>> t = np.linspace(0, 1, 501)
>>> for freq in frequencies:
>>>     y = np.sin(2 * np.pi * freq * t)
>>>     trace_dict = Labber.getTraceDict(y, x0=t[0], x1=t[-1])
>>>     data = {'Signal': trace_dict}
>>>     f.addEntry(data)
```

A.3.4 LogFile class

```
class Labber.LogFile(file_name, instrument_units=False)
```

Bases: `object`

The class handles reading and writing data to and from Labber log files.

Parameters

- **file_name** (*str*) – Labber hdf5 file with log data.
- **instrument_units** (*bool, optional*) – If True, data from the log file is returned in instrument units instead of physical units. Default is False.

addEntry(data)

Add one entry to log file.

Parameters

data (*dict*) – Dictionary with data. The keys should match the channel names, and the values should be 1D numpy arrays or dicts with Labber (x,y) trace data created with the “getTraceDict”-function.

For scalar channels, the length of the array must match the size of the innermost step loop.

If the log contains channels with both scalar and vector data, the dict value for channels that contain vector data should be an iterable with numpy arrays or trace dicts.

getChannelValue(channel_name)

Get value of a channel at the end of the measurement.

Parameters

channel_name (*str*) – Name of channel for getting value.

Returns

value – Channel value as recorded after finishing the measurement.

Return type

float, string, or dict

getChannelValuesAsDict(include_all_quantities=False)

Get value of all channels at the end of the measurement.

Parameters

include_all_quantities (*bool*) – If False, only channels defined in the Measurement dialog are returned. Otherwise, all quantities of all instruments are included

Returns

channels – Dict with all channel values. The key is the channel name.

Return type

dict

getComment(log=-1)

Get comment from log file.

Parameters

log (*int, optional*) – Log number within the log file. Default is -1 (last log)

Returns

comment – String with comment

Return type

str

getData(name=None, entry=None, inner=None, log=-1)

Retrieve data from the log file and return it as a numpy array.

Parameters

- **name** (*str, optional*) – Name or index of the channel with data. If None, data for the first log channel will be returned.
- **entry** (*int or iterable, optional*) – Entry number within log to retrieve. If None, all elements will be returned.
- **inner** (*int or iterable, optional*) – Index of the inner-most loop values to retrieve. If None, all elements will be returned.
- **log** (*int, optional*) – Log number within the log file. Default is -1 (last log)

Returns

data – Depending on the input arguments, the output data will be a floating point number or a one- or two-dimensional numpy array.

Return type

float or np.array

getEntry(entry=-1)

Retrieve an entry from the log file and return a dict with values.

Parameters

entry (*int, optional*) – Entry number to retrieve, as shown in the Log Viewer. Default is -1, which will get the last trace in the file.

Returns

d – Dictionary with entry data. Keys are the channel names, the values are floats, numpy arrays or dicts with vector data. In addition, the dictionary contains a key “timestamp”, which contains a timestamp (seconds since epoch) for the entry.

Return type

dict

getFilePath(tags)

Get path of hdf5 log file.

Returns

path – Full path and name of hdf5 log file.

Return type

str

getLogChannels()

Get log channels in the log file.

Returns

log_channels – List of dicts representing log channels. The dictionaries contain the following keys:

namestr

Name of channel.

unitstr

Unit of channel.

complexbool

If True, the channel contains complex data.

vectorbool

If True, the channel contains vector data.

Return type

list of dict

Examples

Get list of log channels from the log file *TestLog*.

```
>>> import Labber
```

```
>>> l = Labber.LogFile('TestLog')
>>> print(l.getLogChannels())
[{'name': 'Signal', 'unit': 'V', 'complex': False, 'vector': False}]
```

getNumberOfEntries(name=None, log=None)

Get number of entries in the log file for the given channel.

Parameters

- **name** (*str, optional*) – Name of channel for data count. Default is first log channel.
- **log** (*int, optional*) – Log configuration number within the log file. Default is None, which will count all entries for all logs.

Returns

n – Number of entries.

Return type

int

getNumberOfLogs()

Get number of individual log configurations in the log file.

Returns

n – Number of log configurations.

Return type

int

getProject()

Get project name from log file.

Returns

project – String with project name.

Return type

str

getStepChannels()

Get step channels in the log file.

Returns

log_channels – List of dicts representing step channels. The dictionary contains the following keys:

namestr

Name of channel.

unitstr

Unit of channel.

values1D numpy array

Step values for step channels.

complexbool

If True, the channel contains complex data. Always False for step channels.

vectorbool

If True, the channel contains vector data. Always False for step channels.

Return type

list of dict

Examples

Get list of step channels from the log file *TestLog*

```
>>> import Labber
>>> l = Labber.LogFile('TestLog')
>>> print(l.getStepChannels())
[{'name': 'Time', 'unit': 's', 'complex': False, 'vector': False,
  'values': array([ 0.    ,  0.002, ..., 0.998, 1.    ])},
 {'name': 'Frequency', 'unit': 'Hz', 'complex': False, 'vector': False,
  'values': array([ 1.,  5., 10.])}]
```

getTags()

Get tag list from log file.

Returns

tags – List of strings with tags.

Return type

list of str

getTraceXY(*y_channel=None, x_channel=None, entry=-1*)

Retrieve a trace with (x,y) data from the log file .

Parameters

- **y_channel** (*str or int, optional*) – Name or log index of the channel with y-data. Default is first log channel.
- **x_channel** (*str or int, optional*) – Name or step index of the channel with x-data. Default is first step channel.
- **entry** (*int, optional*) – Entry number to retrieve, as shown in the Log Viewer. Default is -1, which will get the last trace in the file.

Returns

(x,y) – A tuple with x and y data as 1-d numpy arrays.

Return type

tuple

getUser()

Get user from log file.

Returns

name – String with user name.

Return type

str

setComment(*comment, log=-1, set_all=True*)

Set comment in log file.

Parameters

- **comment** (*str*) – String with comment.
- **log** (*int, optional*) – Log number within the log file. Default is -1 (last log).
- **set_all** (*bool, optional*) – Set comment of all log numbers within the log file. Default is True.

setProject(*project*)

Set project name in the log file.

Parameters

project (*str*) – String with project name.

setTags(*tags*)

Set list of tags in the log file.

Parameters

tags (*list of str*) – List of string with tags.

setUser(*name*)

Set user name in the log file.

Parameters

name (*str*) – String with user name.

A.4 Script tools

The helper functions in the *ScriptTools* module are designed for repeatedly performing Measurements that each contain one-dimensional sweeps, and where one or multiple parameters of the Measurement configurations are updated between each measurement.

A.4.1 Initialization

The *ScriptTools* functions call the *Measurement* executable for performing the measurements. Before the tools can be used, the path to the executable must be set using the function `setExePath()` defined below.

Note that in Labber version 1.1 and later, the *ScriptTools* module has been moved into the `Labber` module. To make scripts written for older versions of Labber work with version 1.1 and later, replace `import ScriptTools` with `from Labber import ScriptTools`.

A.4.2 Example

The *ScriptTools* functions are best explained by an example, which we'll take from the domain of superconducting qubits. For the purpose of this example, we can view the qubit as a slightly anharmonic oscillator. The qubit is read out by coupling it to a microwave resonator. This coupling is arranged in a way that any drifts or changes in the qubit frequency will cause a shift of the readout resonator frequency.

As an example, we run a mock calibration experiment where we first perform cavity spectroscopy and then use that resonance frequency as the readout LO for qubit spectroscopy. We can leverage the existing fitting features to accomplish this.

We begin by showing the code needed to run the experiments. This code is also available as a Jupyter notebook and can be made available upon request.

```
import Labber
import numpy as np

# # Create a measurement object for resonator spectroscopy
res_spec_meas_obj = Labber.ScriptTools.MeasurementObject(
    'Resonator Spectroscopy.labber',
    'Resonator Spectroscopy.hdf5'
)
# # Use addFit to include the log channel name as well as fit function type
```

```

res_spec_meas_obj.addFit('Resonator Spectroscopy - Signal Frequency',
                        'Lorentzian')

# # Perform measurement set perform_fit = True
res_spec_meas_obj.performMeasurement(return_data = True, use_scheduler = False,
                                    perform_fit = True
                                    );

# # Create a measurement object for qubit spectroscopy
qubit_spec_meas_obj = Labber.ScriptTools.MeasurementObject(
    'Qubit Spectroscopy.labber',
    'Qubit Spectroscopy.hdf5'
)

# # Update the readout resonator LO based on the fit
qubit_spec_meas_obj.updateValueFromFit('Readout LO - Frequency',
    res_spec_meas_obj.sCfgFileFits, 'f0'
)

# # Use addFit to include the log channel name as well as fit function type
qubit_spec_meas_obj.addFit('Qubit Spectroscopy - Signal Frequency',
    'Lorentzian'
)

# # Perform measurement set perform_fit = True
qubit_spec_meas_obj.performMeasurement(return_data = True, use_scheduler = False,
    perform_fit = True);

```

A.4.3 Function definitions

Labber.ScriptTools.setExePath(path)

Set path to the Measurement program, must be done before running scripts

Parameters

path (str) – Path to Measurement.exe program. On Windows, the path is typically 'C:\Program Files\LabberProgram'.

Labber.ScriptTools.load_scenario_as_dict(file_name)

Load Labber measurement scenario from binary .labber or .json file

Parameters

file_name (str) – Path to Labber measurement scenario file (.labber or .json).

Returns

d – Python dict describing measurement scenario.

Return type

dict

Labber.ScriptTools.save_scenario_as_binary(*config*, *file_name*)

Save Labber measurement scenario as binary .labber file

Parameters

- **config** (*dict*) – Python dict describing Labber measurement scenario.
- **file_name** (*str*) – Path to output file.

Labber.ScriptTools.save_scenario_as_json(*config*, *file_name*)

Save Labber measurement scenario as .json file

Parameters

- **config** (*dict*) – Python dict describing Labber measurement scenario.
- **file_name** (*str*) – Path to output file.

A.4.4 MeasurementObject class

class **Labber.ScriptTools.MeasurementObject**(*sCfgFileIn*, *sCfgFileOut*)

Bases: `object`

Class for updating measurement objects and running experiments

Parameters

- **sCfgFileIn** (*str*) – Path of template config file that defines the Measurement.
- **sCfgFileOut** (*str*) – Path to output file that will be created when running the Measurement. This should typically be different from the template file, since the dimensionality of the configuration may change as data is added.
- **sCfgFileFits** (*str*) – Path to fit parameters output file that will be created when running the Measurement. The autogenerated name is `sCfgFileOut + '_fits'`. This is only created if `perform_fit` is set to `True` in `performMeasurement`.

addFit(*log_channel*, *fit_function*, *guess_params*)

A list of dictionaries added to the `fitList` property of the `MeasurementObject` that stores settings for automated fitting.

Parameters

- **log_channel** (*str*) – log channel for the particular fit.
- **fit_function** (*str*) – Function type for fitting, not case sensitive. Options are: 'Lorentzian', 'Gaussian', 'Exponential', 'Exponential W/ Sinusoid' (not sensitive to spaces between words in the last case).
- **guess_params** (*dict, optional*) – Dictionary of keyword: value pairs for initial guess on parameters. If not provided than guess parameters will automatically be generated from data. Users do not need to provide guesses for all fit parameters for a given fit function. Supported keys (not case sensitive) are:
 - Lorentzian: 'f0', 'amplitude', 'width', 'offset'
 - Gaussian: 'mean', 'amplitude', 'sigma', 'offset'
 - Exponential: 'amplitude', 'tau', 'offset'
 - Exponential W/ Sinusoid: 'amplitude', 'tau', 'offset', 'detuning', 'phaseoffset'

Returns

None

Return type

None

performFit()

Performs a fit defined by the entries in fitList to the log file (sCfgFileOut) from the completion of the performMeasurement function. Will save all fit parameters in a new log file with the same name as the data log file but with an appended "_fits" (sCfgFileFits). The log file is created via the LogFile class method createLogFile_ForData with use_database = False meaning that the fit log file will be stored in the same database as the data or the user can specify a different path by using the measurement object property sCfgFileFits. The following parameters are saved in each "_fits" log file.

Parameters

- **Input Waveform** (*trace dict*) – Log channel data waveform used for fitting.
- **Fit Waveform** (*trace dict*) – Waveform based off of the scalar values from the fit as well as the fit function type. To be compared against the input data waveform for comparison.

- **Residuals** (*trace dict*) – A waveform that is the difference between the input data and the fit waveform for visual inspection on goodness of fit.
- **Scalar Values** (*step channels*) – Each scalar value is saved as a step channel in the “_fits” logfile. Supported returns are:
 - Lorentzian: 'f0', 'f0 std', 'amplitude', 'amplitude std', 'width', 'width std', 'offset', 'offset std', 'Q' (loaded quality factor defined as 'f0' / 'width')
 - Gaussian: 'mean', 'mean std', 'amplitude', 'amplitude std', 'sigma', 'sigma std', 'offset', 'offset std'
 - Exponential: 'amplitude', 'amplitude std', 'tau', 'tau std', 'offset', 'offset std'
 - Exponential W/ Sinusoid: 'amplitude', 'amplitude std', 'tau', 'tau std', 'offset', 'offset std', 'detuning', 'detuning std', 'phaseoffset', 'phaseoffset std'

Returns

None

Return type

None

performMeasurement(*return_data=True, use_scheduler=True, perform_fit=False*)

Perform measurement and return (x,y)-tuple.

The function will start the application Measurement.exe.

Parameters

- **return_data** (*bool, optional*) – If True, the function will return a tuple with (x,y) data (see below). If False, the function will return the actual path of the output data file. Default is True.
- **use_scheduler** (*bool, optional*) – If True, the measurement will be executed using Labber’s internal scheduler. If False, a separate instance of the Measurement program will be launched to run the measurement. Default is False.
- **perform_fit** (*bool, optional*) – If True, at the completion of the measurement this will call the performFit function to fit the newly acquired data per the user specified log channel and fit function in the fitList. Before executing performMeasurement the user must use addFit to include fit settings. After

the completion of the fit, values are stored in a separate log file that is the same name as the results for performMeasurement but with an appended "_fits". Default is False.

Returns

(**x,y**) – A tuple with x and y data as 1-d numpy arrays. The x-data is taken from the first step channel, the y-data is taken from the first log channel.

Return type

tuple

rearrangeLog(channel_name, *extra_arg)

Re-arrange a log with N entries of length M to a 2D log with dim (N, M)

The "channel_name" determines which data to use when defining the second dimension. It is also possible to rearrange into a multi- dimensional log by specifying multiple channels, but if so, lists of step values for each dimension need to be specified as well. For example, to rearrange a log with 6 entries into a multi-dimensional log with 3*2 entries, use: rearrangeLog("Channel 1", [1.0, 2.0, 3.0], "Channel 2", [1.0, 2.0])

Parameters

- **channel_name** (*str*) – Path to log file.
- **values** (*list of float, optional*) – Step value of channel_name. If not specified, the values will be taken from log file.

setPrimaryChannel(channel_name)

Specify the primary channel name.

Values of all other updated channels will be defined by look-up tables relative to the primary channel values.

Parameters

channel_name (*str*) – Name of primary channel.

setOutputFile(filename)

Set output file when performing the measurement

Parameters

filename (*str*) – Path to output file.

updateValue(channel_name, value, itemType='SINGLE')

Update a single value in the config file.

The values are kept track of internally until the Measurement.exe program is called.

Parameters

- **channel_name** (*str*) – Name of channel to update.
- **value** (*float*) – New value to set to channel.
- **itemType** (*str, optional*) – Step item parameter to set, must be one of { `single`, `start`, `stop`, `center`, `span`, `step`, `n_pts` }. Default is `single`.

updateValueFromFit(channel_name, fit_logfile, fit_parameter, function=None, index=-1)

Update a single value from a fit logfile in the config file of a measurement.

The values are kept track of internally until the Measurement.exe program is called.

Parameters

- **channel_name** (*str*) – Name of channel to update.
- **fit_logfile** (*str*) – Path of log file containing the scalar fit parameters.
- **fit_parameter** (*str*) – Name of the fit parameter from the fits log file to be used in updating the value.
- **function** (*callable, optional*) – User defined function with the input variable being the fit_parameter value.
- **Index** (*int, optional*) – Default is -1. Index for scalar fit values stored in step channels.

A.5 Configurations

The classes and the function in the *config* module allow Labber Measurement *scenarios* to be modified or created from scratch. In most cases, the recommended workflow is to create a template scenario in the Measurement program with all the instruments and signal connections used in the setup, then load the scenario into the API and use the functions *add_step* and *add_log* which channels to step over or log.

A.5.1 Example

A Labber measurement scenario is represented by an instance of the class `Scenario`, which is described in more detail in the [Scenario class](#) section below. The `Scenario` object contains lists of `Instrument`, `Channel`, `StepItems` and log channel objects, as well as a number of settings and other configuration parameters.

We illustrate the process of creating a scenario with an example. The goal is to create a measurement that will generate a sine waveform with the *Simple Signal Generator* driver, send it to the *Signal demodulation* driver over a signal connection, then perform the demodulation and run the experiment for a few different values of the signal and demodulation frequency.

A.5.1.1 Example - Full code

The code used to generate the example scenario is shown below:

```
from Labber import Scenario
import numpy as np

# create and add instruments
s = Scenario()
instr_signal = s.add_instrument('Simple Signal Generator', name='Sine')
instr_demod = s.add_instrument('Signal Demodulation', name='Demod')

# set a few instrument settings
instr_demod.values['Use phase reference signal'] = False
instr_demod.values['Length'] = 1.0

# add signal connections between channels
s.add_connection('Sine - Signal', 'Demod - Input data')

# add step items, values can be defined with np array or keywords
s.add_step('Sine - Frequency', np.linspace(0, 10, 51))
s.add_step('Demod - Modulation frequency', start=1, stop=9, step=4)

# add Log channels
s.add_log('Demod - Value')

# set metadata
s.comment = 'Comment for log'
s.tags.project = 'My project'
s.tags.user = 'John Doe'
s.tags.tags = ['Tag 1', 'Tag 2/Subtag']

# set timing info
s.wait_between = 0.01
```

```
# set log name and save to disk
s.log_name = 'Test signal demodulation'
s.save('demodulation_scenario')
```

The example will output a file *demodulation_scenario.labber*, which can then be opened in the Measurement program or executed using the ScriptTools API.

A.5.1.2 Example - Detailed description

To describe the various function in more detail, we go through the example line-by-line. We start by creating an empty Scenario and printing the resulting object:

```
>>> from Labber import Scenario
>>> s = Scenario()
>>> print(s)
Scenario:
  instruments: [<Instrument>], #0 items
  channels: [<Channel>], #0 items
  step_items: [<StepItem>], #0 items
  log_channels: [<str>], #0 items
  tags: Tags:
    project:
    user:
    tags: [<str>], #0 items
  settings: Settings:
    send_in_parallel: True
    log_parallel: True
    arm_trig_mode: False
    trig_channel:
    hardware_loop: False
    limit_hardware_looping: False
    n_items_hardware_loop: 1
    update_instruments_if_unchanged: True
    only_send_signal_if_updated: True
    data_compression: 4
    logger_mode: False
  optimizer: Optimizer:
    method: Nelder-Mead
    max_evaluations: 200
    minimization_function: y[0]
    target_value: -inf
    relative_tolerance: inf
    method_settings: {}
  log_name:
  comment:
  wait_between: 0.0
  time_per_point: 0.1
  version: 1.8
```

The print statement lists the properties of the *Scenario* object. These properties fully configure the scenario, and are further described in the *Scenario* class description in the [Scenario class](#) section below.

The properties can be directly modified using standard *Python* notation. For example, the following lines will modify the log comment and the delay setting step channels and measuring log channels in a measurement.

```
>>> s.comment = 'This is a log comment'
>>> s.wait_between = 0.01
```

The first thing we want to do is to add a few instruments to the scenario. This can be done by creating an `Instrument` objects and directly setting it to the `instruments` property of the `Scenario` object. However, it is easier to use the helper function `add_instrument()`:

```
>>> instr_signal = s.add_instrument('Simple Signal Generator', name='Sine')
>>> instr_demod = s.add_instrument('Signal Demodulation', name='Demod')
```

This will create the two instruments and add them to the scenario. To modify the settings of the instruments, we directly update the `values` property of the `Instrument` object:

```
>>> instr_demod.values['Use phase reference signal'] = False
>>> instr_demod.values['Length'] = 1.0
```

Note that the key must match the instrument quantity as defined in the driver definition file. Also note that any undefined quantity values will be initiated to the default values as given in the driver definition file.

The next step is to set the signal connection between the sine waveform and the demodulation input. We do this with the helper function `add_connection()`:

```
>>> s.add_connection('Sine - Signal', 'Demod - Input data')
```

We haven't explicitly defined the channels `Sine - Signal` and `Demod - Input data` used in the signal connection above. The default name for channels follow the convention `<instrument name> - <quantity>`, but it is straightforward to change the name by retrieving a channel with the `get_channel()`-function and then changing its `name` property.

At this point, we are ready to set up the channels to step, and the log channels to measure. This is done with the `add_step()` and `add_log()` functions:

```
>>> s.add_step('Sine - Frequency', np.linspace(0, 10, 51))
>>> s.add_step('Demod - Modulation frequency', start=1, stop=9, step=4)
>>> s.add_log('Demod - Value')
```

Note that the step values can be defined either as a numpy array, or using the keywords `single`, `start`, `stop`, `step`, `n_pts` of the *StepItem* as defined in Section [Scenario class](#) below.

The final thing we need to do is to set the log name and save the scenario to disk:

```
>>> s.log_name = 'Test signal demodulation'
>>> s.save('demodulation_scenario')
```

The resulting file can then be opened in the *Measurement* program or executed using the *ScriptTools* API.

A.5.2 Scenario class

The *Scenario* class contains both properties and helper functions for modifying the configuration.

A.5.2.1 *Labber.Scenario*

```
class labber.config.scenario.Scenario(file_name=None)
```

Class representing a Labber scenario.

The class can be instantiated either as an empty scenario or by loading the Labber scenario provided in the `file_name` input parameter.

Parameters

- **instruments** (*Instrument*, list of) – Configuration of instruments in use in the scenario.
- **channels** (*Channel*, list of) – Channels used in the scenario.
- **step_items** (*StepItem*, list of) – Step items defining channels and values to step or sweep over.
- **log_channels** (*str*, list of) – List of channels to be measured at each step.
- **tags** (*Tags*) – Tags associated with the Labber scenario.
- **settings** (*Settings*) – Measurement settings specific to the scenario.
- **optimizer** (*Optimizer*) – Optimizer settings.
- **log_name** (*str*) – Name of log
- **comment** (*str*) – Comment for scenario
- **wait_between** (*float*) – Time to wait between setting step items and measuring log channels

- **time_per_point** (*float*) – Estimate for time per point, used to calculate duration
- **version** (*str*) – Version of Labber used to create scenario.

__init__(*file_name=None*)

Initialize scenario

Parameters

file_name (*str, optional*) – File with scenario to load, either in .json or .labber format.

add_connection(*source, target*)

Add signal connection between two channels in the scenario.

Parameters

- **source** (*str or Channel*) – Source channel for connection.
- **target** (*str or Channel*) – Target channel for connection.

add_instrument(*driver_name, **kwargs*)

Add instrument to scenario.

Optional keyword arguments are passed on to the Communication object constructor.

Parameters

driver_name (*str*) – Name of driver, must match name in driver database.

Returns

Newly created instrument.

Return type

Instrument

add_log(*channel, index=None*)

Add log item to scenario.

Parameters

- **channel** (*str och Channel*) – Channel for log item. The channel doesn't need to be defined.
- **index** (*int*) – Index of new log item in list. If not given, item is added to end.

add_step(channel, values=None, index=None, **kwargs)

Add step item to scenario.

If the parameter 'values' is not given, additional keywords arguments can be used to initialize the range defining the step item.

Parameters

- **channel** (*str och Channel*) – Channel for step item. The channel doesn't need to be defined.
- **values** (*numpy array, list of float, or float*) – Values for step item.
- **index** (*int*) – Index of new step item in list. If not given, item is added to end.

Returns

Newly create step item

Return type

StepItem

channel_names()

Get list of channels added to the scenario.

The function only returns channels that are active or have been manually added to the configuration. An active channel is used as a step item, log item, or used in a signal connection.

Returns

List of channel names.

Return type

List[str]

get_channel(name)

Get channel by name.

The function can be used to retrieve both active channels and unnamed channels that have not yet been added to the scenario.

For unnamed channels, the name must be of the format “Instrument name - Quantity”. If the instrument/quantity combination is present in the configuration, a new channel will be created and automatically added to the scenario.

Parameters

name (*str*) – Name of channel.

Returns

Channel from scenario.

Return type

Channel

get_config_as_dict()

Create a dict containing the scenario configuration.

Returns

Configuration of scenario.

Return type

dict

get_instrument(name)

Get instrument by name.

Parameters

name (*str*) – Name of instrument to retrieve.

Returns

Instrument from scenario.

Return type

Instrument

get_step(name)

Get step item by name.

Parameters

name (*str*) – Name of step item to retrieve.

Returns

Step item from scenario.

Return type

StepItem

instrument_names()

Get list of instruments present in scenario.

Returns

List of instrument names.

Return type

List[str]

load(file_name)

Load scenario from file.

Parameters

file_name (*str*) – File with scenario to load, either in .json or .labber format.

log_names()

Get list of channel names used as log items.

Returns

List of log names.

Return type

List[str]

remove_channel(name)

Remove channel from scenario.

Note that the function will only remove the channel - the corresponding instrument quantity will still be part of the scenario.

Parameters

name (*str*) – Name of channel to remove.

remove_connection(channel)

Remove signal connection scenario.

Parameters

channel (*str* or *Channel*) – Channel for which to remove connection, can be source or target.

remove_instrument(name)

Remove instrument from scenario.

Parameters

name (*str*) – Name of instrument to remove.

remove_log(channel)

Remove log channel from scenario.

Parameters

channel (*str* or *Channel*) – Log channel to remove.

remove_step(name)

Remove step item from scenario.

Parameters

name (*str*) – Name of step item to remove.

save(file_name, save_as_json=False)

Save Labber scenario to file, either as .labber or .json format.

Parameters

- **file_name** (*str*) – Path to output file.
- **save_as_json** (*bool, optional*) – If True, save to json if no extension is given, by default False

Returns

Final file name, with correct extension

Return type

str

set_log_position(channel, index)

Set position of log item.

Parameters

- **channel** (str or *Channel*) – Channel for log item to move
- **index** (int) – New position for log item in log list

set_step_position(channel, index)

Set position of step item tied to channel.

Parameters

- **channel** (str or *Channel* or *StepItem*) – Channel for step item to move
- **index** (int) – New position for step item in step list

signal_connections()

Get a list of signal connections active in scenario.

Returns

Signal connections, given as list of (source name, target name).

Return type

list of tuple

step_names()

Get list of channel names used as step items.

Returns

List of step names.

Return type

List[str]

A.5.3 Scenario module

This module contains functions and classes for generating Labber scenarios.

A.5.3.1 Enumerations

`class labber.config.scenario.LimitAction`

Enumeration class for actions when channel exceeds limit.

CONTINUE= *'Continue to next step item'*

Continue to next step item

NOTHING= *'Nothing'*

Do nothing

STOP= *'Stop, stay at current values'*

Stop, stay at current values

STOP_RESET= *'Stop, go to init/final configuration'*

Stop, go to init/final configuration

A.5.3.2 Channel

`class labber.config.scenario.Channel(**kwargs)`

Class representing a channel in a Labber scenario.

Parameters

- **name** (*str*) – Channel name.
- **instrument** (*str*) – Instrument used for channel.
- **quantity** (*str*) – Instrument quantity represented by channel.
- **unit_physical** (*str*) – Physical unit of channel
- **unit_instrument** (*str*) – Instrument unit of channel
- **gain** (*float*) – Channel gain, where Instr. value = (Phys. value * Gain + Offset) * Amp
- **offset** (*float*) – Channel offset, where Instr. value = (Phys. value * Gain + Offset) * Amp
- **amp** (*float*) – Channel amplification, where Instr. value = (Phys. value * Gain + Offset) * Amp
- **limit_high** (*float*) – High limit for channel values
- **limit_low** (*float*) – Low limit for channel values
- **limit_action** (*LimitAction*) – Action to take when log channel value exceeds limits

- **signal_source** (*str*) – Channel used as source in signal connections.

get_config_as_dict()

Return the configuration as a dict.

Note that the class variable `_parameter_names` define the list of attributes to include in the dict.

Return type

dict

get_name()

Get name of channel.

If no name is given, the name will be created from the instrument in the form “Instrument - Quantity”.

Returns

Name of channel.

Return type

str

set_signal_source(channel_source=None)

Set channel used as source in signal connection for this channel.

Parameters

channel_source (*str or Channel*) – Channel to be set as source signal. If None, the current signal connection will be removed.

A.5.3.3 Settings

class labber.config.scenario.Settings(kwargs)**

Class representing the settings of a Labber scenario.

Parameters

- **send_in_parallel** (*bool*) – Send values in parallel to multiple instruments.
- **log_parallel** (*bool*) – If True, all channels are measured in parallel

- **arm_trig_mode** (*bool*) – Turn arm/trig mode on/off
- **trig_channel** (*str*) – Trig channel used in arm/trig mode
- **hardware_loop** (*bool*) – Turn hardware loop mode on/off
- **limit_hardware_looping** (*bool*) – Limit hardware looping to first step item.
- **n_items_hardware_loop** (*int*) – Number of step items in hardware loop.
- **update_instruments_if_unchanged** (*bool*) – Update instruments at start even if values are unchanged.
- **only_send_signal_if_updated** (*bool*) – Only send signal if source instrument has been updated.
- **data_compression** (*int*) – Value ranges from 0 (no compression) to 9 (max compression)
- **logger_mode** (*bool*) – If True, object represents a Logger instead of Labber scenario

A.5.3.4 Optimizer

```
class labber.config.scenario.Optimizer(**kwargs)
```

Class representing optimizing settings of a Labber scenario.

Parameters

- **method** (*str*) – Algorithm used for optimization.
- **max_evaluations** (*int*) – Maximum number of function evaluations/measurements before terminating.
- **minimization_function** (*str*) – Function for optimizer to minimize.
- **target_value** (*float*) – Absolute value of minimization function that will terminate optimization.
- **relative_tolerance** (*float*) – Change in value between iterations that is acceptable for convergence.
- **method_settings** (*dict*) – Specific settings for the various optimizer methods.

A.5.3.5 Tags

```
class labber.config.scenario.Tags(**kwargs)
```

Class representing tags of a Labber scenario.

Parameters

- **project** (*str*) – Project name associated with scenario.
- **user** (*str*) – User name associated with scenario
- **tags** (*str, list of*) – List of tags registered to the scenario.

A.5.4 Instrument module

A.5.4.1 Enumerations

`classlabber.config.instrument.Interface`

Enumeration class for defining the communication interface.

ASRL= 'Serial'

Serial - address refers to com port on computer

GPIB= 'GPIB'

GPIB - specify board number in advanced settings

NONE= 'None'

No instrument communication, address is not used

OTHER= 'Other'

Other - address depends on implementation

PXI= 'PXI'

PXI – The address is the PXI slot number. The default chassis number is 1, visible in "Show advanced interface settings." For Keysight products, communication must be established through the Keysight Connection Expert first, and the chassis number should match the ones displayed there. Note that when settings these addresses in the Labber API, the slot number is a string and the chassis number is an integer.

TCPIP= 'TCPIP'

TCPIP - address is TCPIP address

USB= 'USB'

USB - address is USB device name

VISA= 'VISA'

VISA - address is full VISA resource name

`classlabber.config.instrument.Parity`

Enumeration class for defining parity for serial interfaces.

```
EVEN_PARITY= 'Even parity'  
NO_PARITY= 'No parity'  
ODD_PARITY= 'Odd parity'
```

`classlabber.config.instrument.Startup`

Enumeration class for defining the startup operation.

```
DO_NOTHING= 'Do nothing'
```

Leave instrument configuration in its current state

```
GET_CONFIG= 'Get config'
```

Read configuration from instrument at start

```
SET_CONFIG= 'Set config'
```

Set instrument configuration at start

`classlabber.config.instrument.Termination`

Enumeration class for defining termination characters.

```
AUTO= 'Auto'
```

Use system default

```
CR= 'CR'
```

Carriage return

```
CRLF= 'CR+LF'
```

Carriage return + line feed

```
LF= 'LF'
```

Line feed

```
NONE= 'None'
```

No termination

A.5.4.2 Communication

`classlabber.config.instrument.Communication(**kwargs)`

Class representing Labber communication settings.

Parameters

- **name** (*str*) – Instrument name, should be unique.
- **interface** (*Interface*) – Interface type for communication.
- **address** (*str*) – Instrument address, format depends on interface type.
- **startup** (*Startup*) – Operation to perform at instrument startup.
- **server** (*str*) – IP address of server at which instrument is located.
- **lock** (*bool*) – If set, instrument is locked from other users during operation.
- **show_advanced** (*bool*) – Show/hide advanced settings in the instrument configuration window
- **timeout** (*float*) – Maximum time to wait for instrument response.
- **term_char** (*Termination*) – Termination character used by the instrument.
- **send_end_on_write** (*bool*) – Assert end during transfer of last byte of the buffer.
- **lock_visa** (*bool*) – Prevent other programs from accessing the VISA resource.
- **suppress_end_on_read** (*bool*) – Suppress end bit termination on read.
- **tcpip_specify_port** (*bool*) – Use specific TCP port.
- **tcpip_port** (*int*) – TCP port number.
- **tcpip_use_vicp** (*bool*) – Use VICP instead of TCPIP protocol for Teledyne/Lecroy instruments.
- **baud_rate** (*float*) – Communication speed for serial communication.
- **data_bits** (*float*) – Number of data bits for serial communication.
- **stop_bits** (*float*) – Number of stop bits for serial communication, can be 1, 1.5 or 2.
- **parity** (*Parity*) – Parity used for serial communication.
- **gpib_board** (*int*) – The GPIB board number enumeration starts from zero.
- **gpib_go_to_local** (*bool*) – Send GTL over GPIB after closing instrument.
- **pxi_chassis** (*int*) – PXI chassis number.
- **use_32bit_mode** (*bool*) – Run driver in 32-bit mode, for backwards compatibility.

A.5.4.3 Instrument

```
classlabber.config.instrument.Instrument(**kwargs)
```

Class representing the configuration of a Labber instrument.

Parameters

- **hardware** (*str*) – Hardware name, must match instrument driver name.

- **model** (*str*) – Instrument model, must match a model supported by the driver.
- **options** (*str, list of*) – Available instrument options, must match options supported by driver.
- **com_config** (*Communication*) – Communication/interface settings of instrument.
- **values** (*dict*) – Instrument value defining the configuration.
- **version** (*str*) – Version of instrument driver for which configuration is valid.

A.5.5 Step module

A.5.5.1 Enumerations

`classlabber.config.step.FinalAction`

Final action after finishing step

GOTO_FIRST= *'Goto first point'*

Goto first point

GOTO_VALUE= *'Goto value...'*

Goto specific value

STAY_FINAL= *'Stay at final'*

Stay at final

`classlabber.config.step.RangeInterpolation`

Enumeration class for interpolation type of a step item.

LINEAR= *'Linear'*

Linear interpolation.

LOG= *'Log'*

Logarithmic interpolation.

LOGDECADE= *'Log, #/decade'*

Logarithmic interpolation, fixed number of points/decade.

RESONATOR= *'Lorentzian'*

Points are calculated to be equidistant in the complex plane.

`classlabber.config.step.RangeStep`

Enumeration class for step type for a step item.

FIXEDSTEP= 'Fixed step'

Set fixed step size.

N_PTS= 'Fixed # of pts'

Set fixed number of points.

`classlabber.config.step.RangeType`

Enumeration class for defining the range type of a step item.

CENTERSPAN= 'Center - Span'

Center and span values.

SINGLE= 'Single'

Single value

STARTSTOP= 'Start - Stop'

Start and stop values

`classlabber.config.step.StepUnit`

Define step units

INSTRUMENT= 'Instrument'

Define value in instrument units

PHYSICAL= 'Physical'

Define value in physical units

`classlabber.config.step.SweepMode`

Define sweep options of step item

BETWEEN_PTS= 'Between points'

Sweep between fixed points

CONTINUOUS= 'Continuous'

Continuous sweeping

NO_SWEEP= 'Off'

No sweeping

A.5.5.2 *RangeItem*

```
classlabber.config.step.RangeItem(init_value=None, **kwargs)
```

Class representing a single Labber step range item.

Parameters

- **range_type** (*RangeType*) – Range type, can be SINGLE, STARTSTOP, or CENTERSPAN.
- **step_type** (*RangeStep*) – Step length definition, can be either NPTS or FIXEDSTEP.
- **single** (*float*) – Single point value.
- **start** (*float*) – Start point of range.
- **stop** (*float*) – End point of range.
- **center** (*float*) – Center point of range.
- **span** (*float*) – Span of range.
- **step** (*float*) – Step length between points.
- **n_pts** (*int*) – Number of points in the range.
- **interp** (*RangeInterpolation*) – Interpolation type for range.
- **sweep_rate** (*float*) – Sweep rate between points in the range.

```
calc_values()
```

Calculate values for step item.

Returns

Values of step item

Return type

numpy array

```
set_config_from_dict(config)
```

Update config and change range type depending on given settings

Parameters

config (*dict*) – Dictionary with updated values.

`update_parameters()`

Update all parameters (start/end/center/width/etc) to be consistent.

A.5.5.3 *RelationParameter*

```
classlabber.config.step.RelationParameter(for_step_values=False, **kwargs)
```

Class representing a Labber step item relation parameter.

Parameters

- **variable** (*str*) – Parameter name.
- **channel_name** (*str*) – Name of channel represented by parameter.
- **use_lookup** (*bool*) – Turn lookup-table on/off for parameter.
- **lookup** (*LookupTable*) – Lookup-table for parameter.

A.5.5.4 *OptimizerItem*

```
classlabber.config.step.OptimizerItem(**kwargs)
```

Class representing a Labber step item optimizer config.

Parameters

- **enabled** (*bool*) – Enable/disable optimization for this step item.
- **start_value** (*float*) – Start value for optimization process.
- **init_step_size** (*float*) – First step size for optimizer.
- **min_value** (*float*) – Lowest allowed value for optimizer parameter.
- **max_value** (*float*) – Highest allowed value for optimizer parameter.
- **precision** (*float*) – Target precision for optimizer parameter value.

A.5.5.5 *StepItem*

```
classlabber.config.step.StepItem(channel=None, **kwargs)
```

Class representing a Labber step config item.

Parameters

- **channel_name** (*str*) – Name of channel.
- **wait_after** (*float*) – Time (in seconds) to wait after each point.

- **final_value** (*float*) – Value to set after last point. Only relevant if after_last = GOTO_VALUE
- **show_advanced** (*bool*) – Determines if advanced step config dialog is shown by default.
- **use_relations** (*bool*) – Turns relation equation on/off.
- **equation** (*str*) – Equation setting channel relations.
- **step_unit** (*StepUnit*) – Units for step values.
- **after_last** (*FinalAction*) – Final action after finishing last step.
- **sweep_mode** (*SweepMode*) – Define sweep options of step item.
- **use_outside_sweep_rate** (*bool*) – If True, outside sweep rate is set separately from rate between points.
- **sweep_rate_outside** (*float*) – Sweep rate outside sweep range, ie before first and after last point.
- **alternate_direction** (*bool*) – If True, every other step item is executed in reverse order.
- **range_items** (*Rangeltem, list of*) – List with range items defining step values.
- **relation_parameters** (*RelationParameter, list of*) – List with parameters defining relations between channels.
- **optimizer_config** (*OptimizerItem*) – Optimizer configuration for step item.

calc_values()

Calculate and return step values.

Note that the output is the list of values from the range items before applying any relations.

Returns

Step values, before applying any channel relations.

Return type

np.ndarray

update_from_values(values)

Update step item with given values.

Parameters

values (numpy array or list or float) – New values for step item.

A.5.6 Lookup module

A.5.6.1 Enumerations

class `labber.config.lookup.Interpolation`

Enumeration class for defining the interpolation type.

CUBIC= 'Cubic'

Cubick interpolation

LINEAR= 'Linear'

Linear interpolation

NEAREST= 'Nearest'

Nearest x-value

QUADRATIC= 'Quadratic'

Quadratic interpolation

ZERO= 'Zero'

Closest lower x-value

A.5.6.2 LookUpTable

class `labber.config.lookup.LookUpTable(xdata=[], ydata=[], interp=<Interpolation.LINEAR: 'Linear'>)`

Class representing a Labber lookup-table.

Parameters

- **interp** (*Interpolation*) – Interpolation function, default is linear.
- **xdata** (*ndarray*) – X-data for lookup table.
- **ydata** (*ndarray*) – Y-data for lookup table.

calc_values(x)

Calculate values y(x)

Property `x_sorted`

Sorted x-data used by the interpolation function

Return type

ndarray

***property* y_sorted**


Y-data sorted by x-values as used by the interpolation function.

Return type

ndarray

Labber

APPENDIX B: Labber Quantum User Guide



B Appendix B: Labber Quantum User Guide

This appendix provides details on using Labber as a part of Keysight's quantum solution. The configuration will vary depending on requirements, but generally includes the M3202A PXIe Arbitrary Waveform Generator and the M3102A PXIe Digitizer, as well as Pathwave Test Sync Executive and the Quantum Library API software packages.

B.1 System requirements

B.1.1 Recommended Standard Hardware Configurations

Multi-Chassis Configurations (2 or more)	Recommended (minimums)	Preferred	Notes
External Controller	HP Z8	HP Z8	
HP Z8 Configuration	Requirements	Preferred	Notes
BIOS	2.60	Latest	All versions within the recommended test list is sufficient but recommend upgrading to latest when possible
RAM	96 GB		
HDD	500 GB		
CPU	Two lower mid-range Xeon processors	Two upper mid to high range Xeon processors	Must use two CPUs to maximize PCIe connectivity and performance. Other considerations: maximizing the performance of a system does not necessarily mean selecting a processor with the highest available core count. This is because higher core counts generally mean lower turbo frequencies. So for software that is not highly parallel, higher core counts may actually lower performance. Additionally, consideration should be given to Windows licensing which currently has a price increment above 16 core processors for Server OS configurations.
Windows 10 Pro, Windows 10 Pro	Version 10.10.18363	Latest	Always start with latest version. In general, Windows 10 Pro for

Multi-Chassis Configurations (2 or more)	Recommended (minimums)	Preferred	Notes
for Workstations, Enterprise and Windows Server 16+	also known as Version 1909		Workstations is recommended for higher end configurations which includes multi-chassis configurations and/or require NVMe RAID configurations for data storage.

Single-Chassis Configurations	Recommended (minimums)	Preferred	Notes
External Controllers	HP Z4	HP Z8	see Tested PC PXI/AXI Chassis Configurations
HP Z4 Configuration	Requirements	Preferred	Notes
BIOS	HP 61 v02.41	Latest	All versions within the recommended test list is sufficient but recommend upgrading to latest when possible
RAM	64 GB		
SSD	256 GB		
CPU	Lower mid-range Xeon processors	Two upper mid to high range Xeon processors	Must use two CPUs with HP Z8. Other considerations: maximizing the performance of a system does not necessarily mean selecting a processor with the highest available core count. This is because higher core counts generally mean lower turbo frequencies. So for software that is not highly parallel, higher core counts may actually lower performance. Additionally, consideration should be given to Windows licensing which currently has a price increment above 16 core processors for Server OS configurations.
Windows 10 Pro, Windows 10 Pro for Workstations, Enterprise and Windows Server 16+	Version 10.0.18363 also known as Version 1909	Latest	Always start with latest version. In general, Windows 10 Pro for Workstations is recommended for higher end configurations and/or require NVMe RAID configurations for data storage.

Single-Chassis Configuration	
Embedded Controller M9037A	Requirements
BIOS	American Megatrends Inc. AG10 (10/12/2016)
RAM	16 GB
SSD	240 GB
Windows 10 Enterprise 2016 LTSC	Version 10.0.14393

Supporting Hardware Products For External Controllers	Recommended Platform	Notes
M9048A PCIe Desktop Adapter	HP Z4, HP Z8	Single-Chassis
M9049A PCIe High Performance Host Adaptor, Dual Port (x16), Gen 3	HP Z8	Multi-Chassis. Requires Star Connection on PCIe for enumeration. One per chassis up to 4 chassis with dual-cable connection (x16).
M9022A PXIe System Module, Single Port (x8), Gen 3	Single-Chassis	
M9023A PXIe High Performance System Module, Dual Port (x16), Gen 3	Multi-Chassis	Requires Star Connection on PCIe, one per chassis, slot 1 M9024A Not Supported

B.1.2 Labber Software Driver Packages and Dependencies

Depending on the Labber drivers of interest, different sets of dependencies and licenses may be required. See [Installation Guide and Troubleshooting](#) for further instructions.

Software Driver Package	Other Licenses Required	Dependencies
Keysight PXI PWTSE Trigger	KS2201A PathWave Test Sync Executive	Common Software Compatibility
Keysight PXI Sequencer	KS2201A PathWave Test Sync Executive	Common Software Compatibility
Keysight PXI Digitizer Demodulation	KS2201A PathWave Test Sync Executive M5400DMOA (see table below)	Common Software Compatibility; M5400PLSA Software Gateway Compatibility and Dependencies; Instrument Hardware and Firmware
Keysight PXI Agile AWG	KS2201A PathWave Test Sync Executive M5400PLSA (see table below)	Common Software Compatibility; M5400PLSA Software Gateway Compatibility and Dependencies; Instrument Hardware and Firmware
Keysight PXI Digitizer	No	Instrument Hardware and Firmware
Keysight PXI AWG	No	Instrument Hardware and Firmware
Multi-Qubit Pulse Generator	No	
Keysight PXI SMU	No	M960x PXIe instrument driver 3.5.1755.0; Instrument Hardware and Firmware
Keysight PXI Digital IO	No	Common Software Compatibility; M5302A PXIe Digital IO Module Driver 1.0.11601.0; Instrument Hardware and Firmware

B.1.3 M5400xxxA Software/Gateway Compatibility and Dependencies

M5400xxxA	Software Package	API Release	FPGA IP Release	License Required	Dependencies
M5400DMOA	Quantum Library API, FPGA IP Catalog	2.1.2	2.1.1	runtime per module	Common Software Compatibility
M5400PLSA	Quantum Library API, FPGA IP Catalog	2.1.2	2.1.1	runtime per module	Common Software Compatibility

B.1.4 Common Software Compatibility

Vendor	Software/Feature	Release Officially Supported	License Required
Keysight	KS2201A PathWave Test Sync Executive	1.4.15 (2020 Update 1.1)	Yes
Keysight	SD1 3.X	3.02.52	No
Keysight	IO Libraries Suite	18.2.26526.0 (2021)	No
Python	Python software programming language	3.7.x	No

B.1.5 Instrument Hardware and Firmware

Vendor	Module	Options	NOT supported	Firmware
Keysight	M3202A PXIe Arbitrary Waveform Generator 1 GSa/s, 14 bit, 400 MHz	-K41, -K32 -HV1, -FP1, -Mxx, CLF	-CLV option	04.02.11 ONLY
Keysight	M3102A PXIe Digitizer 500MSa/s, 14 bit, 200 MHz	-K41, -K32 -HV1, -FP1, -Mxx, CLF	-CLV option	02.02.06 ONLY
Keysight	M9019A PXIe Chassis: 18-slot, 3U, 24 GB/s, Gen 3	chassis driver M9019A 1.7.402.1		M9019A Chassis Firmware Version 2019EnhTrig
Keysight	M9018B PXIe Chassis: 18-slot, 3U, 8 GB/s, Gen 2	chassis driver M9019A 1.7.402.1	M9018A PXIe Chassis	
Keysight	M9010A PXIe Chassis: 10-slot, 3U, 24 GB/s, Gen 3	chassis driver M9019A 1.7.402.1		M9010A Chassis Firmware Version 2019EnhTrig
Keysight	M9615A PXIe 5-channel Precision Source/Measure Unit			Core: 2.26 SMU: 4.56
Keysight	M5302A PXIe Digital IO Module			5.7.30.0

B.2 Installation Guide and Troubleshooting

B.2.1 Quantum IP Library

The Quantum IP Library is required to run any of the following drivers:

- Keysight PXI Agile AWG (to be used with Keysight PXI Sequencer)
- Keysight PXI Digitizer Demod

Like Labber, licenses for the Quantum IP Library are through the PathWave License Manager. See [Installing and Labber License](#) for instructions on installing a license through this system.

A license for the Quantum IP Library is required on a per PXIe instrument module basis.

- Obtain license file, depending on the count (number of modules):
 - **Single count:** This can be obtained by requesting a trial evaluation license or purchasing license. Once the license request is approved, a license file with a .lic extension is sent as an email attachment. Save this file on your computer. Trial Licenses can be obtained from the Quantum Library website [here](#).
 - **Multi count:** Purchase a multi count license or request a trial license to your Application Engineer specifying the specific product number, number of modules to be used in the system and MAC address.
- Once you obtained the licenses, install the Quantum IP Library. This includes a software programming Quantum IP Library API, PathWave FPGA Quantum IP Library and Keysight PathWave License Manager.

B.2.2 Configuring a Python Environment for Labber Drivers

Labber ships with a built-in Python 3.7 environment for running instrument drivers. It is located in the installation directory - on Windows this is typically 'C:\Program Files\Labber\python-labber'. When a driver is started, a separate Python process is launched using the executable pythonw.exe located in this directory.

The built-in Python environment contains all of the required packages for running the drivers that ship with Labber. However, custom drivers may require additional packages not included by default. One solution is to create a user-defined Python environment which contains all the required packages. Labber can be configured to use this environment instead of its own. In the Edit menu, select 'Preferences...', then in the 'Advanced' section, enter the path the pythonw.exe executable in your environment of choice; for example, 'C:\Users\Administrator\Anaconda3\envs\py37-driver\pythonw.exe'. If this field is left blank (the default), the Labber Python environment is used for all new driver processes. The recommended package manager for creating and maintaining custom Python distributions is Anaconda.

The driver updates included in this release of Labber include a few new dependencies. To ensure that all drivers run properly, the following packages must be installed into your environment:

- `keysight_hvi`
 - The API for Keysight PathWave Test Sync Executive
 - After installing Test Sync Executive
 - navigate to the installation directory (default: 'C:\Program Files\Keysight\PathWave Test Sync Executive 2020 Update 1.1\api\python')
 - In your Python environment, run 'pip install .'
- `keysightSD1`
 - The API for Keysight SD1 3.X
 - After installing Keysight SD1 3.X
 - navigate to the installation directory (default: 'C:\Program Files\Keysight\SD1\Libraries\Python')
 - In your Python environment, run 'pip install keysightSD1-3.X.Y.tar.gz'
- `PyQuLibrary`
 - The API for the Keysight Quantum IP Library (PyQuLibrary)
 - After installing the Quantum IP Library
 - navigate to the installation directory (default: 'C:\Program Files\Keysight\QuantumLibApi\2.1.2')
 - In your Python environment, run 'pip install .'
- `filelock`
 - Can be installed with 'conda install filelock'

B.2.2.1 Additional dependences

The following dependences are also required by some or all of the drivers that ship with Labber. These were present in older versions of Labber, so a previously-working environment likely already has these installed. Most of these packages can be installed with conda. If otherwise noted, install with pip.

- `comtypes`
- `numpy`
- `scipy`
- `h5py`
- `pycrypto`
- `future`
- `scikit-learn`

- dill
- pyvisa (pip)
- qtpy (pip)
- scikit-optimize (pip)
- zhinst (pip)

B.2.3 Multi-chassis configuration

The Labber driver Keysight PXI PWTSE Trigger can be used to operate up to four PXIe chassis (M9010A/M9019A) with PathWave Test Sync Executive for synchronized triggering of up to 62 M3202A Arbitrary Waveform Generators and M3102A Digitizers.

In order to ensure proper operation of a multi-chassis system, there are a few configuration requirements.

B.2.3.1 Firmware and software


Ensure that the correct version of PathWave Test Sync Executive, the Keysight SD1 instrument driver software and firmware, the Keysight IO Libraries, and the chassis firmware, as described in System requirements have been installed.

Instructions for updating chassis firmware can be found here: [PXIe Chassis Firmware Update Guide](#) starting on page 17.

B2.5.2 PCIe connections and chassis/module discovery

PCIe connections are used to connect and transfer data between PXIe chassis and a host PC, where a PCIe adapter is placed in each chassis (M9022A, M9023A, or M9024A). These can be connected in a star-configuration with two dual-slot host PC PCIe adapters (M9049A), or in a cascade configuration with dual-slot adapter in the chassis (M9023A, or M9024A). More information on supported configurations including supported host PCs can be found at <http://www.keysight.com/find/PXIAXIeTestedPC>.


After configuring the PCIe connections, chassis and modules should be visible in the Keysight Connection Expert program:



For troubleshooting information: <https://literature.cdn.keysight.com/litweb/pdf/M9019-90005.pdf?id=3024867>

B.2.3.2 Chassis enumeration

Chassis may not appear numbered as expected. Ensure that the numbering of chassis as seen in Connection Expert is the same as desired for the PXIe trigger sharing cascade (see next section). Chassis can be assigned new numbers in the 'Chassis Numbers' page:



Numbering should be serial and begin at 1.

B.2.3.3 PXI Trigger Line Sharing


To share real-time trigger information between chassis, Keysight M9031 cards are used to connect share PXI lines between chassis. More information is available in the TSE User Manual.

Important: to avoid line delays and trigger distortion, M9031 cards should be connected to one another with the shortest possible SMB cables. Use 390 mm or 12 inch SMB cables. This will require all chassis to be stacked atop one another in the same equipment rack. Ensure all chassis have rubber feet, leaving a gap between chassis for proper airflow.


B.2.4PXI Trigger Reservations

When updating a system from M3601A (HVI) to KS2201A (PathWave Test Sync Executive), the situation is sometimes encountered that some of the PXI trigger lines in the M90XX chassis remain reserved by M3601A or some other program. This will cause an error when turning on the Output of the Keysight PXI PWTSE Trigger driver, with an error message about being unable to reserve trigger lines.

To check if the fault is old trigger reservations, open the Keysight Connection Expert program, choose "PXI/AXIe Chassis" → "Chassis Triggers."



This tool will display the existing trigger reservations for each of the eight PXI lines in each of the two or three segments:



All PXI_TRIG lines should be vacant. If any are reserved by "KtlviDriver_Default," these will need to be cleared. If you have trouble clearing the lines, contact Keysight customer support.

B.2.5 Common Problems When Getting Started with TSE/Quantum IP Library


When getting started with PathWave Test Sync Executive and the PathWave Quantum IP Library it is essential to ensure that one has compatible versions of the following:

- Labber and Labber drivers used for instance:
 - a. Keysight PXI Agile AWG
 - b. Keysight PXI AWG
 - c. Multi-Qubit Pulse Generator
 - d. Keysight PXI Digitizer
 - e. Keysight PXI Digitizer Demod
 - f. Keysight PXI Sequencer
- PathWave Test Sync Executive (HVI)
- Keysight SD1
- Quantum IP Library (PyQuLibrary)


Keysight Quantum Engineering Solutions is committed to continuous improvements to its software solutions. As such it is quite easy to become a few versions behind the latest software that may not be backwards compatible. Please check the earlier part of this Users' Guide to see minimum requirements for compatibility.

B.2.5.1 Checking Labber Version

To check for the version of Labber being used, open the Instrument Server.




From there click on Help and then About:




B.2.5.2 Checking Labber Python Version

Labber installs with its own python distribution and uses that by default. The minimum compatible version is python 3.7. To see what python distribution Labber is using, go to Edit, then Preferences:




Then in the Preferences window go to the Advanced Tab:



If the python distribution is left blank, then Labber is using the distribution that ships with Labber which is located in the installation directory. This can be overwritten by the user if you have configured your path variables differently. As such be careful on the environment being used to run Labber and if one is observing that the default is not being used it may be because of settings in your path.

B.2.5.3 Checking PathWave Test Sync Executive (HVI) Version

When starting the Keysight PXI Sequencer or Keysight PXI PWTSE Trigger the installed version and its installation location of PathWave Test Sync Executive are displayed (when in Basic Log level) in the Instrument Server Instrument Log. To change the log level in the Instrument Log. Go to the Instrument Server Preferences. Then on the Server tab, change the Instrument log level to 'Basic'.



From the Instrument Log when starting the Keysight PXI Sequencer:

```
16:30:25,451: Keysight Pathwave TSE (HVI) Python API Version: 1.4.7
16:30:25,452: Keysight Pathwave TSE (HVI) Python API Install Location: C:\Program Files\Labber\python-labber\lib\site-packages\keysight_hvi\__init__.py
16:30:25,452: Keysight SD1 Python API Install Location: C:\Program Files\Keysight\SD1\Libraries\Python\keysightSD1.py
```

For help on installing, uninstalling, etc please reference the relevant installation guide.


B.2.5.4 Checking Quantum IP Library Version

When starting the Keysight Agile AWG or the Keysight Digitizer Demod driver the Quantum IP Library Python API version and its installation location are displayed (when in Basic Log level) in the Instrument Server Instrument Log.

```
09:24:16,532: Agile AWG: Quantum IP Library Version: 1.0.0
09:24:16,532: Agile AWG: Quantum IP Library Location: C:\ProgramData\Anaconda3\envs\py37-driver\lib\site-packages\PyQULibrary\__init__.py
```

For help on installing, uninstalling, or other issues please reference the relevant installation guide.

B.2.5.5 Common PathWave Quantum IP Library Installation Issue



If one encounters the above error when trying to install the Keysight PathWave Quantum IP Library from the GUI then that means your python environment has not been configured in an amenable way for the installer.

However, it is easy to do a "manual installation" of the Python API. One simply needs to:

- Go to the new installation directory selected for the Quantum IP Library. In this example it was 'C:\Program Files\Keysight\QuantumLibApi\1.0.0'. Note this will install into a location set by your environment path variable.
- Type 'pip install .' which will install the Quantum IP Library
 - More information on this is available in the PathWave Quantum IP Library Labber Users' Guide.

B.2.5.6 Checking the Keysight SD1 Version for PXI Instruments

We have also included, when in Basic Log Level, information on the individual AWG or Digitizer Firmware and Hardware SD1 is using:

Example Keysight PXI AWG:

```
16:35:21,027: AWG Hardware Version: 04.13.00  
16:35:21,028: AWG Firmware Version: 04.01.20
```


Example Keysight PXI Digitizer:

```
16:31:53,960: Digitizer Hardware Version: 04.09.00  
16:31:53,960: Digitizer Firmware Version: 02.01.60
```

B.2.5.7 Wrong Version of Dependency Still Being Used Even After Install

If you have installed a new version of PWTSE (HVI) or SD1 but the Instrument Log still shows the older version there are a few things to try.

- Stop and restart the Instruments. This should reload the python definitions used by the instrument driver.
 - a. In rare circumstances this may fail and one would need to first stop the instrument to prevent any hanging connections. Then, remove the instrument from the Instrument Server. Then, add the Instrument back to the Instrument Server.
- If after doing the above the version is still incorrect. Use the Windows Registry Editor to display the installed versions of TSE/SD1/QIP.
 - To use the Registry Editor, go to the search bar in the bottom left and type: Registry Editor



- If nothing is found, use the Windows command prompt and type 'regedit'

Registry Edit Examples:

Quantum Library API


	Name	Type	Data
ab	(Default)	REG_SZ	(value not set)
ab	Location	REG_SZ	C:\Program Files\Keysight\QuantumLibApi\1.0.0
ab	Timestamp	REG_SZ	***unknown variable timestamp***
ab	Version	REG_SZ	1.0.0

HVI Core/Pathwave Test Sync Executive

	Name	Type	Data
ab	(Default)	REG_SZ	(value not set)
ab	CurrentVersion	REG_SZ	1.0.18+0
ab	InstallLocation	REG_SZ	C:\Program Files\Common Files\Keysight\PathWave Test Sync Executive\HVIcore\1.0

	Name	Type	Data
ab	(Default)	REG_SZ	(value not set)
ab	CurrentVersion	REG_SZ	1.0.18+0
ab	InstallLocation	REG_SZ	C:\Program Files\Keysight\PathWave Test Sync Executive 2020

Keysight SD1



Name	Type	Data
(Default)	REG_SZ	(value not set)
baseDir	REG_SZ	C:\Program Files (x86)\Keysight\SD1
baseDir64	REG_SZ	C:\Program Files\Keysight\SD1
driver_version	REG_SZ	1.40.000000
qt_version	REG_SZ	5.12.8
version	REG_SZ	3.00.95

This will show if there are multiple installations/versions of the same software. If there are multiple installations, one should uninstall the software the undesired software by first using the Keysight installation managers. If that fails, use the Windows Control Panel. If that also fails, as a last resort, delete it from the Registry Editor.

B.2.5.8 Python API Version in Instrument Server Log isn't Updated to Latest Download

If after downloading and installing a new version of either PathWave TSE (HVI) or Quantum IP Library and the version displayed in the Instrument Log is the older version the following methods should be tried to resolve the issue:

- Make sure to run 'pip install .' in the correct location to create the new Python API for PathWave TSE (HVI) or the Quantum IP Library.
- If the problem persists, stop and restart the Instrument.
- If the problem persists, right click and reload the Instrument definitions.
- If the problem persists, stop and remove the Driver from the Instrument Server. Add the Driver back to the Instrument server.
- If the problem persists, make sure one is using the latest released version of Labber. Check www.keysight.com/find/labber.

B.2.5.9 Digitizer not added to Sequencer

When one is configuring an Agile Sequence, it is essential to make sure that all components are included. For instance, if one does not enter in the Chassis/Slot information for the digitizer an error will be thrown as shown below.

Throws Error message:




An error occurred when sending a value to an instrument.

Instrument name: Keysight PXI Sequencer (HVI)
Address: PXI:
Quantity: Pulse Sequence
Value: [0.00e+00 1.01e-07, ..., nan nan]

Error message:
Must include at least one digitizer in the sequence

Unfiltered error message:
File "C:\Program Files\Labber\python-labber\multiproc-include\py37\BaseDriver.py", line 451, in processOperation
self.performOperation(dOp, callId)
File "C:\Program Files\Labber\python-labber\multiproc-include\py37\BaseDriver.py", line 493, in performOperation
bWaitForSweepInFunction=dOp['wait_for_sweep'], options=dOp)
File "C:\Program Files\Labber\python-labber\multiproc-include\py37\BaseDriver.py", line 800, in _performSetValue
return self._performSetValueNoNone(quant, value, 0.0, options)
File "C:\Program Files\Labber\python-labber\multiproc-include\py37\BaseDriver.py", line 790, in _performSetValueNoNone
quant, value, sweepRate=sweepRate, options=options)
File "C:\Users\Administrator\Labber\Drivers\Keysight_PXI_Sequencer\Keysight_PXI_Sequencer.py", line 303, in performSetValue
self.program_hvi()
File "C:\Users\Administrator\Labber\Drivers\Keysight_PXI_Sequencer\Keysight_PXI_Sequencer.py", line 568, in program_hvi
self.program_hvi_sequence()
File "C:\Users\Administrator\Labber\Drivers\Keysight_PXI_Sequencer\Keysight_PXI_Sequencer.py", line 813, in program_hvi_sequence
self.define_registers()
File "C:\Users\Administrator\Labber\Drivers\Keysight_PXI_Sequencer\Keysight_PXI_Sequencer.py", line 724, in define_registers
raise Exception('Must include at least one digitizer in the sequence')

If one opens the Instrument Driver for the Keysight PXI Sequencer, it becomes apparent what the issue is when inspecting the Readout Configuration.



B.3 Keysight PXI Labber Drivers

B.3.1 Keysight PXI AWG

B.3.1.1 Hardware looping for DC offsets

The Keysight PXI AWG driver supports one-dimensional hardware loop scans for DC offsets, similar to hardware looping for waveforms. A loop of N DC offsets is accomplished by

- creating a series of N short, constant waveforms with amplitude equal to the desired offset

- setting the AWG to hold on the last sample of the waveform until the next trigger
- uploading and queueing the N waveforms to play in a loop

To turn an ordinary DC scan into a hardware loop scan:


- Enable Hardware Loop mode
- Ensure the AWG channels used are in AWG mode, not DC mode

This mode of operation supports sweeping the channel offset directly or through a scalar signal connection, which can be useful when scanning relative to a calibration.

B.3.1.2 Multi-dimensional hardware loop sweeps

If scanning two DC offsets in a 2D grid, open the 'Show Settings' pane, select the option 'Limit hardware looping to first step item,' and set the 'Number of step items in hardware loop' to 1. See image below. This will ensure that all the offsets for the inner loop (first step channel) are played in a loop for each step of the outer loop (second step channel).


Multi-dimensional sweeps will only work properly if the hardware looping is limited to the first step item.



B.3.1.3 Example: 2D DC scan

In this example we show how to set up a 2D DC offset scan with inner hardware loop. In this example, the measurement time is reduced from 75 seconds to 12 seconds by using hardware loop.

This example uses two channels on the same AWG, but the setup is the same if the channels are on separate AWGs.




Important notes:

- In this example, a waveform from the Single-Qubit Pulse Generator ('Pulses') is connected to each of the AWG channels. This is optional, but serves to show that this combination is valid. In this case the DC offsets will simply be added to the single supplied waveform.
- The Amplitude of the AWG channels must be set to be at least as large as the range of DC offsets used, since the Amplitude sets how large the waveform amplitude can be. In this case the Amplitude for Channel 1 is set to 1.5 V because the scan ranges from -1.2 V to +1.2 V, plus an additional small-amplitude waveform.
- The AWG trigger modes must be set to 'Software / HVI', and the trigger channel must be the 'Output' channel of the Keysight PXI PWTSE Trigger driver.

Results:


Channel 1 inner loop for first outer loop step (but the output is independent of outer step)



We see a few features here:

- The step up at around 300 ns is the start of the waveform which implements the DC offset. On the first Offset value, the value before the step is zero because there is no waveform playing yet. On successive steps, the value holds from the last waveform. This example uses only one average – on successive averages the value before the first step would be the *last* value, +1.2 V.
- The oscillatory pattern is a 20 MHz square pulse added to the DC offset, with a delay of 200 ns chosen to demonstrate the timing more clearly.

Channel 2 inner loop for outer loop step 10 (+250 mV):



The output is held at +250 mV for all values of the inner loop, except for the oscillatory waveform.

B.3.2 Keysight PXI PWTSE Trigger

B.3.2.1 Driver Usage


The new Keysight PXI PWTSE Trigger driver for Labber is designed to support software synchronized experiments using the PathWave Test Sync Executive product line.

This page describes the use of the PathWave Test Sync Executive Trigger using the Labber interface, including multi-chassis triggering.

B.3.2.1.1 Loading and setting a HVI Instrument

- In the Instrument Server, add a Keysight PXI PWTSE Trigger driver
- Double click to open the driver settings panel. In the Communication section:
 - Select the address of the PXI chassis
 - You can specify a name for the instrument (here we choose *HVI*)
 - Select the model of the PXI chassis (current options: M9010 for 10 slots, M9019 for 18 slots)
- In the Settings sections:
 - Trigger:

- Output (default = True): when on, the HVI is running. Usually, that is the desired mode. Labber will toggle the output when running a measurement in arm/trig mode (see below)
- Trig period (s): specify the trigger interval in steps of 10 ns
- Digitizer delay (s): specify an optional delay between AWG and digitizer triggers
- Configuration:
 - Multi-chassis: when on, a panel to set the inter-chassis links will be displayed (see below)
 - Auto-detect: when on, it will scan the chassis for all connected modules.
 - Scan for devices: Scan to get a list of connected modules (only shown when Auto-detect is off)
 - Slot #-#: if auto-detect, or scanned for devices, this list will be populated with the detected devices (either AWG or Digitizer). Alternatively, one can select the desired modules directly from this menu.






Fig. 2.4.1: Instrument Settings window for Keysight PXI PWTSE Trigger, with a M9010 chassis (10 slots). The slot configuration can be set on section '1st chassis' (highlighted in red), or detected automatically with the 'Auto-detect' option or by pushing the 'Scan' button.

B.3.2.1.2 Using HVI in a Measurement

Once the HVI instrument is configured as explained above, one can add set up a measurement as usual using the Measurement Console. The HVI is added from server, in addition to the all the PXI modules and other instruments needed in the experiment.

The variable quantities (Trig period and Digitizer delay) can be set or stepped like any other quantity in Labber.

The only supported measurement mode is the Arm/trig mode (with or without the hardware loop option). This mode prevents potential conflicts between HVI and AWG

(or Digitizer) drivers, by stopping the HVI sequence before a measurement starts, and at every step of a sweep.

To activate the arm/trigger mode and the optional hardware looping (see section 6.11 of the Labber Manual), one needs to add a HVI quantity in the step sequence (see Fig.2.4.2). This will enable the HVI - Output option in the Arm/trig mode menu, which is required to start the HVI only once the digitizer is armed. This, in combination with the hardware loop option, allows for the acquisition to be synchronized with the looping hardware, such as different steps in a sequence played by the AWGs.




Fig.2.4.2: Example of a Measurement Setup for HVI use in hardware loop mode. The stepped quantity (here Single Qubit Pulse Generator - Width #1) will be looped through at an interval set by Trig period, until the set Digitizer - Number of averages is reached. If the hardware loop box is not checked, instead, the program will acquire the desired number of averages at every step, before moving to the next pulse width.

B.3.2.1.3 Multi-Chassis Usage

This driver supports up to four interconnected chassis. For a description of the setup that is necessary to connect multiple chassis, refer to the [KS2201A PWTSE User Manual](#).

When checking the 'Multi-chassis' box, the driver scans for all the available chassis, overriding the address provided in the Communication section. When combined with the 'Auto-detect' option, or by pushing the 'Scan' button, the corresponding Chassis sections are populated with all the connected modules (AWGs or Digitizers). Alternatively, one can enable/disable modules manually, provided that the conditions for the minimum number of modules to ensure communication between chassis is maintained (see same link above).

In addition, a new set of 'Interchassis link' dialogue boxes appears, allowing the user to input the chassis and slot number of each pair of connected M9031A modules, mapping to the hardware configuration. By default, these modules are configured to connect the chassis in ascending order of address. If the system detects fewer chassis than those supported by the driver, or if the user wants to disable some of them manually, the chassis numbers for the extra Interchassis links are set to 0 (indicating that they are disabled).

Once the Interchassis setup is complete, the PWTSE driver runs as described above for a single chassis.




Fig.2.4.3: Example configuration showing two Interchassis links, sufficient to connect 3 chassis with each other. Each link is specified by a pair of source and destination M9031A modules, each identified by a chassis and slot number.

B.3.3 Keysight PXI Digitizer Demod

A driver which uses Keysight FPGA IP and Quantum Library API to replace and extend the functionality provided by the Keysight PXI Digitizer FPGA Demod driver.

Supports up to 10-tone simultaneous demodulation with rectangular or custom integration windows. This driver requires an HVI sequence to trigger acquisition, and should be used with the Keysight PXI PathWave Test Sync Executive Trigger or the Keysight PXI Sequencer driver.

B.3.3.1 Requirements

This driver requires that the Quantum Library Python API (PyQuLibrary) be installed in the Python environment used by Labber drivers. For more information, see [Configuring a Python Environment for Labber Drivers](#).

B.3.3.2 Input

2 channels of IF signal (that is, I and Q outputs of an IQ mixer). Only input channels 1 and 2 are used by this driver

B.3.3.3 Output

Up to 10 demodulated IQ values, both average and single-shot

B.3.3.4 Driver parameters

B.3.3.4.1 General parameters

- **Number of samples:** DOUBLE
Number of samples per acquisition
- **Number of records:** DOUBLE
Number of distinct records per step in a Labber sweep
- **Number of averages:** DOUBLE
Number of repetitions per record and per step in a Labber sweep
- **Input range:** COMBO
Input voltage range, in Volts. Applied to channels 1&2. One of 4, 2, 1, 500m, 250m, 125m, 62.5m

Note: In the current version of the Quantum Library API, the coupling, impedance, and range are set automatically. The coupling and impedance are hard-coded to be set to 50 Ohm and DC, respectively, so the driver enforces this.

B.3.3.4.2 Trigger settings

- **Trig mode:** COMBO
Set the trigger mode of the demodulators. One of "Software/HVI" or "Digital trigger." "Digital trigger" mode uses the trigger input on the front panel of the module

- **External Trig Config:** COMBO
Trigger mode for digital front panel trigger. One of "Rise," "Fall", "High," or "Low"
- **Trig Sync Mode:** COMBO
Trigger clock synchronization setting. One of "None" or "Ext/PXI 10 MHz." Sets whether the trigger is sampled using the internal 100 MHz clock or the PXI backplane clock, respectively.

B.3.3.4.3 Demod settings

- **Number of qubits:** COMBO
the number of frequency-domain multiplexed (FDM) channels. From 1 to 10
Per channel demod settings. # refers to the index of the FDM channel, ranging from 1 to **Number of qubits**.
- **Window type #:** COMBO
Type of window function. 'Rectangular' (created by driver) or 'Custom' (loaded from file)
- **Demodulation frequency #:** DOUBLE
center frequency of demod window in Hz. When **Window type** == 'Custom' this is an optional additional frequency shift applied to window from file (see **Shift custom window**)
- **Demodulation phase #:** DOUBLE
phase of demod window in degrees. When **Window type** == 'Custom' this is an optional additional phase shift applied to window from file (see **Shift custom window**)
- **Window file # :** STRING
full file name of file specifying the window. Either .npz (single complex array) or .npz (dict w/ multiple windows, keyed as 'window_#')
- **Shift custom window #:** BOOLEAN
If True, custom window is frequency and phase shifted using **Demodulation frequency #** and **Demodulation phase #**. If False, unaltered array is used.
- **Window gain #:** DOUBLE
When set to 1, the demodulated amplitude should be representative of the RMS amplitude of the incoming signal. For a constant pulse, the demodulated amplitude will be equal to the amplitude of the signal. Increase the gain or reduce the input range if encountering noticeable signal discretization. See **Output resolution and window gain**.
- **Threshold:** DOUBLE
the value (in V) of the threshold used to discriminate state 1 (for shots above the threshold) from 0 (below). See also **State Assignment**.
- **Discrimination quadrature:** COMBO
the chosen quadrature (I or Q) to be converted to assigned qubit states.

- **Invert sign:** BOOLEAN
invert state assignment if selected. (0 if above threshold).

B.3.3.4.4 Driver outputs

- **IQ_# - Single shot:** VECTOR_COMPLEX
Array of length **Number of averages * Number of records** containing demodulated complex IQ point for each record and repetition
- **IQ_#:** VECTOR_COMPLEX
Array of length **Number of records**, containing the average demodulated complex IQ points.
- **Mag_#:** VECTOR
Average demodulated absolute magnitude: $\text{mean}(\text{abs}(\text{IQ_}\# - \text{Single shot}))$. Not equal to $\text{abs}(\text{IQ_}\#)$.
- **window_#:** VECTOR_COMPLEX
The window function loaded to the FPGA. Can be accessed and logged for visualization and verification
- **Assigned state_# - Single shot:** VECTOR
Array containing the assignments of each value in IQ_# - Single shot to 0 or 1, if that value is below or above the set Threshold, respectively.
- **Assigned state_#:** VECTOR
Average of Assigned state_# - Single-shot

B.3.3.5 Use model within Labber

Labber supports multiple operation modes, defined by the selection of "Hardware trigger mode" and "Hardware loop mode". These use cases are summarized here.

Non-hardware trigger mode

Not supported. This would require configuring the FPGA IP blocks while the HVI sequence is running, which is not supported by the hardware. This means that data acquired in the Instrument Server interface is not guaranteed to be configured properly. Likewise, with a front panel trigger, the trigger source must be stopped while configuring the FPGA.

Hardware trigger mode (w/out hardware loop)

This is the most flexible use mode, though it is sub-optimal for long scans where many waveforms must be loaded to the AWG(s) (see Hardware loop mode, below).

In "Software/HVI" trigger mode with the trigger channel set to "HVI Trigger - Output", the HVI trigger sequence is stopped and unloaded from hardware (Output set to "Off"), all instrument settings are applied, the HVI sequence is loaded and started, then data is read. In this case, this happens at each step of the experimental scan. Averaging is done at each point, with the number of averages N set by the Digitizer Demod driver. This

allows the HVI sequence to be started after the IP blocks are registered and configured, which ensures that the data acquired are the first N shots of the sequence.

If using an external digital trigger, the trigger source should correspond to a setting which toggles the output of that source.

This mode supports up to 3 million averages per step, irrespective of trigger source.

Hardware loop mode

For sweeps with many unique waveforms loaded to the AWGs, this is recommended, as it is by far the most efficient mode of operation.


The demodulation parameters cannot be scanned in hardware loop mode, because there is no way to pre-load them - they must be sent to the FPGA using the Quantum Library API, which cannot be done in real-time.

However, if the scan involves only AWG waveform parameters, this is allowed, because the IP blocks need only be configured once. In this case, the Digitizer Demod driver expects to acquire all of the data for the scan at the end of the measurement. The data is reshaped and returned in the way Labber expects.

A limitation in this mode is the maximum number of shots allowed by the DataWrapppper (3,121,257). With a large scan and many averages, this can be hit easily. For now, use Hardware trigger mode to get around this (which is unfortunately slower). See **Max number of samples**.

B.3.3.6 Demodulator layout

The layout of the FPGA includes up to 10 demodulator blocks, each of which is configured independently via the settings described above. Each demodulator block is composed of 4 downsamplers (shown below) which perform point-wise multiplication and integration of one of the input channels against the real or imaginary part of the complex integration window $w(t)$, and are used to produce I and Q of the output demodulated point.



**NB: Qubit State Decoder not included in this version. All "+/-" blocks perform addition, as the sign is handled by the complex window creation.*

Since the FPGA performs multiplication with real numbers, this is accomplished with four simultaneous multiplication steps, added pairwise as shown in the above diagram to produce the complex output. From top to bottom, the multiplication and integration performed is

$\sum Re(w(t)) \times Ch_1(t)$
$\sum -Im(w(t)) \times Ch_2(t)$
$\sum Im(w(t)) \times Ch_1(t)$
$\sum Re(w(t)) \times Ch_2(t)$

Rows 1 and 2 are added together to produce output *I*, and rows 3 and 4 are added together to produce output *Q*.


B.3.3.7 Switching between raw and demod modes

There are two Labber drivers used to control M3102A digitizer modules:

- Keysight PXI Digitizer
- Keysight PXI Digitizer Demod

The first is used to acquire raw time trace signals ("raw mode"). The second performs real-time demodulation of input signals with up to 10 independent integration windows ("demod mode"), used to demodulate components of a frequency-multiplexed signal. If using both modes in different measurements, the Instrument Server will likely be running

both of these drivers for each physical digitizer. Below is an example using three digitizers in a four-chassis configuration:



B.3.3.8 Switching from demod to raw signal acquisition

The two digitizer drivers require different FPGA sandbox images. The loading and unloading of the sandbox region is handled by the Digitizer Demod driver. The bitfile is loaded when the instrument is started, and unloaded when the instrument is stopped. This means the raw mode driver will not work properly when the demod mode driver is running. To run an experiment with the raw mode Digitizer driver, make sure the corresponding Digitizer Demod driver is stopped. In the above image, the driver "Digitizer1-15" will not function as expected.

The opposite situation is allowed: the Digitizer Demod driver will work properly when the Digitizer driver is running. **However, it is recommended to set the configuration of the Digitizer driver after closing the Digitizer Demod driver.** This is because the Digitizer Demod driver can set parameters like the input range, which cannot be queried. Stopping and starting the Digitizer driver, or running a measurement which sets the configuration, is sufficient to ensure the actual instrument settings match the Labber driver parameters.

B.3.3.9 Errors to watch for

Before some stability improvements were made in Labber, there were a few common errors. Certain resource locks are in place which should prevent these errors. Please report them to the Labber development team if they are encountered.

- TimeoutError of the form "Driver timeout waiting 10 seconds for another process to release resource pxi_module_1-3. Try increasing this driver\'s timeout in the "advanced interface settings" section of the "Communication" settings." This would indicate that one Labber driver is waiting for another to finish accessing the PXI module in chassis 1, slot 3. Try increasing the timeout of this driver. If the error persists, or occurs after an abrupt shutdown of Labber, there may be uncleared locks preventing access to the PXI module. Check for .lock files like pxi_module_1-3.lock in the Labber "Temporary items" folder (by default, ~\AppData\Roaming\Labber). If found, stop all running Labber PXI drivers and delete all such .lock files.
- Windows Error: previously, this error would sometimes occur at the start of a measurement when the PXI Digitizer Demod and PXI PWTSE Trigger try to access the same physical resources on the M3102A module. Try stopping all Labber drivers and then starting them again, one by one.
- "operation "HVI port reset" failed:" This error was sometimes found after the previous Windows Error, and prevented further access to the M3102A module. If this error occurs and the user is then unable to connect to the digitizer by opening the Labber driver, a power cycle of the chassis is required.
- PC freezes: previously, this error would sometimes occur when the PXI Digitizer Demod and PXI PWTSE are both started at the same time. If encountered, a power cycle of the host PC and chassis is required.

B.3.3.10 Output resolution and window gain

The Keysight PXI Digitizer Demod Labber driver allows for easy control of demodulation windows to optimize the performance for multiplexed signals. The relevant parameters are the "Input range" in the Settings tab, and the "Window gain" in the "Demod params" tab. There is a separate window gain for each of the 10 frequency channels.

The DownsamplerX5N produces a complex number which consists of two unsigned 16-bit integers. The Quantum Library API converts these integers to 16-bit *signed* doubles, and multiplies them by the input range (between 62.5 mV and 4V), which is set in the Labber driver. Consequently, each quadrature of the resulting IQ point has 16 bits of resolution. For example, with a 1V input range, the full scale is 2 V (-1 to +1), and the resolution is $2 \text{ V} / 2^{16} = 30.5 \text{ uV}$. For convenience, the resolution in each range is summarized below:


Input range	Output resolution
4 V	122 uV
2 V	61 uV
1 V	30.5 uV
500 mV	15.2 uV
250 mV	7.6 uV
125 mV	3.8 uV

Input range	Output resolution
62.5 mV	1.9 μ V

The Input range parameter sets both the input range of the ADCs and the output range of the demodulator. This is usually convenient, since the demodulator aims to yield an amplitude (in Volts) which is consistent with that of the input signal. More specifically, if the input signal and integration window are both of constant amplitude (modulo an IF modulation) and the window gain is set to 1, then the output value will have the same amplitude of the input signal. For signals with time-varying amplitudes, one often uses a match integration window which also varies in time. If the input signal is matched to the integration window and the gain is set to 1, then the output amplitude will be equal to the RMS average of the input signal. Thus, for optimal resolution, the input range should be chosen to be the minimum which accommodates the largest expected signal amplitude.


However, if the input signal contains many frequency components, as is the case when demodulating a multiplexed readout signal, then the peak input voltage can be many times larger than any individual signal. In this case, the input range must be chosen to be larger to avoid clipping of the input signal. This leads to a coarsening of the resolution, which can lead to unwanted discretization of the demodulated output.

In the example below, we demodulate a 1 μ s constant pulse at 5 MHz with a Rectangular window. The input signal has an amplitude of 20 mV, and we set the input range to 250 mV, which would accommodate 10 such tones. However, the resultant signal suffers from discretization error



The resolution of 7.6 uV is clearly visible in each quadrature. The amplitude is $(14.9^2 + 13.1^2)^{.5} = 20$ mV, as expected.

The "Window gain" parameters in the "Demod params" tab is used to compensate for this. Since the input range is 10 times larger than would be necessary for a single carrier tone, one should increase the window gain by a factor of 10. This factor increases the amplitude of the integration window and more efficiently utilizes the available output range:



Note that the signal amplitude has been multiplied by 10 in each quadrature, but the phase is unaffected.

B.3.3.11 Max Number of Samples

The number of datapoints which can be recorded is finite, due to limitations in the Data Acquisition (DAQ) blocks on the digitizer. The maximum number of times data can be pushed to a given DAQ is slightly more than one million. The DataWrapperX5 packages demodulated results in groups of 3×5 (3 IQ points per frequency-domain channel). This means the max number of IQ points which can be acquired per frequency channel is thus $3 \times 1\text{M}$ (actually 3,121,257). Asking more than this many samples will force an error in the configuration of the acquisition, and the Labber driver will close.

Ordinarily, this means the number of averages must be kept less than 3 million. However, in Hardware Loop mode, the digitizer acquires a sample for each step of the loop before averaging. This means the total number of datapoints is $M \times N$, where M is the size of the loop and N is the number of averages. In this case, the max number of averages must be less than $3,121,257/M$. For instance, if sweeping a signal amplitude over 100 points, the max number of averages is about 30,000.

B.3.3.12 Example: Multiple Digitizers

This example demonstrates the use of the Digitizer Demod driver in a multi-chassis configuration.

In this example, we have three M3102A digitizers and five M3202A AWGs installed across four chassis. Below, the AWG channels are arranged in IQ pairs where I and Q signal originate from different modules. The IQ pairs are processed by the Demod driver. Note, none of the modules used in this example are in chassis 3.

```
# 3 digitizers
dig_addresses = [
    (1, 15),
    (2, 16),
    (4, 18)
]

# AWG channels operate in 3 pairs (chassis, slot, channel number)
pair_addresses = [
    ((1,4,1), (2,8,1)), # -> Digitizer 1-15
    ((2,5,1), (1,7,1)), # -> Digitizer 2-16
    ((4,7,1), (1,4,2)) # -> Digitizer 4-18
]
```

We will play the same pairs of waveforms (a frequency-modulated IQ signal) on each pair of AWG channels, and analyze them with the digitizers. We will sweep the frequency of the signal.

B.3.3.12.1 Instruments

The example Labber configuration, name `freq_sweep_mixed_reference.hdf5`, is included in the reference examples.

This measurement uses 5 AWGs and 3 digitizers. Signals can be acquired and demodulated on the M3102A using the Keysight PXI Digitizer Demod driver ("demod mode"), or just acquired using the Keysight PXI Digitizer driver ("raw mode"). Digitizers 2-16 and 4-18 are in demod mode, and Digitizer 1-15 is in raw mode.

Channels				
Instrument	Name/Address	Instr. value	Phys. value	Server
> FDMReadout	readout			localhost
> Keysight PXI AWG	AWG2-5			localhost
> Keysight PXI AWG	AWG2-8			localhost
> Keysight PXI AWG	AWG1-4			localhost
> Keysight PXI AWG	AWG4-7			localhost
> Keysight PXI AWG	AWG1-7			localhost
> Keysight PXI Digitizer Demod	Digitizer Demod 2-16			localhost
> Keysight PXI Digitizer Demod	Digitizer Demod 4-18			localhost
> Keysight PXI Digitizer	Digitizer1-15			localhost
> Keysight PXI HVI Trigger	trigger			localhost

The HVI sequence which synchronizes the instruments is created and executed by the **Keysight PXI PWTSE Trigger driver**. The slots with M3XXX modules installed are populated in the interface:

Keysight PXI HVI Trigger - PXI: 1 - trigger

At startup:

Sections

- Communication
- Settings
- 1st chassis
- 2nd chassis**
- 3rd chassis
- 4th chassis

Instrument settings

Slot 2-1:	-
Slot 2-2:	-
Slot 2-3:	-
Slot 2-4:	-
Slot 2-5:	AWG
Slot 2-6:	-
Slot 2-7:	-
Slot 2-8:	AWG
Slot 2-9:	-
Slot 2-10:	-
Slot 2-11:	-
Slot 2-12:	-
Slot 2-13:	-
Slot 2-14:	-
Slot 2-15:	-
Slot 2-16:	Digitizer
Slot 2-17:	-
Slot 2-18:	-

☐ Send/receive
 ☐ Active

Some modules are represented in the PWTSE Trigger driver slot options, but not in the Labber Measurement configuration. This is because the HVI sequence needs to use a minimum set of modules for multi-chassis communication, but this particular Labber measurement does not send to data to certain unused AWGs. For more information, see section "Multi-chassis Usage" in the **Keysight PXI PWTSE Trigger** documentation.

B.3.3.12.2 Signal source

The signal is generated by an FDMReadout driver, which can create a frequency-multiplexed tone. In this case, we only use a single carrier frequency, and play a square pulse. The relevant FDMReadout instrument settings are

- Pulse length: 10^{-6} s (the length of the square pulse)
- Acquisition length: 10^{-6} s (the total length of the waveform, here chosen to match Pulse length)
- Pulse amplitude: 0.5 V

The Signal Connections between the FDMReadout driver (called "readout") and the AWG channels are as shown below, reflecting the connections between the AWGs and the digitizers above.

Signal connections		
Target	Source	Channel
▼ AWG2-5		
Ch1 - Waveform	readout	Output I
▼ AWG2-8		
Ch1 - Waveform	readout	Output Q
▼ AWG1-4		
Ch1 - Waveform	readout	Output I
Ch2 - Waveform	readout	Output Q
▼ AWG4-7		
Ch1 - Waveform	readout	Output I
▼ AWG1-7		
Ch1 - Waveform	readout	Output Q

B.3.3.12.3 Acquisition and demodulation

Relevant settings for the digitizer instruments are:

- Number of samples: the length of the actual acquisition in samples. The sampling rate of the M3102A is 500 MS/s, so the length in nanoseconds is twice the number of samples. This is chosen separately for each digitizer
 - Digitizer 4-18 (raw mode): 750 samples, to show the full trace
 - Digitizer 1-15 (demod mode): 500 samples, the full pulse length
 - Digitizer 2-16 (demod mode): 200 samples, only part of the pulse. This will show the different frequency responses of the two different integration windows.
- Input range: 1 V for all channels
- For demod mode only,
 - Number of qubits: The number of independent frequency domain demodulation channels used. Can range from 1 to 10. Here we choose 2.

- Demod channel settings (see image below). Integer k goes from 1 to "Number of qubits"
 - Demodulation frequency k: the center frequency of the demodulator for frequency channel k We choose two different frequencies, **-10 MHz** and **+25 MHz**
 - Window type k: here, "Rectangular," meaning the complex integration window is of the form $\exp(2\pi j f)$, for demodulation frequency f. Labber generates this window automatically
 - Window gain k: The gain of the demodulation stage (typically set to 1). See [Output resolution and window gain](#) for more detail.

The image shows a software interface with two sections, 'Qubit 1' and 'Qubit 2'. Each section contains four settings: 'Demodulation frequency [Hz]', 'Demodulation phase [deg]', 'Window gain', and 'Window type'. For Qubit 1, the frequency is -10E6, phase is 0, gain is 1, and type is Rectangular. For Qubit 2, the frequency is 25E6, phase is 0, gain is 1, and type is Rectangular. Each setting is in a text box with a small up/down arrow icon.

B.3.3.12.4 Measurement settings and steps

The step and log quantities are shown here:

Step sequence

Channel	# pts.	Step list	Output range
readout - Modulation frequency 1	241	-30 MHz - 30 MHz	-30 MHz - 30 MHz
trigger - Trig period	1	20 us	20 us
trigger - Digitizer delay	1	300 ns	300 ns

☒ Arm/trig mode: trigger - Output
☒ Hardware loop
Edit...
Remove

Log channels

Channel	Instrument	Address
Ch1 - Signal	Keysight PXI Digitizer	PXI: 15
Ch2 - Signal	Keysight PXI Digitizer	PXI: 15
IQ_1	Keysight PXI Digitizer Demod	PXI: 2:16
IQ_2	Keysight PXI Digitizer Demod	PXI: 2:16
IQ_1	Keysight PXI Digitizer Demod	PXI: 4:18
IQ_2	Keysight PXI Digitizer Demod	PXI: 4:18

☒ Log in parallel
Edit...
Remove

In this example, we perform a one-dimensional sweep of the signal frequency from -30 MHz to +30 MHz with a 250 kHz step size.


As explained in Use model, we must operate in Hardware arm/trigger mode, triggering on the start of the HVI sequence created by a **PXI PWTSE Trigger driver**. We further choose to operate in Hardware Loop mode to speed up the measurement.

We log the raw signals on channels 1 and 2 of the digitizer in raw mode, and log the average IQ results for the two frequency domain channels for the digitizers in demod mode. The PXI Digitizer Demod driver also provides the single-shot IQ data (for example, "IQ_1 - Single shot"), which we do not log here.


The quantity "trigger - Digitizer delay" sets the relative delay between the HVI triggers for the AWGs and for the digitizers. This is chosen to approximately align the start of the demodulation window with the arrival of the pulse from the AWGs.

B.3.3.12.5 Results

The acquired data is shown in the plot below. First, we see the raw signal measured on channel 2 (Q quadrature) of digitizer 1-15, demonstrating the frequency sweep (time axis zoomed in for clarity):



Next, we show the response of the demodulators. We plot the magnitude of the two demodulation channels versus the signal frequency. Digitizer 2-16 is the top pane, digitizer 4-18 is on the bottom. Frequency channel 1 (-10 MHz) is in orange, frequency channel 2 (+20 MHz) is in black.



We see the characteristic sine-like response, with each channel centered at its demodulation frequency. Note that the response of digitizer 4-18 is wider, due to its shorter time-domain acquisition. Note that in all cases, the peak amplitude is 0.5 V, the amplitude of the input signal.

B.3.3.13 Example: Window Calibration

This example demonstrates the calibration of custom demodulation windows.

The general procedure is as follows:

1. Acquire reference trajectories for different qubit states using the raw signals digitizer in a calibration scenario.
2. Construct window functions by using the complex conjugate of the difference of these two trajectories.
3. Apply the computed custom windows in the Demod driver and it use in a verification scenario.

In a real-world calibration, this example can be adapted by replacing the FDMReadout driver with a Multi-Qubit Pulse Generator in a "Readout training" sequence, and replacing the "State" steps with a single step of "Training, input state".

This section walks through this example, which is provided as a reference example called `calibrate_envelopes.py`. The example uses one digitizer and one pair of AWG channels.

B.3.3.13.1 Hardware configuration

We must specify the mapping between the physical in/out channels and the logical IQ pair. We also specify inter-chassis links if using multiple chassis. To run this example, the user must edit the file `test_config.py`, the contents of which are shown below:

```

# -----
#      user-defined system configuration
# -----

# list of ( source chassis, source slot, destination chassis, destination slot ) tuples
interchassis_links = [(1, 12, 2, 12), (2, 13, 3, 12), (3, 13, 4, 12)] # list of source chassis
- source module - destination chassis - destination module tuples
multichassis = len(interchassis_links) > 0

# this example uses 3 digitizers
dig_addresses = [
    (1, 15),
    (2, 16),
    (4, 18)
]
n_dig = len(dig_addresses)

# AWG channels operate in 3 pairs
awg_addresses_all = [(1, 4), (1, 7), (2, 5), (2, 8), (3, 6), (3, 9), (3, 17), (4, 7), (4, 10)]
pair_addresses = [
    ((1, 4, 1), (2, 8, 1)), # -> Dig1-15
    ((2, 5, 1), (1, 7, 1)), # -> Dig2-16
    ((4, 7, 1), (1, 4, 2)) # -> Dig4-18
]

# need all the modules in the HVI sequence so the interchassis links operate
properly
# these will not all be added to Labber config
awg_addresses_all = [(1, 4), (1, 7), (2, 5), (2, 8), (3, 6), (3, 9), (3, 17), (4, 7), (4, 10)]
dig_addresses_all = [(1, 15), (2, 16), (4, 18)]

```

This example will calibrate envelopes for a particular IQ pair used for readout, chosen from one of the above pairs. This choice is made at the top of the script:

```
# -----  
#      digitizer selection  
# -----  
  
# of available signal pairs in test_config, which one to use  
i_pair = 0  
# digitizer address used for readout - ( chassis, slot )  
dig_addr = dig_addresses[i_pair]  
  
# AWG channels used for readout operate in an IQ pair - ( chassis, slot, channel  
number )  
pair_addresses = pair_addresses[i_pair]
```

B.3.3.13.2 General settings

There are a number of global settings used by both the calibration and verification measurements (time and frequency units like ns, us, MHz, etc. are defined earlier in the script to conform with Labber conventions):

```

# -----
#      measurement parameters
# -----
n_qubits = 5 # can go up to 10
n_reps = 20e3
n_reps_cal = 40e3
pulse_len = 1*us
acq_len = 2*us
trig_period = 10*us
sample_rate_awg = 1*GHz
dig_delay = 300*ns
n_samples = round(acq_len / DT)

# signal amplitude and input range
sig_amp = 0.005
range_str = '62.5 mV'
range_val = 0.0625
window_gain = range_val / sig_amp # >1 b/c the range is much larger than the signal,
for multiplexed tests

# Labber measurement settings
hardware_trig = True # only supported when True
hardware_loop = False

```

These parameters are collected into the relevant instrument settings, not shown in this walk-through.

B.3.3.13.3 Readout channel settings

This example uses the FDMReadout driver for signal generation. Here we include the settings for the readout channels, including frequency and dispersive readout parameters. Parameters for all channels have different frequencies, phases, and resonator linewidths

```

freqs = np.linspace(-85e6, 95e6, 10)[:n_qubits]
# add some arbitrary frequency staggering, so they're not necessarily equally spaced
df = 3e6
dfs = (2*np.random.rand(n_qubits) - 1) * df
freqs = freqs + dfs

# resonator linewidths
kappa_min, kappa_max = 0.2e6, 3e6
kappas = np.random.rand(n_qubits) * (kappa_max - kappa_min) + kappa_min

# input signal phases
phases = np.random.rand(n_qubits) * 360 - 180

# parameters for the FDMReadout object are defined like this
ch_params1 = {
    'Modulation frequency 1' : freqs[0],
    'Signal phase 1' : phases[0],
    'kappa 1' : kappas[0],
    'chi 1' : -.5e6,
    'delta0 1' : -.25e6,
    'Pulse length 1' : pulse_len,
    'Amplitude 1' : sig_amp
}
# we define the other channels likewise

```

B.3.3.13.4 Connect to Instrument Server

To switch between raw and demod mode, we need access to some of the instruments used in the measurements. See note [Switching between raw and demod modes](#).

```

# -----
#      connect to server and get/create instruments
# -----
print('Connecting to server instruments')

server = Labber.connectToServer('localhost')
# create new instruments if they do not exist on the server

# digitizer demod driver
inst_dig_demod = server.connectToInstrument('Keysight PXI Digitizer Demod',
{
    'name' : demod_name,
    'interface' : 'PXI',
    'PXI chassis' : chassis,
    'address' : str(slot),
    }, bCreateNew=True)

```

B.3.3.13.5 Creating Measurement Scenarios

Most of the instrument settings and signal connections are the same in the calibration and verification scenarios. The function 'create_generic_scenario' prepares a Labber measurement Scenario with those general settings, signal connections, and log channels, as well as some specific one which depend on the type of measurements. 'create_generic_scenario' has a Boolean keyword argument 'demod_mode,' which controls the mode of the digitizer. The calibration scenario sets this to be False, and the verification scenario sets it True. Here, we walk through this function.

First, we add the Instruments:

```

s = Scenario()
fdmr = s.add_instrument('FDMReadout', name=readout_name)
hvi = s.add_instrument('Keysight PXI PWTSE Trigger', name=trig_name,
    interface='PXI', address=1)

for addr in awg_addresses_exp:
    name = awg_name(*addr)
    this_awg = s.add_instrument('Keysight PXI AWG', name=name,
        interface='PXI', pxi_chassis=addr[0], address=addr[1])

# add digitizer to Measurement scenario
# choice of driver depends on demod_mode
chassis, slot = dig_addr
if demod_mode:
    demod = s.add_instrument('Keysight PXI Digitizer Demod', name=demod_name,
        interface='PXI', pxi_chassis=chassis, address=slot)
else:
    dig = s.add_instrument('Keysight PXI Digitizer', name=dig_name,
        interface='PXI', pxi_chassis=chassis, address=slot)

```

Settings for each of these instruments are also added here (not shown).

The measurement includes the following signal connections between the readout pulse generator and the AWG(s):

```

# add signal connections between signal generator and AWGs
addrI, addrQ = pair_addresses
nameI = awg_name(*addrI[:2])
s.add_connection(readout_name+' - Output I', nameI+' - Ch{} -
Waveform'.format(addrI[-1]))
nameQ = awg_name(*addrQ[:2])
s.add_connection(readout_name+' - Output Q', nameQ+' - Ch{} -
Waveform'.format(addrQ[-1]))

```

We set the log channels, which depends on the type of measurement.

```

if demod_mode:
    # log demodulated signals
    for n in range(n_qubits):
        s.add_log(demod_name+' - IQ_{} - Single shot'.format(n+1))
else:
    # log raw signals
    s.add_log(dig_name+' - Ch1 - Signal')
    s.add_log(dig_name+' - Ch2 - Signal')

```

Finally, we set the measurement step sequence, which will flip different qubit states in the FDMReadout signal generator. This is done differently in the calibration and verification measurements, and we show them below.

We also provide a function 'run' for executing a measurement, and 'collect_data' for getting the results.

B.3.3.13.6 Calibration scenario

We create the raw calibration scenario and run it:

```

s = create_generic_scenario('calib', demod_mode=False)

# save scenario
cal_scenario_file = os.path.join(output_folder, fn_raw + '.labber')
s.save(cal_scenario_file)

# run scenario
output_fn = os.path.join(output_folder, fn_raw + '.h5')
run(cal_scenario_file, output_fn)

```

Since we wish to know the change in average signal due to each qubit individually, we flip them one at a time. This can be done in Labber using a dummy variable as an index for the qubits and using channel relations (also accessible in the Measurement Editor under 'Show advanced settings' in the Step configuration dialog):

```

dummy_steps = list(range(0, n_qubits+1))
s.add_step(dummy_variable_name, dummy_steps)
for n in range(1, n_qubits+1):
    step = s.add_step(readout_name+' - State {}'.format(n))
    step.set_config_from_dict({
        'use_relations' : True,
        'equation' : 'd1 == {}'.format(n) # only when dummy variable is equal to n
    })
    relation = step.relation_parameters[0]
    relation.set_config_from_dict({
        'variable' : 'd1',
        'channel_name' : dummy_variable_name
    })

```

Thus, our step sequence has $n_{\text{qubits}}+1$ steps. In this example, we use the Timer driver for the dummy variable.

B.3.3.13.7 Construct custom windows

We load the result of the calibration experiment, and compare the response for the two qubit states individually for each qubit:

```

# get raw signals from the calibration measurement
data_names = [dig_name+' - Ch1 - Signal', dig_name+' - Ch2 - Signal']
[dataI, dataQ] = collect_data(output_fn, data_names)
dataIQ = dataI + 1j*dataQ

# indicies of data traces used from above sweep to prepare reference trajectories
idx = [ [0, n] for n in range(1, n_qubits+1) ]

window_dict = {}
for i, (i_g, i_e) in enumerate(idx):
    traj_g = dataIQ[i_g]
    traj_e = dataIQ[i_e]
    # conjugate of difference creates matched filter at correct frequency
    window = np.conj( traj_e - traj_g )
    window_name = 'window_{}'.format(i+1)
    window_dict[window_name] = window

# save window file
np.savez(window_fn, **window_dict)

```

These windows are also plotted and saved as windows.png:

B.3.3.13.8 Create verification scenario

Now we create a scenario with 'demod_mode=True' to be used to verify our window functions.

```
s = create_generic_scenario(fn_demod, demod_mode=True, window_type='Custom')

# auxiliary settings
demod = s.get_instrument(demod_name)
for i in range(n_qubits):
    demod.values['Window type {}'.format(i+1)] = 'Custom'
    demod.values['Window file {}'.format(i+1)] = window_fn
    demod.values['Shift custom window {}'.format(i+1)] = False
```

Since the demodulation channels are independent, we can flip all of the qubits at the same time, which makes this verification measurement more efficient and demonstrates the independence:

```
s.add_step(readout_name+' - State 1', ['0', '1'])
for n in range(2, n_qubits+1):
    # step in relation to State 1
    step = s.add_step(readout_name+' - State {}'.format(n))
    step.set_config_from_dict({
        'use_relations': True,
        'equation': 's1'
    })
    relation = step.relation_parameters[0]
    relation.set_config_from_dict({
        'variable': 's1',
        'channel_name': readout_name+' - State 1'
    })
```


We also add another step which scans the demodulation window type from 'Custom' to 'Rectangular' to see the difference that using a correctly configured custom window has:

```

# step window type to see both rectangular and custom
s.add_step(demod_name+' - Window type 1', [1, 0])
for n in range(2, n_qubits+1):
    step = s.add_step(demod_name+' - Window type {}'.format(n))
    step.set_config_from_dict({
        'use_relations' : True,
        'equation' : 'w1'
    })
    relation = step.relation_parameters[0]
    relation.set_config_from_dict({
        'variable' : 'w1',
        'channel_name' : demod_name+' - Window type 1'
    })

```


We also plot the measured windows:




B.3.3.13.9 Run verification scenario

We save and run the above scenario and compare the results of the custom and rectangular windows. Below are plotted the 2D histogram of IQ points for one frequency channel, with rectangular and custom windows.

Rectangular:



Custom:



The custom window orients the signal along the correct axis, and improves signal to noise (not obvious in this example).

B.3.3.14 Known Issue: Phase Stability

The phase of a demodulated signal is sensitive to relative timing jitter or instability between modules. Though the modules all share a common 10 MHz reference, the modules themselves use internal phase-locked loops to boost this to the 100 MHz clock

for HVI instruction execution, and other higher-frequency clocks for sampling. Thus these higher-frequency clocks are not necessarily perfectly synchronized.

The M3102A and/or M3202A modules appear to exhibit oscillation of the timing of instructions. This manifests as oscillation of the phase of the demodulated signal from shot to shot. The magnitude of the phase shift depends on the IF frequency. For instance, a time shift of 100 ps corresponds to a phase shift at 100 MHz IF of $(0.1 \text{ ns}) \cdot (0.1 \text{ GHz}) = 0.01 \text{ cycles}$, or 3.6 degrees.

The magnitude of the oscillation varies from system to system.


- On a single-chassis system, we have measured the peak-to-peak variation to be on the order of 20 ps, or 0.7 degrees at 100 MHz, peak-to-peak.
- On a four-chassis system, the variation is at most 100 ps, or 3.6 degrees at 100 MHz, peak-to-peak.

The period of the oscillation is roughly the same on the two systems, and is about 0.3 ms, or 3 kHz.

B.3.3.14.1 Example measurement on multiple chassis

We can see evidence of this oscillation with a simple measurement. We play a single-tone pulse on the AWGs, which are directly connected to the digitizers. By repeatedly measuring the phase from pulse to pulse, we can detect drifts or oscillations in the timing of the pulse and/or demodulator.

Below we see an example of phase oscillation with three digitizer across three chassis, all acquiring the same pulse at 85 MHz IF, sourced from different AWGs. The vertical axes are phase in degrees, while the horizontal axis is in samples. The measurement cadence is was 10 μ s, so 1000 samples is 10 ms, and we see about 30 oscillations per 10 ms, a 0.3 ms period. The magnitude of the oscillation is largest on the second digitizer, at most about 3 degrees.




The oscillations are not phase-locked between digitizers. This measurement does not explain whether the origin of the oscillation is the AWG, the digitizer, or both; nor does it determine whether the absolute timing of triggers is changing (phase modulation), or if the relative sample rates are changing (frequency modulation). Further tests suggest the jitter is in the pulse/demodulation timing, not the real sample rate.

B.3.3.15 State Assignment

The demodulated and integrated I/Q outputs can be digitized in software according to selected thresholds.

By enabling the Assign state option below, each single-shot result will be assigned to a qubit state. Each qubit has the following set of options:


- **Threshold** : the value (in V) of the threshold used to discriminate state 1 (for shots above the threshold) from 0 (below).
- **Discrimination quadrature**: the chosen quadrature (I or Q) to be converted to assigned qubit states. You can use a custom integration window (see **Window type**) to ensure that the signal is predominantly in one of these two quadratures.
- **Invert sign**: invert state assignment if selected. (0 if above threshold).




Selecting the Assign state option enables a number of "Soft reset" options (Post-processing section) to to post-process the binary outcomes. These options are:

No reset. No post-processing after state assignment.

Post-select. Discard outcomes when the first measurement in a sequence returns 1. This requires a sequence with at least an additional measurement in the beginning, and Number of records > 1 (see [Creating a custom sequence](#)). Both single-shot and averaged Assigned state traces ignore the records following the 1st, when the 1st record is assigned to 1.




Post-process. Similar to post-select, but the assigned state for each record after the 1st is inverted (0 to 1 and 1 to 0), instead of discarded.



Restless. Invert state assignments based on the last measurement record. This mode does not require additional measurement pulses in the sequence. It always uses the most recent measurement record, whether it results from a different measurement pulse in the sequence, or from the previous iteration of the same sequence. The very first record in each continuous acquisition is left unchanged (that is, it assumes that the qubit begins in the ground state). Note that this causes a different behavior whether a sweep is in hardware-loop mode or not.

- In a hardware-loop sweep, the acquisition is continuous: the measurement repeatedly loops through all the steps until the desired Number of averages is reached. The sequences are equally spaced from each other, with the end-to-start time approximately set by Init time in the Keysight PXI Sequencer. Thus, all measurements occur at deterministic times, and it is always possible to correlate a result with the previous one.
- In a sweep without hardware-loop, the program repeats each step Number of averages times, then interrupts the acquisition, and rearms instruments before proceeding to the next step. Because of this interruption, the first measurement in the next step is not conditioned on the result of the last measurement in the previous step.



Example, with Number of records = 3, Number of averages = 1, for a sweep consisting of 4 steps.

Assigned State_1 - Single shot results, for different selections of Soft reset:

	Step 1	Step 2	Step 3	Step 4
No reset	0, 1, 0	0, 1, 1	1, 0, 1	1, 1, 0

	Step 1	Step 2	Step 3	Step 4
Post-select	0, 1, 0	0, 1, 1	1, - , -	1, - , -
Post-process	0, 1, 0	0, 1, 1	1, 1, 0	1, 0, 1
Restless, hw-loop ON	0, 1, 0	0, 1, 0	0, 1, 1	0, 0, 1
Restless, hw-loop OFF	0, 1, 0	0, 1, 0	1, 1, 1	1, 0, 1

B.3.4 Multi-Qubit Pulse Generator for Agile sequences

This guide illustrates the usage of the Labber Multi-Qubit Pulse Generator driver (MQPG) to generate sequences of individually triggered pulses. This is made possible by the combination of M5400PLSA and KS2201A PathWave Test Sync Executive (PWTSE) - see requirements below - which are both integrated in the Keysight PXI Agile AWG and Keysight PXI Sequencer Labber drivers, respectively.

B.3.4.1 Agile sequencing - Driver configuration

B.3.4.1.1 Multi-Qubit Pulse Generator

The MQPG driver can produce sequences to support two modes of operation with PTSE. In the first mode, the MQPG produces the full waveforms for each AWG in the sequence, comprising all the pulses with defined timing and shape properties. These waveforms are simultaneously triggered at a specified trigger interval using the Keysight PXI PWTSE Trigger. In this mode, the user connects the output of each waveform channel of the MQPG to a waveform channel of the Keysight PXI AWG driver.

In the second mode, which we refer to as Agile sequencing, the MQPG defines pulse sequences in a similar way, but each pulse in the sequence can now be treated as a separate entity. This allows for the control of the following pulse properties at runtime: amplitude, frequency, phase, and timing (the latter only for selected sequences). Also, an initialization time, defined as the interval between the end of a sequence and the beginning of a new one, can now be set in place of a fixed trigger period. In Agile sequencing mode, the MQPG is used with the Keysight PXI Agile AWG driver (see below).




Fig.B3.4.1: Example of a CPMG sequence for 1 AWG using the waveform and the agile sequencing modes. In the waveform mode, the full sequence is stored and triggered as a single waveform. In the agile mode, a single waveform wf0, corresponding to an individual pulse, is stored in the AWG memory. The sequence is constructed by setting timing, amplitude and phase of each instance of that pulse.

To enable the Agile mode, select the **Generate waveform primitives** option in the Primitives section:




Fig. B3.4.2: Setting Agile sequencing mode

The Qubit primitives and Readout primitives groups specify which pulse properties are pre-compiled into waveform memory (checked **Apply phase/amplitude modulation to qubit primitives** and **Apply frequency/phase/amplitude modulation to readout primitives**) or set at runtime (unchecked). By default, frequency, phase, and amplitude are all set at runtime, which means the value of these parameters will be encoded in the HVI sequence. This choice allows for a tradeoff between the real-time setting of those values (see below), which more-efficiently utilizes waveform memory, and the minimum pulse distance achievable (enforced by the compiler). Pre-compiling the qubit control frequency is not supported. Instead, frequency modulation is always applied and updated at runtime to maintain a constant rotating reference frame.

The onboard FPGA waveform memory of the Keysight PXI Agile AWG driver can store up to 64 unique waveforms per output channel, which can be played and reused in arbitrary order in a sequence. The option to set parameters at runtime provides more efficient use of this limited memory, as a single primitive can be used to effect different operations on a qubit by appropriate selection of amplitude and phase.

Hardware-dynamic parameters

When waveform parameters are set at runtime, there is a further option to specify the values of these parameters in the HVI sequence when it is compiled (default), or to use the hardware-dynamic option, which allows this value to be modified after compile time. Activating a hardware-dynamic option has two effects:

- 1) The specified quantity can be swept linearly in hardware-loop mode. This allows for efficient scans where the set values are stepped through at runtime in a single HVI sequence, removing the need of uploading new waveforms at every step. See Hardware-loop dynamic sweeps for an example.
- 2) The specified quantity can be set to a different value between measurement runs without uploading new waveforms or programming a new HVI sequence.

This option is enabled in the Hardware-dynamic parameters group. By selecting the generic **Enable hardware-dynamic amplitude** and/or **Enable hardware-dynamic frequency** options, XY or Z pulses with dynamic properties (amplitude and/or frequency for XY, amplitude for Z) can be specified in the 1-QB gates XY or 1-QB gates Z section, respectively. In the following example, the amplitude of all XY gates for Qubit #1 and the modulation frequency of Qubit #2 can be swept in hardware by selecting **Hardware-dynamic amplitude** in group Pulse #1 and **Hardware-dynamic frequency** in group Pulse #2, respectively.




Fig. B3.4.3: Selecting qubit-specific dynamic quantities

In addition to pulse- and qubit-specific dynamic settings, built-in sequences can also have dynamic properties. For example, selecting the 'CP/CPMG' sequence exposes a few new options for sequence-specific hardware-dynamic parameters:




Fig. B3.4.4: Selecting sequence-specific dynamic quantities

Dynamic last pulse phase. If the phase is not pre-compiled in the pulse primitives (corresponding box is unchecked in Fig. B3.4.2), this option enables the fast sweep of the phase of the second $\pi/2$ pulse in the sequence.

Dynamic duration. This option enables the fast sweep of the total sequence time.

Finally, one can define custom sequences with an arbitrary choice of dynamic parameters, as outlined in Section B.3.4.6.2.

B.3.4.1.2 Keysight PXI Sequencer

This driver contains the graphical interface used to set up the mapping between logical and physical channels.



The Settings section specifies:

- **Number of qubits and Number of readout buses.** The buses are defined as a way to group multiple measurement channels into the same physical IQ pair (i.e. for readout multiplexing).
- **Number of chassis:** if > 1, enable the configuration of links between multiple chassis. See Keysight PXI PWTSE Trigger.
- **Output:** Set whether the sequence is running. This is typically controlled by the Measurement program in hardware loop mode.
- **Init time (s):** Minimum time between the end of the last pulse in a sequence and the beginning of a new iteration.
- **Digitizer delay (s):** Delay of the digitizer trigger with respect to the measurement pulse. If positive, the digitizer is triggered later. Can be negative, up to the time of the readout pulse trigger.

The Qubit Configuration section includes chassis, slot, and channel numbers for the control AWG of each qubit. I/Q channels are specified as pairs, whereas Z and 'gate' outputs as single channels. If the configuration has an M5302A, one may use a single SMB channel as the 'gate' instead. To do so, simply select 'Use Digital I/O for Gate' for the desired qubit and specify the chassis, slot, and channel to be used. Each qubit is linked to a readout bus.

The Readout Configuration includes a similar configuration for the readout AWGs, as well as the corresponding Digitizer number. This assignment will trigger acquisition on all four channels of an M3102A digitizer. The configuration of that digitizer can be controlled in Labber with a Keysight PXI Digitizer or Keysight PXI Digitizer Demod driver.

Note: When using the Keysight PXI Sequencer in the Labber Measurement program, **Arm/trig mode** must always be enabled with the trigger channel set to Keysight PXI Sequencer - Output.






Fig. B3.4.5: Mapping qubits and modules in the Keysight PXI Sequencer. The middle screen capture shows the ability to use the M5302 as the 'Gate'.

B.3.4.1.3 Keysight PXI Agile AWG

See [Keysight PXI Agile AWG](#)

B.3.4.1.4 Sequence limitations

Delays. The maximum static delay between pulses is 83.8 ms for a M3202 AWG. The maximum delay between pulses and the maximum pulse length is 27.9 ms for a M5302 Digital I/O. The limitation on delays (but not on pulse length) does not apply to sequences with Dynamic duration or with user-defined pulses with dynamic timing (see Section B.3.4.6.2). If a delay or length exceeds the maximum value, an error of the following type is displayed:

Delay [...] exceeds the maximum of XXX ns for static delays [...]

Waveforms. The maximum number of unique waveforms per AWG channel is 64. When using multiple tones per channel, each tone has independent waveform memory. The exception is when using more than 2 tones per channel, in which case the number of unique waveforms is reduced to 8 per tone. If your sequence is exceeding the maximum allowed number of waveforms, make sure "Apply amplitude/phase to waveform primitives" is set to False. See [Multi-Qubit Pulse Generator](#) for more detail.

Number of pulses. The maximum number of compiled instructions (including those that trigger pulses) is 1024 per module. This is an upper bound on the number of pulses in a sequence, as some of the instructions are used to set the parameters that are not pre-compiled in the waveforms (see Section B.3.4.1).

B.3.4.2 Set up driver connections in a Labber measurement

To set up a measurement with the Agile mode, the following channels connections must be established (see Fig. 2.6.7):

- MQPG - Pulse Sequence — Keysight PXI Sequencer - Pulse Sequence
 - This connection carries the pulse order, timing, and runtime-parameter values like amplitude and phase.
- MQPG - Waveform primitives — Agile AWG - Input waveforms
 - This connection carries the pulse primitives, which are uploaded to hardware by the Agile AWG driver
- Keysight PXI Sequencer - PXI config — Agile AWG - PXI config
 - This connection carries the mapping of qubits to physical control and readout channels, so the Agile AWG driver uploads the correct waveforms.

Note that all the AWG drivers are connected to the same MQPG and PXI Sequencer outputs.




Fig. B3.4.6: Schematics of the relations between the Labber Drivers at play in the Agile Sequencing mode.

B.3.4.3 Hardware-loop dynamic sweeps

Regardless of which quantities are selected as dynamic (see Agile sequencing section above), a fast sweep is always set up using the hardware-loop mode. The trigger channel must be Keysight PXI Sequencer - Output. In the following example, we set up the fast sweep of the pulse amplitude for Qubit#1, following the configuration of amplitude as a dynamic property as shown in Figure 3. Note that, in contrast to the non-Agile sequence mode, in this case only a single waveform is uploaded on the AWG for that pulse, with its amplitude cycling from 0 to 1 V with a 10 mV step for every loop iteration.





Fig. 2.6.7: Measurement Setup for a dynamic amplitude hardware-loop sweep

B.3.4.4 2D sweeps using hardware-dynamic variables

An efficient way to sweep multiple variables in a gate sequence is to combine the hardware-loop sweep of one of these variables (inner loop) with additional sweeps for the other ones (outer loops). If all these variables are set as hardware-dynamic, the corresponding pulse parameters are updated during the measurement execution without the overhead introduced by uploading new waveforms or compiling a new HVI sequence. To set up such a scan, first check the Limit hardware looping to first item option in the Show Settings panel. Then add the desired sweep parameters as usual in the Step sequence, with the top sweep being the innermost loop (Last pulse phase in the example below).



B.3.4.5 Multi-tone control with shared AWG channels

If multiple qubits are associated to the same I/Q pair, there are 2 different behaviors, depending on whether the channels are used for control or readout:

- **Multiplexed readout.** If multiple qubits share the same Readout bus, the measurement waveforms are combined into a single one by the Agile AWG driver. To enable this mode, the three Apply frequency/phase/amplitude modulation to readout primitives options need to be enabled. Currently, this mode supports only one of the following two possibilities in the same sequence: 1) simultaneous readout of all qubits sharing the same bus; 2) readout of a single qubit at a time.
- **Multiplexed control.** If multiple qubits share the same control channels, the PXI Sequencer allows for the coherent control of those qubits, either by alternating between their control frequencies, or by simultaneously applying any combinations thereof. Differently from multiplexed readout, phase is maintained when switching between frequencies by use of multiple oscillators inside the AWG, and it is still possible to set all the pulse parameters dynamically (see above). Note the following limitations:
 - The maximum number of multiplexed control frequencies is 4 for channels 1/2 and 2 for channels 3/4.
 - If 1 or 2 qubits are assigned to Ch 1/2, there are no modifications to the Sequencer behavior (with the exception in bold below)
 - If 3 or 4 qubits are assigned to Ch 1/2, the Sequencer will introduce additional instructions to switch between the control of different qubits. This may require increasing the start-start distance between pulses, as indicated by a corresponding error message of the Sequencer driver. Furthermore, the maximum number of waveform primitives on Ch 1/2 is then 8 (instead of 64).

B.3.4.6 Creating a custom sequence

This guide can be used to create an Agile sequence (Multi-Qubit Pulse Generator for Agile sequences) that goes beyond the built-in sequences in the MQPG driver. By selecting the Custom sequence type in the MQPG (see Fig. 2.6.1), one can point to a file containing the desired sequence. Here, we explain how to create such a sequence and provide some examples making use of the dynamic features of the sequencer. A few examples are also provided with the Keysight PXI Sequencer driver (in the sequence directory).




Fig.2.6.1: Interface to select a custom sequence in the Multi-Qubit Pulse Generator.

B.3.4.6.1 How to define hardware-dynamic parameters

The first step in creating a sequence with arbitrary dynamic parameters is to identify which parameters can be changed in hardware. This is done using the variables `dynamic_QUANT`, where `QUANT` can be one of: amplitude, frequency, phase or timing. For example, let's define a Rabi sequence containing a pulse with variable amplitude. Such a sequence can be defined by creating a `.py` file containing the following object:

```

class CustomSequence(Sequence):
    def generate_sequence(self, config):
        """Generate sequence by adding gates/pulses to waveforms."""
        # get parameters
        amp = config['Parameter #1'] # in V

        # define the variable-amplitude pulse
        pulse_var = RabiGate(amplitude = amp, dynamic_amplitude = 1)

        # add some gates
        # ...
        # add a gate to qubit 0 with a variable amplitude pulse
        self.add_gate(0, pulse_var, dt = delay)
        # add other gates
        # ...

```

In the defined pulse_var, the amplitude amp is set through the MQPG quantity Parameter #1. **Note:** this requires the option Apply amplitude modulation to qubit primitives in the Primitives section to be unchecked.

Additionally, setting the pulse property dynamic_amplitude = 1 (or any other integer, see below) allows for the amplitude to be modified at runtime using the hardware loop mode (see [Multi-Qubit Pulse Generator for Agile sequences](#)).

The same pulse pulse_var can be used again in the same sequence, either on the same qubit or on a different one. In this case, the amplitude of all those gates will be modified together and by the same amount, as they are all linked to the same variable dynamic_amplitude = 1.

To create a sequence with multiple, distinct amplitudes, one needs to define a second pulse, with a different variable for both amplitude (referring to a different Parameter #) and dynamic_amplitude. Whereas, only 1D sweeps are currently supported in hardware-loop mode, the two variables can be linked to be linearly swept together as described in the Labber docs.

Other quantities that can be similarly looped through are frequency and phase.

Warning: switching a channel to a second frequency and then back to the first does not preserve phase relations between pulses at the first frequency.

Note that, differently from built-in sequences, custom sequences with dynamic parameters do not require selecting the corresponding checkboxes (such as, dynamic timing, enable hardware-dynamic amplitudes, and so on). Those parameters must explicitly be defined as dynamic in the sequences themselves (see examples in this page and in the sequences directory).

B.3.4.6.2 Dynamic timing

Differently from the quantities above, which are properties of the pulse itself, the pulse timing is defined in relation to the other pulses in a sequence. There are two options to specify the timing of a pulse, either the relative time dt or the absolute time $t0$.

```
delay = config['Parameter #1'] # in s

pulse_var = SingleQubitXYRotation(dynamic_timing = 1)

self.add_gate(0, X2p)
self.add_gate(0, pulse_var, dt = delay)
self.add_gate(0, X2p, dt = 1e-6)
```

Similar to the description above, we define a `dynamic_timing` value in the pulse to make the hardware-loop sweep possible. However, the timing value is only defined later when the pulse is added to the sequence. In this case, the second pulse will be delayed compared to the first pulse by a value equal to `delay` (which must be positive).

In the above example, the third pulse will track the second and always start 1 microsecond after it, even if it doesn't have a `dynamic_timing` value defined. This is because the pulse timings are compiled by the Keysight PXI Sequencer into the start-start distance between consecutive pulses, which is in this case unchanged for the third pulse. If instead we want to keep the position of the third pulse fixed, we need to make its relative timing variable compared to the previous pulse:

```
delay = config['Parameter #1'] # in s

pulse_var1 = SingleQubitXYRotation(dynamic_timing = 1)
pulse_var2 = SingleQubitXYRotation(dynamic_timing = 2)

self.add_gate(0, X2p)
self.add_gate(0, pulse_var1, dt = delay)
self.add_gate(0, pulse_var2, dt = 1e-6 - delay)
```

Since the relative delays for pulse 2 and 3 are different, they each need to be defined with a distinct value of `dynamic_timing`.

Even in cases where the distance between pulses is nominally the same and repeated in a sequence (for example, in a CP/CPMG sequence with variable duration), distinct values of `dynamic_timing` may still be required for each pulse. This is due to the fact that the timing of the compiled instructions may differ depending on its position in the sequence. If this is the case, the PXI Sequencer driver will throw an error of the type:

"Multiply-defined register value `dynamic_timing_1` on `awg[...]`"

This issue can be solved by assigning different values of `dynamic_timing` for each of those pulses, as in the example above.

It is not possible to dynamically sweep the timing of one pulse with respect to another if: a) the two pulses are emitted by the same AWG; *and* b) the beginning of the two pulses coincides for one of chosen timings. This is because generating simultaneous pulses on the same AWG requires the compiler to condense two trigger instructions into one, which cannot be done dynamically.

B.3.4.6.3 Multiple dynamic quantities

For each quantity, the value of the `dynamic_QUANT` variable is arbitrary, but pulses referring to the same variable **and for the same QUANT** use the same register and will be modified together. For example, in the script below,


```
# get parameters
amp = config['Parameter #1'] # in V
freq = config['Parameter #2'] # in Hz
phase = config['Parameter #3'] # in rad
last_phase = config['Parameter #4'] # in rad

pi2_first = SingleQubitXYRotation(phase = phase, dynamic_phase = 1)
pulse_var = RabiGate(phase = phase, amplitude = amp, frequency = freq,
dynamic_phase = 1, dynamic_amplitude = 1, dynamic_frequency = 1)
pi2_last = SingleQubitXYRotation(phase = last_phase, dynamic_phase = 2)

self.add_gate(0, pi2_first)
self.add_gate(0, pulse_var, dt = delay)
self.add_gate(0, pi2_last)
```

there is no relation between the different quantities for pulse_var (phase, amplitude, frequency), but the phase for pulse_var and pi2_first are linked together. An independent phase can be set for pi2_last by referring to a second dynamic variable $\text{dynamic_phase} = 2$.

B.3.4.7 Using an External Gate



Example of using an external gate in the Multi-Qubit Pulse Generator from the M5302A

To use a 'gate' channel from the M5302 rather than AWG channel do the following: Select 'Generate waveform primitives' in the 'Primitives' section as well as select 'External gate' in the 'Output' section of the Multi-Qubit Pulse Generator. The remaining logic of the gate following the XY pulse operates the same as if it were coming from the same AWG.

B.3.4.8 Intermediate measurements

By default, a readout pulse is added at the end of every custom sequence. In addition, one can insert readout pulses at any point during the sequence by adding one or more ReadoutGate objects at the specified times. When the sequence is executed, the digitizer will be triggered at the start of each readout pulse (with an offset equal to the Digitizer delay on the PXI Sequencer front panel). For proper averaging and visualization of the data, the Number of Records on the PXI Digitizer or PXI Digitizer Demod driver must be set equal to the number of measurements in the sequence.

B.3.4.9 Creating Custom Multi-Qubit Pulse Generator and Sequencer Drivers

The Multi-Qubit Pulse Generator driver is configured for 16 qubits by default. However, it can be overwritten to drive an arbitrary number of qubits between 1 and 64 in

waveform sequencing mode (see definition). The default nearest neighbor lattice for two qubit gate pairs can also be overwritten. Finally, the number of readout channels allowed has been expanded from two to the number of qubits in the driver, up to 64. This guide describes how to create a custom Multi-Qubit_Pulse Generator driver to support an arbitrary number of qubits for an arbitrary set of two qubit gate pairs. To use a custom MQPG with more than 16 qubits for agile sequencing a custom Sequencer must also be generated (see Section 3.10).


N.B. Custom drivers should be created in the Local Drivers directory (by default C:\Users\<user_name>\Labber\Drivers), rather than overwriting the ones in the installation directory (by default C:\Program Files\Keysight\Labber\Drivers). This is to avoid Labber updates overwriting custom drivers, and because writing files in the installation directory typically requires administrator permissions. To create modified versions of the Multi-Qubit Pulse Generator and Keysight PXI Sequencer drivers, first copy the top-level folder of each, including all Python files, into the local Drivers directory.

B.3.4.9.1 *Creating a new Multi-Qubit Pulse Generator driver*

In the **MultiQubit_PulseGenerator** driver folder, there is a file called **generate_mqpg.py**. In order to create a new driver .ini file this file can simply be called from the command line. The number of qubits in the driver can be specified with the **n_qubits** input argument. If no **n_qubits** is specified, the default is 16 qubits. For example, to create a MultiQubit_PulseGenerator driver for 64 qubits cd to **drivers/MultiQubit_PulseGenerator** and run:

```
python generate_mqpg.py --n_qubits 64
```

This will replace your local **MultiQubit_PulseGenerator.ini** file with a new one. The newly created driver can be added as usual from the Instrument Server and run as usual from the Measurement Editor or the Labber API.



B.3.4.9.2 Adding a custom list of two-qubit gates

By default, the Multi-Qubit Pulse Generator driver assumes nearest-neighbor connectivity between qubits in a line for two-qubit gates. For each qubit n , gates are possible with qubit $n+1$ and qubit $n-1$. This can be overwritten with a custom connectivity by supplying a json file consisting of a list of two-qubit pairs labeled "gates". For example, for 4 qubits with all-to-all connectivity the file would look like this:


connected_qubits.json

```
{"gates": [[1, 2], [1, 3], [2, 3], [2, 4], [3, 4], [4, 1]]}
```

The path to this file can then be supplied as the command line argument **gates_file** for **write_ini.py**. Creating a Multi-Qubit Pulse Generator driver for 4 qubits with the above connectivity would be done by running:

```
python generate_mqpg.py --n_qubits 4 --gates_file connected_qubits.json
```

Once this is run, the new MultiQubit Pulse Generator can be added to the Instrument Server and will have the two-qubit gate pairs specified in **connected_qubits.json**:



B.3.4.9.3 Creating a new Sequencer Driver

In the folder containing the Keysight PXI Sequencer driver, there is a file named *generate_sequencer.py*. In order to create a new sequencer driver, this file can be called from the command line. To do this, navigate to the Keysight PXI Sequencer driver folder and run

```
python generate_sequencer.py
```

There are 3 possible input arguments:

n_qubits (Integer): required, maximum 64. Specifies the number of qubits the driver will contain

n_chassis (Integer): optional, default 1, maximum 4. Specifies the number of chassis the driver will contain.

mqqpg (Boolean): optional, default False. If true, overwrites the existing .ini file for the MultiQubit Pulse Generator with one containing the number of qubits as the new sequencer.

(Note: A custom MultiQubit Pulse Generator driver can also be created from the file `generate_mqqpg.py` in the MultiQubit Pulse Generator driver folder. If a custom connectivity for two-qubit gates is desired the MultiQubit Pulse Generator driver must be created using `generate_mqqpg.py`, not `generate_sequencer.py`. Documentation on how to do that is [here](#).)

For example, to create a sequencer to control 64 qubits on 4 chassis with a matching MultiQubit Pulse Generator driver, the syntax would be:


```
python generate_sequencer.py --n_qubits 64 --n_chassis 4 --mqqpg True
```

Once this is run, the new sequencer can be added in the Instrument Server and run as normal.

B.3.5 Keysight PXI Agile AWG

The Keysight PXI Agile AWG driver uses the Quantum IP Library Gapless Agile AWG to play pulse sequences from waveform primitives. It is used in conjunction with a pulse sequence generator driver (like the Multi-Qubit Pulse Generator) which supplies Waveform-Primitive-type data, and with the Keysight PXI Sequencer, which provides the mapping between logical qubit primitives and physical AWG channels.

The driver allows AWG channels to be operated individually or in IQ pairs. In individual mode, each channel has a distinct amplitude scale and phase. In IQ mode, channels can be grouped into pairs (1/2 and 3/4), and a pair has a single-phase offset, as well as an amplitude ratio and phase skew (referenced to 90 degrees) for IQ mixer imbalance. These settings are implemented by setting amplitude and phase registers in a DUC FPGA block, and are applied to all tones on the same channel. In either mode, channels have individual DC offsets, which are persistent. As shown in the interface below, the two IQ pairs can be grouped or ungrouped independently, allowing for one IQ pair and two independent channels:



B.3.5.1 Requirements

This driver requires that the Quantum Library Python API (PyQuLibrary) be installed in the Python environment used by Labber drivers. For more information, see [Configuring a Python Environment for Labber Drivers](#). For information on usage, see Multi-Qubit Pulse Generator for Agile sequences.

B.3.5.2 Instrument settings

In individual or IQ mode, all channels n have the following setting available

- **Ch n DC Offset [V]: DOUBLE**
Static DC offset. Ranges from -1.5 to +1.5 V

In individual mode, each channel n has the following output settings:

- **Ch n Amplitude Scale: DOUBLE**
Scale factor used to adjust the amplitude of input waveforms.
- **Ch n Phase [deg]: COMBO**
Additional phase offset applied to all waveforms.

In IQ mode, each pair $n/n+1$ has the following output settings, for $n=1, 3$:

- **Ch n/n+1 IQ Amplitude Scale: DOUBLE**
Scale factor used to adjust the amplitude of input waveforms. Applied to both channels.
- **Ch n/n+1 IQ Amplitude Ratio: DOUBLE**
Scale factor used to adjust the amplitude ratio of input waveforms. If <1 then channel n+1 is multiplied by this factor, if >1 then channel n is divided by this factor.
- **Ch n IQ Phase Skew [deg]: COMBO**
Differential phase shift applied to channel n+1. This is referenced to 90 degrees, so a value of zero will result in signals in quadrature.
- **Ch n IQ Phase Offset [deg]: COMBO**
Additional phase offset applied to all waveforms on channels n and n+1.

B.3.5.3 Input channels

- **Input Waveforms: VECTOR**
Waveforms to be uploaded to AWGs. Drivers supplying this data should use the WaveformPrimitives class.
- **PXI Config: VECTOR**
System mapping from qubit control and readout to physical AWG channels. Drivers supplying this data should use the PXIRouting class.


B.3.5.4 Use with HVI sequences

Since the Agile AWG driver controls channel amplitudes and phases by modifying the waveform memory, these parameters are in addition to any controls performed in real-time with an HVI sequence. The waveform modulation frequencies are exclusively set in the Multi-Qubit Pulse Generator driver. Pulse-dependent amplitudes and phases are applied in addition to the Agile AWG settings. See Multi-Qubit Pulse Generator for Agile sequences for more details on the interworking of these drivers.

B.3.5.5 IQ Sideband Modulation

The FPGA design which defines the Agile AWG uses a combination of phase and amplitude modulation with an AWG which outputs phase and magnitude data, and a digital up-converter (DUC) which specifies a frequency, phase, and modulation gain.

The layout of the FPGA includes a Dual AWG which outputs an amplitude and phase, which are used to modulate a function generator:



In total the modulated output of the system is

$$(A + G_m \times MAG(t)) \times \cos(2\pi ft + \phi + G_p \times PHS(t))$$

Where

$MAG(t)$ and $PHS(t)$ are the GAWG outputs defined on $(-1,1)$. They are the magnitude and phase of the uploaded complex waveform, scaled appropriately.

G_m is the mod gain which is set to 1.5 V,

G_p is the phase gain which is set to 180 degrees, and

$A = 0$ (this would otherwise output a constant signal at f)

Consider a complex waveform $1.5 \text{ V} \times (I(t) + j \times Q(t))$ with which we want to drive an IQ mixer using two output channels.

We set up the AWGs as follows:

AWG1:

$$z_1(t) = I(t) + j \times Q(t)$$

AWG2:

$$z_2(t) = r \times (I(t) + j \times Q(t)) \times \exp(j\pi \times (90 + \phi_{sk})/180)$$

using the instrument settings

ϕ_{sk} is the phase skew in degrees, and

r is the amplitude imbalance ratio

For convenience, we define the mag and phase of the input waveform $I(t) + j \times Q(t)$ as

$$|z(t)| = |I(t) + j \times Q(t)|$$

$$\phi(t) = \arg(I(t) + j \times Q(t)) \text{ in degrees}$$

For a baseband pulse, $\phi(t)$ is constant.

For a baseband pulse with DRAG, $\phi(t)$ is some small variation on constant.

For a detuned pulse, $\phi(t)$ is linear in time.

Then the waveforms uploaded and played are

$$MAG_1(t) = |z(t)|$$

$$PHS_1(t) = \phi(t) / 180$$

$$MAG_2(t) = r \times |z(t)|$$

$$PHS_2(t) = (\phi(t) + 90 + \phi_{sk}) / 180$$

The DUCs are both configured with the same frequency f and the same phase ϕ_0 .

The output of DAC channels 1 and 2 are then:

$$V_1(t) = G_m \times |z(t)| \times \cos(2\pi ft + \phi_0 + \phi(t))$$

$$V_2(t) = r \times G_m \times |z(t)| \times \sin(2\pi ft + \phi_0 + \phi(t) + \phi_{sk})$$

For a perfect IQ mixer we would set $r = 1$ and $\phi_{sk} = 0$, and have

$$V_1(t) = G_m \times |z(t)| \times \cos(2\pi ft + \phi_0 + \phi(t))$$

$$V_2(t) = G_m \times |z(t)| \times \sin(2\pi ft + \phi_0 + \phi(t))$$

B.3.6 Keysight PXI SMU


B.3.6.1.1 Summary

This driver controls the Keysight M9614A or M9615A SMUs, which are 5-channel high-precision source/measure units. This driver supports simple voltage and current sourcing, I/V measurements, and resistance measurements. All source and measurement functions can be controlled independently for each of the 5 channels.

B.3.6.1.2 Python IVI-C Wrapper


This instrument does not have support for typical VISA/SCPI communication, and thus the structure of the driver is somewhat different from typical Labber drivers. Translation between python and IVI-C is currently encapsulated in the pyviktm960x package, which is included in the driver folder. This package must be findable in order for the driver to work.

B.3.6.1.3 Source Functions



- **Output State: BOOLEAN**
Enables or disables the output on a given channel.
- **Source Function: COMBO**
Switches the output between sourcing current and sourcing voltage.
- **Voltage Amplitude: DOUBLE**
Sets the DC output voltage amplitude in Volts. Limit $\pm 30V$.
- **Current Amplitude: DOUBLE**
Sets the DC output current in Amps. Limit $\pm 500mA$.
- **Sweep Rate: DOUBLE**
Sets the maximum sweep rate for current or voltage, in V/s or A/s, respectively.

B.3.6.1.4 Measurement Functions:



- **Measurement Type: COMBO**
Determines whether the channel is measuring current, voltage, or resistance.
- **Acquisition Duration: DOUBLE**
The duration of each acquisition in seconds.
- **Measured Voltage: DOUBLE**
The measured voltage in Volts. (Read-only)
- **Measured Current: DOUBLE**
The measured current in Amps. (Read-only)
- **Measured Resistance: DOUBLE**
The measured resistance in Ohms. If the channel is sourcing voltage, this is calculated by the source voltage by the measured current. If the channel is sourcing current this is calculated by dividing the measured voltage by the source current. (Read-only)

B.3.6.1.5 Sample Code

This code will sweep the voltage on Channel 1 from -1 to 1 Volt in steps of 0.1V, measuring and printing the resistance at each step

```

import Labber
import numpy as np

"""Replace 'localhost' and the instrument address with the values for your setup. The
instrument server must be running
and the SMU must be added for this to connect"""
client = Labber.connectToServer('localhost')
smu = client.connectToInstrument('Keysight PXI SMU', dict(interface='PXI',
address='6'))
smu.startInstrument()

# Configure measurement on Channel 1
smu.setValue('Ch1 Output State', 'On')
smu.setValue('Ch1 Source Function', 'Voltage')
smu.setValue('Ch1 Measurement Type', 'Resistance')

# Ramp voltage from -1 to 1V in steps of 0.1V, printing measured resistance at each
step
for v_amp in np.arange(-1.0, 1.1, 0.1):
    smu.setValue('Ch1 Voltage Amplitude', v_amp)
    res = smu.getValue('Ch1 Measured Resistance')
    print(res)

smu.stopInstrument()
client.close()

```